



PKU Test Agar Base

M282

PKU Test Agar Base is recommended for estimation of phenylalanine in blood for detection of Phenylketonuria (PKU).

Composition**

Ingredients	Gms / Litre
L-Glutamic acid	0.500
DL-Alanine	0.500
Asparagine	0.500
Dextrose	10.000
Dipotassium phosphate	15.000
Monopotassium phosphate	5.000
Ammonium chloride	2.500
Ammonium nitrate	0.500
Sodium sulphate	0.500
Magnesium sulphate	0.050
Manganese chloride	0.005
Ferric chloride	0.005
Calcium chloride	0.0025
Agar	15.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 50.06 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE OR OVERHEAT. Cool to 50°C and add *Bacillus subtilis* spores and 3.3 mg β-2-thienylalanine. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Phenylketonuria is a congenital defect caused due to absence of phenylalanine hydroxylase. As a result of this, phenylalanine accumulates in the blood, which is excreted via urine hence it is called as phenylketonuria. Subsequently this deficiency may cause brain damage resulting in mental retardation. Guthrie and Tiekemann (3) devised a modified inhibition assay for early detection of PKU using blood / urine samples of newborn infants having low levels of phenylalanine by determining the serumphenylalanine levels or the level of phenylpyruvic acid in urine. PKU Test Agar Base, developed by Demain (1) is a chemically defined medium, which is supplemented with β-2-thienylalanine. It provides the agar bed for testing blood samples.

The Guthrie test (1-4) was developed on the observation that *Bacillus subtilis* is normally inhibited in presence of b-2-thienylalanine but grows well when L-phenylalanine is added to the medium. Phenylalanine neutralizes the b-2-thienylalanine and allows bacteria to grow. The phenylalanine level can be read to determine the level of amino acid in blood. Other than phenylalanine, proline, phenylpyruvic acid or phenyllactic acid can be used. Small filter paper discs saturated with patients blood are placed on PKU Test Agar with b-2-thienylalanine inoculated with *Bacillus subtilis*. Control discs impregnated with different levels such as 2, 4, 6, 8,10,12 and 20 mg% of L-phenylalanine are also placed on the medium. After overnight incubation, zones of growth around the paper discs are observed and compared with zones around control discs. A response comparable to 4 mg% control disc is considered as presumptive positive. The results can be repeated using a duplicate test disc and a chemical or spectrofluorometric procedure (5, 6).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light yellow coloured clear to slightly opalescent gel forms in Petriplates

Reaction

Reaction of 5% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

M282: Cultural characteristics observed with added B-2-thienylalanine, after an incubation at 35-37°C for 12-16 hours

Organism	Growth w/ 2% Phenylalanine	Growth w/ 4% Phenylalanine	Growth w/ 6% Phenylalanine	Growth w/ 8% Phenylalanine	Growth w/ 10% Phenylalanine	Growth w/ 12% Phenylalanine
----------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------	-----------------------------

Cultural Response

Bacillus subtilis ATCC 6633 none-poor luxuriant luxuriant luxuriant luxuriant luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry date on the label.

Reference

1. Demain A. L., 1958, J. Bacteriol., 75:517.
2. Guthrie R., 1961, J. Am. Med. Assoc., 178:863.
3. Guthrie R. and Tiekemann H., 1960, London Conference on the Scientific study of Mental Deficiency, London.
4. Guthrie R. and Susi A., 1963, Pediatrics, 32:338.
5. Ambrose J. A., Ingerson A., Gorretson L. G., Chung L. W., 1967, Clin. Chem. Acta., 15:493.
6. Ambrose J. A., 1969, Clin. Chem., 15:15.

Revision : 1 / 2011

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.