

## MB580 HiPurA® Forensic Multi Sample DNA Purification Kit

### Kit Contents

Product Code	Reagents provided	MB580		
		20 preps	50 Preps	100 Preps
DS0015	Lysis Solution (AL)	48 ml	120 ml	240 ml
DS0010	Lysis Solution (C1)	14 ml	35 ml	70 ml
DS0011	Prewash Solution Concentrate (PW)	8 ml	20 ml	40 ml
DS0012	Wash Solution Concentrate (WS)	4 ml	10 ml	20 ml
DS0040	Elution Buffer (ET) [10 mM Tris-Cl, pH 8.5]	1.6 ml	4 ml	8 ml
DS2280	Proteinase K	10 mg	25 mg	50 mg
DS0075	Bone Extraction Buffer	28 ml	70 ml	140 ml
DS0077	Wash Solution Concentrate (BW)	8 ml	20 ml	40 ml
DBCA03	HiElute Miniprep Spin Column (Capped) [in DBCA016 Collection Tube]	20 nos	50 nos	100 nos
DBCA016	Collection Tube (Uncapped), Polypropylene (2.0 ml)	40 nos	100 nos	200 nos
DBCA017	HiPer® Lock Micro- Centrifuge Tube 2ml	40 nos	100 nos	200 nos

### Intended Use

Recommended for isolation of DNA from human blood/animal blood, plasma, serum and other samples such as cells, tissues, stool sample and body fluids.

### Introduction

HiPurA® Forensic Multi Sample DNA Purification Kit provides a fast and easy method for purification of total DNA for reliable applications in PCR and Southern blotting technique. The DNA purification procedure using the miniprep spin columns comprises of three steps viz. adsorption of DNA to the membrane, removal of residual contaminants and elution of pure genomic DNA. HiMedia's HiElute Miniprep Spin Column (Capped) format allows rapid processing of multiple samples. The columns have a high binding capacity and high quality genomic DNA is obtained from various species. The genomic DNA obtained is compatible with downstream applications such as restriction enzyme digestion, PCR and Southern blotting.

### HiPurA® Forensic Multi Sample DNA Purification Kit

This kit simplifies isolation of DNA from cigarette butts, envelopes, hair, nail clippings, stamps, or blood, saliva or semen stains, Bone and tooth samples with spin-column procedure. Following lysis, DNA binds to the silica gel membrane of the HiElute Miniprep Spin Column (Capped) to yield purified DNA. Two rapid wash steps remove trace salt and protein contaminants resulting in the elution of high quality DNA in the Elution Buffer (ET) provided with the kit.

### **HiElute Miniprep Spin Column (Capped) [DBCA03]**

HiElute Miniprep Spin Column (Capped) is based on the advanced silica binding principle presented in a microspin format. The system efficiently couples the reversible nucleic acid-binding properties of the advanced silica gel membrane and the speed plus versatility of spin column technology to yield high quantity of DNA. The use of spin column facilitates the binding, washing and elution steps thus enabling multiple samples to be processed simultaneously. This column eliminates the need for alcohol precipitation, expensive resins, and harmful organic compounds such as phenol and chloroform, otherwise employed in traditional DNA isolation techniques. DNA binds specifically to the advanced silica-gel membrane while contaminants pass through. PCR inhibitors such as divalent cations and proteins are completely removed in two efficient wash steps, leaving pure nucleic acid to be eluted in the buffer provided with the kit.

#### **Elution**

A single elution with 50 µl of Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) will provide sufficient DNA to carry out multiple amplification reaction. Elution with volume less than 50 µl will increase the final DNA concentration, but will reduce the overall DNA yield.

#### **Concentration, yield and purity of DNA**

Spectrophotometric analysis and agarose gel electrophoresis will reveal the concentration and the purity of the genomic DNA. Use Elution Buffer to dilute samples and to calibrate the spectrophotometer, measure the absorbance at 260 nm, 280 nm, and 320 nm using a quartz microcuvette. Absorbance readings at 260 nm should fall between 0.1 and 1.0. The 320 nm absorbance is used to correct for background absorbance. An absorbance of 1.0 at 260 nm corresponds to approximately 50 µg/ml of DNA. The  $A_{260}-A_{320}/A_{280}-A_{320}$  ratio should be 1.6 – 1.9. Purity is determined by calculating the ratio of absorbance at 260 nm to absorbance at 280 nm. DNA purified by HiPurA® Forensic Multi Sample DNA Purification Kit is free of protein and other contaminants that can inhibit PCR or other enzymatic reactions.

Concentration of DNA sample (µg/ml) = 50 x  $A_{260}$  x dilution factor.

#### **Materials needed but not provided**

- Heating Block or thermomixer at 55°C
- 70°C water bath or heating block
- Motor and pestle
- Liquid nitrogen
- Tabletop Microcentrifuge (with rotor for 2.0 ml tubes)
- 1M DTT (Dithiothreitol) solution
- RNase A Solution (20mg/ml) (Product Code: DS0003)
- 1X PBS (ML116)
- Ethanol (96 – 100%)
- Molecular Biology Grade Water (Product code: ML024)

#### **Storage**

Store the HiPurA® Forensic Multi Sample DNA Purification Kit between 15-25°C except certain components as specified on each labels. Under recommended condition kit is stable for 18 months.

#### **General Preparation Instructions**

1. Set the heating block or thermomixer at 55°C.

2. Preheat a water bath or heating block to 70°C.
3. **Thoroughly mix reagents**  
Examine the reagents for precipitation. If any kit reagent forms a precipitate (other than enzymes), warm at 55- 65°C until the precipitate dissolves and allow cooling to room temperature (15-25°C) before use.
4. Ensure that clean & dry tubes and tips are used for the procedure.
5. **Dilute Prewash Solution Concentrate (PW) (DS0011) as follows:**

Number of Preps	Prewash Solution Concentrate (PW)	Ethanol (96-100%)
20	8 ml	12 ml
50	20 ml	30 ml
100	40 ml	60 ml

6. **Dilute Wash Solution Concentrate (WS) (DS0012) as follows:**

Number of Preps	Wash Solution Concentrate (WS)	Ethanol (96-100%)
20	4 ml	12 ml
50	10 ml	30 ml
100	20 ml	60 ml

7. **Dilute Wash Solution Concentrate (BW) (DS0077) as follows:**

Number of Preps	Wash Solution Concentrate (BW)	Ethanol (96-100 %)
20	8 ml	8 ml
50	20 ml	20 ml
100	40 ml	40 ml

8. **Reconstitute DTT (DL-Dithiothreitol) (MB070) (NOT PROVIDED)**

Resuspend the DTT (DL-Dithiothreitol) powder in Molecular Biology Grade Water (ML024) to obtain a 1M stock solution.

Number of Preps	DTT (DL-Dithiothreitol)	Molecular Biology Grade Water
20	30.85 mg	0.2 ml
50	77.12 mg	0.5 ml
100	154.24 mg	1 ml

The product as supplied is stable at room temperature (15-25°C), upon reconstitution store at -20°C in dark (wrapped in an aluminium foil), as mentioned in storage instructions.

**NOTE:** The 1M DTT solution must be added directly to each sample preparation every time.

9. **Reconstitute Proteinase K (DS2280)**

The HiPurA® Forensic Multi Sample DNA Purification Kit contains Proteinase K. Intensive research has shown that it is the optimal enzyme for use with the Lysis Solution provided in the kit. It is completely free of DNase and RNase activity. Proteinase K is the enzyme of choice for use with an SDS containing Lysis Solution. The specific activity of Proteinase K is 33.5 units/mg dry weight.

Resuspend the Proteinase K powder in Molecular Biology Grade Water (ML024) to obtain a 20 mg/ml stock solution.

Number of Preps	Proteinase K	Molecular Biology Grade Water
20	10 mg	0.5 ml
50	25 mg	1.25 ml
100	50 mg	2.5 ml

The product as supplied is stable at room temperature, upon reconstitution store at -20°C as mentioned in storage instructions.

**NOTE:** The Proteinase K solution must be added directly to each sample preparation every time. Do not combine the Proteinase K and Lysis solutions for storage.

### Centrifugation

All centrifugation steps are carried out in conventional laboratory centrifuge e.g. Beckman CS-6KR, Heraeus Varifuge 3.0R, or Sigma 6k10 with fixed angle rotor. The tubes provided with the kit are compatible with almost all laboratory centrifuges and rotors. All centrifugation steps are performed at room temperature and are given in g, the correct rpm can be calculated using the formula:

$$RPM = \sqrt{RCF/1.118} \times 10^5 r$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary g-force.

### Specimen Handling and Collection

Collect whole blood in an anticoagulant tube (an EDTA tube is preferred) under sterile conditions (if to be used for future) and store the samples at 2-8°C for short term storage or -20°C for long term storage.

Collect cells, tissues, Serum, Plasma, stool sample and body fluids in a clean sterile container and store the samples at 2-8°C for short term storage or -20°C for long term storage. Repeated freeze- thaw of samples should be avoided.

Ensure that the blood/ tissue/ cells/stool/bacterial culture/yeast Culture, Serum, Plasma and body fluids sample is at room temperature (15- 25°C) before beginning the protocol. After use, contaminated material must be sterilized by autoclaving before discarding.

### Types of Specimen

Clinical samples: Whole blood, tissues, stool and body fluids.

### Procedure

#### 1. Preparation of Forensic Sample

Prepare the following forensic samples according to the steps given below:

- A. Cigarette butts
- B. Envelopes, stamps, chewing gum or epithelial cells on fabric
- C. Hair roots
- D. Hair shafts without roots
- E. Nail clippings or soft tissue
- F. Material stained with blood, saliva or semen
- G. Bone and tooth
- H. Liquid blood

**A. To lyse cigarette butts**

1. Cut a 1cm<sup>2</sup> piece of outer paper from the end of the cigarette or filter. Cut it further into 6 smaller pieces and transfer the pieces into a 2.0ml capped microcentrifuge tube.

2. **Lysis**

Add 300 µl Lysis Solution AL (DS0015), 20 µl Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) and 20 µl of 1M DTT (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 55°C for at least 1- 1½ hours with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Add 300 µl of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

5. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

6. Centrifuge at 12,000-16,000 x g (≈13,000-16,000 rpm) for 1 minute at room temperature.

7. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 6, into a new 2.0 ml collection tube (not provided) and add 300 µl of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

8. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650 µl of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g (≈8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

9. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

10. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g (≈10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

11. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700 µl of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g (≈13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

12. Centrifuge the column for 2 minutes at 20,000 x g (≈14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

### 13. DNA Elution

Pipette 20 - 50 µl of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at ≥6,500 x g (≈10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

14. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

### B. To lyse envelopes, stamps, chewing gum, epithelial cells on fabric (DNA extraction from sexual assault samples)

1. Cut a 0.5 - 2.5 cm<sup>2</sup> sample from the envelope, stamp, chewing gum or fabric into small pieces. Cut it further into smaller pieces and transfer the pieces into a 2.0ml capped microcentrifuge tube.

**NOTE:** Surface contamination can be reduced by using a swab moistened with distilled water, before cutting the envelope or stamp.


#### 2. Lysis

Add 500 µl Lysis Solution AL (DS0015), 20 µl Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 56°C for at least 1 hour with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Remove the solid particles from the lysate by centrifuging the tube at 14,000 rpm for 5 minutes. Follow the flowchart.



Incubate at 56°C for 1 hour with shaking at 900 rpm in a thermomixer. (Continue with step 5)

5. Add 300 µl of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

6. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

7. Centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature.

8. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 7, into a new 2.0 ml collection tube (not provided) and add 150 µl of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

9. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650 µl of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g ( $\approx$ 8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

10. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

11. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 500 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

12. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add 500 µl of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

13. Add 500 µl of Ethanol (96-100%) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

14. Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

15. **DNA Elution**

Pipette 20 - 50  $\mu$ l of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50  $\mu$ l increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

16. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

**C. To lyse hair roots**

1. Cut a 0.5 - 1.0 cm piece starting from the hair bulb. Transfer the piece into a 2.0ml capped microcentrifuge tube.

2. **Lysis**

Add 300  $\mu$ l Lysis Solution AL (DS0015), 20  $\mu$ l Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) and 20  $\mu$ l of 1M DTT (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 55°C for at least 1- 1½ hours with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Add 300  $\mu$ l of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

5. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

6. Centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature.

7. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 6, into a new 2.0 ml collection tube (not provided) and add 300  $\mu$ l of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

8. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650  $\mu$ l of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g ( $\approx$ 8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

9. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

10. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

11. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700 µl of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

12. Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

13. **DNA Elution**

Pipette 20 - 50 µl of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

14. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

**D. To lyse hair shafts (without roots)**

1. Cut the hair shaft into 0.5-1.0 cm pieces. Transfer the pieces into a 2.0ml capped microcentrifuge tube.

2. **Lysis**

Add 300 µl Lysis Solution AL (DS0015), 20 µl Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) and 20 µl of 1M DTT (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 55°C for at least 1- 1½ hours with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Add 300 µl of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

5. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

6. Centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature.

7. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 6, into a new 2.0 ml collection tube (not provided) and add 300  $\mu$ l of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

8. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650  $\mu$ l of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g ( $\approx$ 8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

9. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

10. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700  $\mu$ l of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

11. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700  $\mu$ l of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

12. Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

13. **DNA Elution**

Pipette 20 - 50  $\mu$ l of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50  $\mu$ l increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

14. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

**E. To lyse nail clippings or soft tissue**

1. Clip the nails into small pieces and chop soft tissues (approximately upto 25mg) to fine pieces. Transfer the pieces into a 2.0ml capped microcentrifuge tube.

2. **Lysis**

Add 300 µl Lysis Solution AL (DS0015), 20 µl Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) and 20 µl of 1M DTT (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 55°C for at least 1- 1½ hours with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Add 300 µl of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

5. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

6. Centrifuge at 12,000-16,000 x g (≈13,000-16,000 rpm) for 1 minute at room temperature.

7. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 6, into a new 2.0 ml collection tube (not provided) and add 300 µl of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

8. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650 µl of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g (≈8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

9. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

10. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g (≈10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

11. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700 µl of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

12. Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

13. **DNA Elution**

Pipette 20 - 50 µl of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

14. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

**F. To lyse material stained with blood, saliva or semen/ buccal swabs**

1. For material stained with blood, saliva or semen- Take 0.5cm<sup>2</sup> area of the stained material and cut into smaller pieces. Transfer the pieces into a microcentrifuge tube. (Continue with Step 2)

For buccal swabs- Cut buccal swab of 0.5cm<sup>2</sup> and transfer the pieces into a microcentrifuge tube. (Continue with Step 2)

2. **Lysis**

Add 300 µl Lysis Solution AL (DS0015), 20 µl Proteinase K (20 mg/ml) (**Refer to General preparation Instructions**) and 20 µl of 1M DTT (**Refer to General preparation Instructions**) to the material in the collection tube. Mix thoroughly by pulse-vortexing for 10-15 seconds.

3. Incubate at 55°C for at least 1- 1½ hours with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 10 minutes for 10-15 seconds to improve lysis of the sample material.

4. Add 300 µl of Lysis Solution C1 (DS0010) and mix thoroughly by pulse-vortexing for 10-15 seconds.

**NOTE:** To ensure efficient lysis, mix the sample and Lysis Solution C1 thoroughly to form a homogenous solution. A white precipitate may form when Lysis Solution C1 is added to Lysis Solution AL. This precipitate will dissolve during incubation in Step 5.

5. Incubate at 70°C for 10 minutes with shaking at 900 rpm in a thermomixer.

**NOTE:** If using a heating block or water bath, vortex the tube after every 3 minutes for 10-15 seconds to improve lysis of the sample material.

6. Centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature.

7. **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 6, into a new 2.0 ml collection tube (not provided) and add 300 µl of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

8. **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650 µl of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g (≈8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

9. Repeat step 8 with the remaining sample. Discard the flow-through liquid.

10. **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g (≈10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

11. **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700 µl of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g (≈13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

12. Centrifuge the column for 2 minutes at 20,000 x g (≈14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

13. **DNA Elution**

Pipette 20 - 50 µl of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at ≥6,500 x g (≈10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

14. Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

**G. To lyse bone and tooth**

**DAY 1**

**Procedure for Decalcification of Bone and Tooth samples**

- Crush the bone/ tooth into small fragments. Grind it to a fine powder using liquid nitrogen. If necessary use a freezer mill.

**NOTE:** It is important to grind the bone fragments into a fine powder, since, finer the powder more is the DNA released. However, if you use a freezer mill do not overgrind, as this may fragment the DNA.

- Weigh 500 mg of bone/ tooth sample and transfer it to a sterile 2.0 ml Collection tube (provided).

**NOTE:** The bone/ tooth sample powder can be stored at room temperature (15-25°C), but we recommend to perform the extraction as soon as possible.

- Add 1.2ml of Bone Extraction Buffer (DS0075) to the sample to decalcify the bone/tooth
- Incubate the sample at 55°C for 16 hours or overnight at approximately 100rpm.

#### **DAY 2**

- Centrifuge the decalcified sample for 2 minutes at 6,500 x g ( $\approx$ 10,000 rpm) at room temperature (15-25°C).

- **Prepare for Binding**

Transfer the supernatant carefully, obtained from step 1, into a new 2.0 ml collection tube (not provided) and add equal volume of ethanol (96- 100%) for preparation of the lysate for binding. Mix thoroughly by vortexing for 5-10 seconds.

- **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650  $\mu$ l of the lysate obtained in the above step, including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6,500 x g ( $\approx$ 10,000 rpm) at room temperature (15-25°C). Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

- Repeat the above step with the remaining sample. Discard the flow-through liquid.

- **Wash**

**(Prepare the Wash Solution (BW) (DS0077) as indicated in General Preparation Instructions)**

Add 500  $\mu$ l of diluted Wash Solution (BW) to the HiElute Miniprep Spin Column (Capped) and centrifuge for 1 minute at 6,500 x g ( $\approx$ 10,000 rpm) at room temperature (15-25°C). Discard the flow-through liquid and reuse the same collection tube.

- Repeat the wash step by adding another 500  $\mu$ l of diluted Wash Solution (BW) (DS0077) to the HiElute Miniprep Spin Column (Capped) and centrifuge for 1 minute at 6,500 x g ( $\approx$ 10,000 rpm) at room temperature (15-25°C). Discard the flow-through liquid and reuse the same collection tube.

- Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

- **DNA Elution**

Pipette 20-50  $\mu$ l of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

- Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

#### H. To lyse liquid blood

- ❖ Add 20 µl of the reconstituted Proteinase K solution (20 mg/ml) (**Refer to General Preparation Instructions**) into 2.0 ml capped collection tube containing 200 µl of the whole blood. Vortex for 10-15 seconds to ensure thorough mixing.

**NOTE:** If the sample is less than 200 µl, add 1X PBS (ML116) to bring the volume upto 200 µl.

##### **Optional RNase A treatment**

If RNA-free genomic DNA is required, add 20 µl of RNase A solution (20 mg/ml) (DS0003). Vortex for 10-15 seconds and incubate for 2 minutes at room temperature (15-25°C).

- ❖ **Lysis reaction**

Add 200 µl of the Lysis Solution (C1) (DS0010) to the sample, vortex thoroughly for a few seconds to obtain a homogenous mixture. Incubate at 55°C for 10 minutes.

**NOTE:** If cell clumps are visible, the sample can be mixed gently by pipetting to obtain a homogenous mixture.

- ❖ **Prepare for Binding**

Add 200 µl of ethanol (96-100%) to the lysate obtained from the above step for preparation of lysate for binding to the spin column. Mix thoroughly by gentle pipetting.

**NOTE:** A homogenous solution is essential.

- ❖ **Load lysate in HiElute Miniprep Spin Column (Capped) [DBCA03]**

Add 650 µl of the mixture from including any precipitate, which may have formed, to the HiElute Miniprep Spin Column (Capped). Centrifuge for 1 minute at 6000 x g (~8000 rpm) at room temperature. Discard the flow-through.

**NOTE:** Use a wide bore pipette tip to reduce shearing of the DNA when transferring contents into the column.

- ❖ Repeat the above step with the remaining sample. Discard the flow-through liquid.

- ❖ **Prewash**

**(Prepare the Prewash Solution (PW) (DS0011) as indicated in General Preparation Instructions)**

Place the column in a same 2.0 ml collection tube and add 700 µl of diluted Prewash Solution (PW) to the column. Centrifuge at 6,500 x g (~10,000 rpm) for 1 minute at room temperature (15-25°C). Discard the flow-through liquid and re-use the same collection tube with the column.

- ❖ **Wash**

**(Prepare the Wash Solution (WS) (DS0012) as indicated in General Preparation Instructions)**

Add another 700  $\mu$ l of diluted Wash Solution (WS) to the HiElute Miniprep Spin Column (Capped) and centrifuge at 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) for 1 minute at room temperature. Discard the flow-through liquid and reuse the same collection tube.

- ❖ Centrifuge the column for 2 minutes at 20,000 x g ( $\approx$ 14,000 rpm) to dry the column membrane to remove the traces of residual ethanol, if observed. Place the column in a new uncapped 2.0 ml collection tube.

#### ❖ **DNA Elution**

Pipette 20 - 50  $\mu$ l of the Elution Buffer (ET) (DS0040) or Molecular Biology Grade Water (ML024) directly onto the center of the HiElute Miniprep Spin Column (Capped) membrane without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at  $\geq$ 6,500 x g ( $\approx$ 10,000 rpm) for 1 minute to elute the DNA.

**NOTE:** To increase the elution efficiency, incubate for 5 minutes at room temperature (15-25°C) after adding the Elution Buffer, then centrifuge. Elution with volumes less than 50  $\mu$ l increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

- ❖ Transfer the eluate to a fresh capped 2ml collection tube for longer DNA storage.

#### **DNA concentration**

Precipitate the DNA by adding 0.1 volumes of 3.0 M sodium acetate, pH 5.2, and 0.7 volumes of isopropanol to the recovered DNA. Mix well by inversion and centrifuge at 10,000 rpm at 4°C for 30 minutes. Decant the supernatant, taking care not to disturb the pellet. Rinse the pellet with 1.5 ml of 70% ethanol and centrifuge at 10,000 rpm at 4°C for 10 minutes. Carefully decant the supernatant and air-dry the pellet until the residual ethanol has evaporated. Resuspend the DNA pellet in 20-50  $\mu$ l (or desired volume) of Elution Buffer or Molecular Biology Grade Water.

**Storage of the eluate with purified DNA:** The eluate contains pure genomic DNA. For short- term storage (24-48 hrs) of the DNA, 2-8°C is recommended. For long-term storage, -20°C or lower temperature (-80°C) is recommended. Avoid repeated freezing and thawing of the sample which may cause denaturing of DNA. The Elution Buffer will help to stabilize the DNA at these temperatures.

#### **References**

1. Sambrook, J., *et al.* Molecular Cloning: A laboratory Manual, 2<sup>nd</sup> ed. (Cold Spring Harbor Laboratory Press, Plainview, NY, 1989)
2. Birren, B. and Lai, E. Pulsed Field Gel Electrophoresis: A practical guide (Academic Press, San Diego, CA, 1993)

#### **Warning and Precautions**

Not for Medicinal Use. Read the procedure carefully before beginning the protocol. Wear protective gloves/protective clothing/eye protection/face protection. Follow good clinical laboratory practices while handling clinical samples. Standard precautions should be followed as per established guidelines. Safety guidelines may be referred in safety data sheets of the product.

#### **Limitations**

The yield of DNA depends upon the type and the volume of starting material used.

## Performance and Evaluation

Each lot of HiMedia's HiPurA® Forensic Multi Sample DNA Purification Kit is tested against predetermined specifications to ensure consistent product quality.

## Quality Control

Type of Sample	DNA Purity
Body Fluids/ Cells/ Tissue	1.6-1.9

## Troubleshooting guide

Sr. No.	Problem	Possible Cause	Solution
1.	Spin column is clogged	Lysate/ethanol mixture is not homogenous	Vortex the tubes for atleast 5-10 seconds in order to obtain homogenous solution before applying it to the column. If minimally sheared DNA is required for downstream applications like PCR, mix with gentle pipetting or inversion until homogenous instead of vortexing.
		DNA elution is improper	Ensure that the DNA elution is in 50 µl of Elution Buffer. To improve the DNA yield incubate for 5 minutes at room temperature after Elution Buffer is added to the column.
		Ethanol was omitted during binding	Ensure that ethanol is added before adding the sample to the HiElute Miniprep Spin Column (Capped)
		Eluate contains residual ethanol from the wash	Remove ethanol from the second wash completely before eluting the DNA. Spin for additional 2 minutes to dry the membrane completely. In order to avoid the interference of ethanol, always use a fresh tube for elution.

		Use of water instead of Elution Buffer for elution of DNA	Elution Buffer is recommended for optimal yields and storage of the genomic DNA. If water is used instead of the Elution Buffer, the pH should be at least 7.0 to avoid acidic conditions which may cause acid hydrolysis of DNA when stored for long periods of time.  <b>NOTE:</b> Only DNase/RNase free water should be used for eluting DNA.
2.	Purity of the DNA is lower than expected; ( $A_{260}/A_{280}$ ratio is low)	Eluate was diluted in water for absorbance measurement.	Use either the Elution Buffer provided or 10 mM Tris-HCl, pH 8.0.
		Background reading is high due to silica fines	Spin the DNA sample at maximum speed for 1 minute, the supernatant can be used to repeat the absorbance readings.
3.	$A_{260}/A_{280}$ ratio is too high	RNA contamination	RNase A treatment can be included in future isolations or the final product can be treated with RNase A Solution and repurified.
4.	Shearing of genomic DNA	Improper handling of genomic DNA	All pipetting steps should be executed as gently as possible. Wide orifice pipette tips are recommended to eliminate shearing of the DNA to a large extent. If the isolated DNA is to be used for PCR, instead of vortexing mix with gentle pipetting or invert until homogenous. This reduces shearing of DNA considerably.
5.	Downstream applications are inhibited	Traces of ethanol present in the final genomic DNA preparation	After the washing steps, the eluate should not come in contact with the column. Spin the column for 2 minutes at maximum speed 12,000-16,000 x g ( $\approx$ 13,000-16,000 rpm) if necessary, after emptying the collection tube.
		Salt is carried over in the final eluate containing DNA	The column should be transferred to a new 2.0 ml collection tube before adding the wash solution in steps 10-11.

Please refer disclaimer Overleaf.

## Safety Information

The HiPurA® Forensic Multi Sample Genomic DNA Purification Kit is for laboratory use only, not for drug, household or other uses. The Lysis Solution (C1) contains chaotropic salts, which are irritants. Take appropriate laboratory safety measures and wear gloves when handling. Not compatible with disinfecting agents containing bleach. Please refer the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.









## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed off in accordance with current laboratory techniques.

## Technical Assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail to [mb@himedialabs.com](mailto:mb@himedialabs.com).

## Symbols

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIMB580  
Rev. No.: 09  
Date of Issue: 2025-06

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)