

MB585 HiPurA® SuperPlant DNA Maxiprep Purification Kit

Kit Contents

Product Code	Reagents provided	MB585	
		10 preps	25 preps
DS0200	SuperPlant Extraction Buffer	50 ml	115 ml
DS0070	Additive-II	5 ml	11.5 ml
DS0071	Additive-III	1 gm	2.5 mg
DS0003	RNase A Solution (20 mg/ml)	1.1 ml	2.6 ml
DS0019	Wash Solution Concentrate (WSP)	18 ml	39 ml
DS0040	Elution Buffer (ET) [10mM Tris-Cl, pH 8.5]	6 ml	14 ml
DBCC01	HiPure Maxiprep Column (in PW143 Collection Tube)	10 No.	25 No.
PW143	Collection Tubes (50 ml conical)	20 No.	50 No.

Introduction

DNA extraction from plant tissues, unlike DNA isolation from mammalian tissues, remains difficult due to the presence of a rigid cell wall surrounding the plant cells. This method can be used both on freeze dried leaves and on fresh leaves. The scale of extraction is dependent on the amount of starting material, for e.g. 1 gm of material requires 4.5 ml of SuperPlant Extraction Buffer and yields up to 150 µg of DNA.

HiPurA® SuperPlant DNA Maxiprep Purification Kit

This kit simplifies isolation of DNA from fresh plant material with spin column procedure. The procedure is optimized for a maximum of 1 gm of wet-weight of the starting material. The sample (fresh) is cut and ground in liquid nitrogen along with SuperPlant Extraction Buffer. SuperPlant Extraction Buffer contains CTAB (Cetyltrimethylammonium bromide), a detergent used to break open plant cells and solubilize the contents. Chlorophyll and some denatured proteins are removed from green plant tissue in an organic chloroform- isoamylalcohol step, and the organic phase is separated by centrifugation. Since the extract contains DNA and RNA, RNA can be removed by the addition of RNase A. The flow-through fraction is then mixed with a solution that enhances the binding of DNA to the column. The solution is then passed through HiPure Maxiprep Column that is followed by washing steps to remove trace contaminants. High quality DNA is eluted in the Elution Buffer (ET) provided in the kit.



Registered Office

HiMedia Laboratories Pvt Ltd.

Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area,
Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797
Tel : 00-91-22-6147 1919, 6903 4800

Fax : 6147 1920
Web : www.himedialabs.com
Email : info@himedialabs.com
mb@himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein

HiPure Maxiprep Column (DBCC01)

HiPure Maxiprep Column is based on the advanced silica binding principle presented in a microspin format. The system efficiently couples the reversible nucleic acid-binding properties of the advanced gel membrane and the speed plus versatility of spin column technology to yield high quantity of DNA. The use of spin column facilitates the binding, washing and elution steps thus enabling multiple samples to be processed simultaneously. This column eliminates the need for alcohol precipitation, expensive resins, and harmful organic compounds such as phenol and chloroform, otherwise employed in traditional DNA isolation techniques. DNA binds specifically to the advanced silica-gel membrane while contaminants pass through. PCR inhibitors such as divalent cations and proteins are completely removed in two efficient wash steps, leaving pure nucleic acid to be eluted in the buffer provided with the kit. The purified DNA is upto 20 - 30 kb in length and can be used for further downstream applications.

Elution

The yield of genomic DNA depends on the sample type and the number of cells in the sample. The eluted DNA ranges in size up to 20-30 kb, and is suitable for direct use in PCR, restriction digestion and Southern blotting applications.

Concentration, yield and purity of DNA

Spectrophotometric analysis and agarose gel electrophoresis will reveal the concentration and the purity of the genomic DNA. Use Elution Buffer (ET) to dilute samples and to calibrate the spectrophotometer, measure the absorbance at 260 nm, 280 nm, and 320 nm using a quartz microcuvette. Absorbance readings at 260 nm should fall between 0.1 and 1.0. The 320 nm absorbance is used to correct for background absorbance. An absorbance of 1.0 at 260 nm corresponds to approximately 50 µg/ml of DNA. The $A_{260}-A_{320}/A_{280}-A_{320}$ ratio should be 1.6-1.9. Purity is determined by calculating the ratio of absorbance at 260 nm to absorbance at 280 nm. DNA purified by HiPurA® SuperPlant DNA Purification Kit is free of protein and other contaminants that can inhibit PCR or other enzymatic reactions.

Concentration of DNA sample (µg/ml) = 50 x A_{260} x dilution factor.

Materials needed but not provided:

- Chloroform: Isoamylalcohol (24:1) (Product Code: MB115)
- Ethanol (96-100%)
- Mortar and pestle
- Liquid nitrogen
- Tabletop Centrifuge with rotor for 50 ml tubes capable of $\geq 3,500 \times g$ (5,000 rpm).
- 65°C water bath or heating block

Storage

Store the HiPurA® SuperPlant DNA Maxiprep Purification Kit between 15-25°C except certain components as specified on each labels. Under recommended condition kit is stable for 1 year.

RNase A enzyme treatment

RNase A is a type of RNase that is commonly used in research. RNase A (e.g., bovine pancreatic ribonuclease A) is one of the sturdiest enzymes in common laboratory usage. It cleaves 3'end of unpaired C and U residues.

Unit Definition for RNase A

One unit of the enzyme causes an increase in absorbance of 1.0 at 260 nm when yeast RNA is hydrolyzed at 37°C and pH 5.0. Fifty units are approximately equivalent to 1 Kunitz unit. It is completely free of DNases and proteases. The specific activity is 90 U/mg.

The product as supplied is stable at room temperature (15-25°C).

Centrifugation

All centrifugation steps are carried out in a conventional laboratory centrifuge e.g. Beckman CS-6KR, Heraeus Varifuge 3.0R, or Sigma 6k10 with fixed angle rotor. The tubes provided with the kit are compatible with almost all laboratory centrifuges and rotors. All centrifugation steps are performed at room temperature and are given in g, the correct rpm can be calculated using the formula:

$$RPM = \sqrt{RCF/1.118} \times 10^5 r$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary *g*-force.

General Preparation Instructions

1. Grinding of the plant material can be done using mortar and pestle. Midrib should be removed from the material before grinding, as midrib is a major source of carbohydrate contamination.
2. **SuperPlant Extraction Buffer:** Immediately prior to use, add 450 µl of Additive-II and 90 mg of Additive-III in 4.5 ml of SuperPlant Extraction Buffer. Preheat the solution to 65°C.
3. **Dilute Wash Solution Concentrate (WSP) (DS0019) as follows:**

Number of Preps	Wash Solution Concentrate (WSP)	Ethanol (96-100 %)
10	18 ml	42ml
25	39 ml	91ml

DNA Isolation Protocol

Sample Preparation

Finely cut the leaf material before grinding. Weigh 1 gm of the finely cut plant material and grind properly using a mortar and pestle in liquid nitrogen to a fine powder. Allow the liquid nitrogen to evaporate. DO NOT ALLOW THE SAMPLE TO THAW (keep samples on ice if needed). Proceed immediately to the DNA isolation protocol.

Protocol

NOTE: Ensure that Additive-II and Additive-III are added to SuperPlant Extraction Buffer (DS0200) as mentioned in General Preparation Instructions.

1. To 1 gm of the ground material add 4.5 ml of SuperPlant Extraction Buffer (DS0200) (preheated to 65°C) (Refer General Preparation Instructions) and transfer the sample to a Collection Tubes (50 ml conical). Mix by vortexing.
2. Incubate the samples for 60-90 minutes with occasional inversion at 65°C.
3. Add 4.5 ml of Chloroform: Isoamylalcohol (24:1) and mix gently by inversion for 5 minutes.
4. Centrifuge the samples at $\geq 3,500 \times g$ ($\geq 5,000$ rpm) for 15 minutes at room temperature (15-25°C).

5. Transfer the top aqueous layer (containing DNA) into a fresh collection tube (50 ml conical) and add 100 μ l of RNase A Solution (20 mg/ml) (DS0003). Incubate for 5-10 minutes at room temperature (15-25°C).
6. Add equal volume of Ethanol (96-100%) to the lysate obtained from the above step and mix by pipetting.
7. **Load lysate on HiPure Maxiprep Column (DBCC01)**
Transfer half (approx. 15-20 ml) of the lysate obtained from step 4 onto HiPure Maxiprep Column provided, placed in a 50 ml collection tube. Centrifuge at $\geq 3,500 \times g$ ($\geq 5,000$ rpm) for 3 minutes. Discard the flow-through liquid and place the column back into the tube.
8. Repeat the above step with the remaining sample. Discard the flow-through liquid and reuse the 50 ml collection tube.
9. **Wash**
(Prepare the diluted Wash Solution (WSP) (DS0019) as indicated in General Preparation Instructions)
Add 2.5 ml of diluted Wash Solution (WSP) and centrifuge at $\geq 3,500 \times g$ ($\geq 5,000$ rpm) for 2 minutes. Discard the flow-through liquid and place the column back into the tube.
10. Add another 2.5 ml of diluted Wash Solution (WSP) to the column and centrifuge for 2 minutes at $\geq 3,500 \times g$ ($\geq 5,000$ rpm). Discard the flow-through and reuse the same collection tube.
11. Centrifuge the empty column for 10 minutes at the same speed if residual ethanol is observed. Discard the collection tube containing the flow-through liquid and place the column in a new 50 ml collection tube.

NOTE: The column must be free of ethanol before eluting the DNA. Incubating the column for 10 minutes at 70°C to evaporate residual ethanol is recommended.

NOTE: If the centrifugal force is below 5,000 rpm, it is strongly recommended to follow the above step.

12. **DNA Elution**
Pipette 500 μ l of the Elution Buffer (ET) (DS0040) directly onto the column without spilling to the sides. Incubate for 5 minutes at room temperature. Centrifuge at $\geq 3,500 \times g$ ($\geq 5,000$ rpm) for 2 minutes to elute the DNA.

NOTE: For highly concentrated DNA, reload the eluate onto the column membrane, incubate at room temperature for 5 minutes. Centrifuge at $\geq 3,500 \times g$ ($\geq 5,000$ rpm) for 2 minutes.

Storage of the eluate with purified DNA: The eluate contains pure genomic DNA. For short-term storage (24-48 hours) of the DNA, 2-8°C is recommended. For long-term storage, -20°C or lower temperature (-80°C) is recommended. Avoid repeated freezing and thawing of the sample which may cause denaturing of DNA. The Elution Buffer (ET) will help to stabilize the DNA at these temperatures.

Precautions

Read the procedure carefully before starting the experiment.

Performance and Evaluation

Each lot of HiMedia's HiPurA® SuperPlant DNA Maxiprep Purification Kit is tested against predetermined specifications to ensure consistent product quality.

Quality Control

Type of Sample	DNA Yield	DNA Purity
Banana leaf	Up to 150 µg of DNA	1.6-1.9

Troubleshooting Guide

Sr. No.	Problem	Probable Cause	Solution
1.	Carbohydrate contamination in the sample	Grinding of the midrib along with the leaf material	Remove the midrib from the leaf before grinding as these plant species have prominent midribs. Removal of the midrib is not important in case of very young leaves.
2.	DNA appears degraded (as a smear running down the gel)	The plant material for freeze-drying is not immediately frozen	When harvesting plant material for freeze drying, ensure that the tissue is immediately frozen, as this reduces DNA degradation.
		DNA appears fragmented or broken	DNA being a large molecule can be broken by shear forces. Therefore, mix the samples gently. To minimize shearing, always use a wide bore pipette tip for mixing.
3.	Difficulty to dissolve DNA in Elution Buffer (ET)	This is due to over-drying of DNA pellet	The DNA should not be allowed to over-dry at any stage during the preparation, as it hinders the resuspension and solubilization in Elution Buffer (ET). Rehydrate the DNA by incubating at 65°C for 1 hour in Elution Buffer (ET).

Safety Information

HiPurA® SuperPlant DNA Maxiprep Purification Kit is for laboratory use only; not for drug, household or other uses. Take appropriate laboratory safety measures and wear gloves when handling. Avoid contact with skin, and use eye protection. In case of contact, wash with large amount of water. Seek medical attention. Not compatible with disinfecting agents containing bleach. Please refer the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.


Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed off in accordance with current laboratory techniques.

Please refer disclaimer Overleaf.

Technical assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail to mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



11/2025

PIMB585_0/1122

MB585-00

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com