

MBLA007

HiPurA® Cervical Specimen Collection Kit
(Liquid Cytology Collection Kit)

Intended Use:

Recommended for the collection of Endocervical samples

Introduction

HiPurA® Cervical Specimen Collection Kit is an essential tool in sample collection for PAP test. PAP test is a procedure in which a small brush is used to gently remove cells from the surface of the cervix and the area around it so they can be checked under a microscope for cell changes that may lead to cervical cancer. A Pap test may also help find other conditions, such as infections or inflammation.

Description

HiMedia's HiPurA® Cervical Specimen Collection Kit contains a Cervical brush which efficiently collects Endo cervical cells. This kit contains a unique proprietary buffer which keeps the cell morphology intact. It also helps to maintain the integrity of nucleic acid present inside the cells.

Kit Contents

S.No	Product Code	Description	Pack size (1NO)	Pack size (5NO)	Pack size (25NO)	Pack size (50NO)	Storage
1	MBLA007A	Sterile Disposable Endocervical Brush	1NO	5NO	25NO	50NO	15-25°C
2	MBLA007B	HiSolu™ Container	1 No	5NO	25NO	50NO	15-25°C

Salient features

- 100% safe sample collection
- Ready to use
- Easy and safe transport

Storage conditions

HiPurA® Cervical Specimen Collection Kit has to be stored at (15-25°C). Under recommended condition, the reagent is stable for 12 months.

Directions to Use

1. Insert the central (longest) bristles of the brush into the endocervical canal.
2. Apply gentle pressure until the shorter side bristles come in contact with the endocervix
3. Rotate the brush through 5 complete (360° clockwise revolutions with constant gentle pressure sufficient to keep the side bristles bent. **DO NOT OVER ROTATE.**

Please refer disclaimer Overleaf.

4. Remove the brush gently from the endocervical canal. Remove the brush head (only bristles) and quickly put in to the tube containing HiSolu™ Container. Cap the container tightly to avoid spillage. Discard the stem.
5. Label the container and store the specimen at Room Temperature or at 2-8°C.
6. It is recommended to send specimen to the laboratory soon after the collection.

Warning and Precautions

Read the SDS carefully before beginning the protocol. Wear protective gloves/protective clothing/eye protection/face protection. Follow good clinical laboratory practices while handling clinical samples. Standard precautions should be followed as per established guidelines. Safety guidelines may be referred in safety data sheets of the product.

Performance and Evaluation


Performance of the product is expected when it is stored at recommended temperature and within the expiry period.

Safety Information

The HiPurA® Cervical Specimen Collection Kit is for laboratory use only, not for drug, household or other uses. Take appropriate laboratory safety measures and wear gloves and safety goggles when handling. Not compatible with disinfecting agents containing bleach. Please refer the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



10/2027

PIMBLA007_O/1024

MBLA007-03

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com