

MBPCR105

Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit

Description:

Human papillomavirus (HPV) infection is caused by a human papillomavirus, a DNA virus from the papillomavirus family, of more than 200 related viruses and is the most common sexually transmitted infection (STI). While most HPV infections are benign causing warts (papillomas) on areas of the body including the hands, feet and genitals, certain strains put a person at a higher risk of developing particular types of cancers. HPV infection is limited to the basal cells of stratified epithelium, the only tissue in which they replicate. They typically spread by sustained direct skin-to-skin contact with oral, vaginal and anal sex being the most common methods. Risk factors for persistent HPV infections include early age at first sexual intercourse, multiple partners, smoking and poor immune function. Occasionally, it can spread from a mother to her baby during pregnancy. The high-risk HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 are known to cause several types of cancer. Among the high-risk HPV types identified, HPV types 16 and 18, are responsible for most HPV-caused cancers. Hence, the focus on these genotypes gives the physician useful details upon which to make effective treatment decisions. Real time PCR technique is considered simple and fast method for diagnosis with respect to the standard PCR technique. This technique has been successfully used for the rapid and accurate method for diagnosis of HPV infection in humans. Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit detects and differentiate HPV 16 and 18 genotypes in a single tube multiplex assay with high sensitivity and specificity.

NOTE: Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is for *in vitro* use only.

Intended Use:

Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is intended for use by qualified laboratory personnel trained in Real-Time PCR. The kit is recommended for sensitive and specific detection of HPV16 and HPV18 in clinical samples. The results from Hi-PCR Human Papilloma virus (HPV) Genotyping (Multiplex) Probe PCR kit must be interpreted within the context of all relevant clinical and laboratory findings.

Principle:

Real-time polymerase chain reaction, also called quantitative Polymerase Chain Reaction (qPCR) or kinetic Polymerase Chain Reaction, is a laboratory technique based on the principle of PCR. This technique is used to amplify a targeted DNA sequence by use of hydrolysis probes that are short oligonucleotides that have a fluorescent reporter dye attached to the 5' end and a quencher dye to the 3' end. Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is designed to detect **the E6 region of HPV16 and HPV18 in FAM channel and Texas Red channel respectively, and an endogenous Internal Control in Cy5 channel** in a single tube reaction. The kit allows sensitive detection and differentiation of HPV16 and HPV18 in a single tube assay.

Internal control (IC): An endogenous internal control is a DNA sequence that is naturally present in the test samples and is amplified in the same tube with a different set of primers and probe which may aid in identification of possible PCR inhibition, DNA purification efficiency and assessing the quality of the sample.

Registered Office

HiMedia Laboratories Pvt Ltd.

Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797
Tel : 00-91-22-6147 1919, 6903 4800


Fax : 6147 1920
Web : www.himedialabs.com
Email : info@himedialabs.com
mb@himedialabs.com



Positive control (PC): A positive control mimics a sample which contains all the target DNA sequences that the PCR is designed to amplify. It is included in a PCR assay to check the proper and intended functioning of all the PCR reagents.

No template control (NTC): A No Template Control is a PCR reaction that contains all the PCR reagents but does not include any DNA template. PCR grade water is used as the template to confirm that any observed amplification in the test samples is not due to contamination or non-specific amplification.

Diagrammatic representation of preferential binding of probe specific to DNA fragments in Real-Time PCR



While the probe is intact, the proximity of the quencher dye greatly reduces the fluorescence emitted by the reporter dye by fluorescence resonance energy transfer (FRET). The probes are designed such that they anneal within a DNA region amplified by a specific set of primers. During PCR amplification, these probes will hybridize to the target sequences located in the amplicon i.e. the DNA. As the *Taq* DNA polymerase replicates the template with the bound probe, the 5'-nuclease activity of the polymerase enzyme cleaves the fluorescent probe. The end result in cleavage of the probe is separation of the reporter dye from the quencher dye and increasing the reporter dye signal. As the probe is removed from the target strand, primer extension continues to the end of the template strand. Hence, fluorescence detected in the quantitative PCR thermal cycler is directly proportional to the fluorophore released and the amount of DNA template present in the PCR. Thus, inclusion of the probe does not inhibit the overall PCR process.

Features

Molecular Features:

- Simultaneous detection and differentiation of HPV16 and HPV18 in a single assay
- Highly sensitive – 2 IU/μL for both HPV16 and HPV 18
- Highly specific – No cross reactivity with pathogens with similar clinical presentation

Technology features:

- Fast and reliable results within 90 minutes.
- Includes all reagents & controls for validity of the test.
- Open system – Compatible with 4-channel and 5-channel qPCR cyclers.
- Wet-lab assays validated on the Bio-Rad CFX Opus 96, Applied Biosystems QuantStudio 5 and Insta Q96® Plus Real Time PCR Systems.

Recommended Sample Types: Cervical swabs, Cervico-vaginal samples in women, Formalin fixed paraffin embedded (FFPE) tissue specimen

Specimen collection and Handling

Follow appropriate techniques for handling specimens; after use, contaminated materials must be sterilized by autoclaving before discarding. Standard precautions as per established guidelines should be followed while handling clinical specimens and items contaminated with blood and other body fluids. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf life

The provided kit has a shelf-life of 12 months when stored between -10°C to -20°C. Repeated thawing and freezing of PCR reagents should be avoided, not more than 5 freeze-thaw cycles, as this may reduce the sensitivity. If the reagents are to be used multiple times, we recommend storing reagents as aliquots to avoid repeated freeze and thaw. Exposure to light, heat or humidity may also affect the shelf life of some of the kit components and should be avoided. Degradation of specimen/ extracted RNA can also reduce the sensitivity of the assay. HiMedia Laboratories does not recommend using the kit after the expiry date stated on pack.

Kit Contents: The provided PCR kit contains:

Components	Product Code	Reagents provided for (reactions)* (µL)		
		25R	50R	100R
HPV16-18 Master Mix	DS2285	135	270	530
HPV16-18 Primer Probe Mix	DS2286	81	162	318
HPV16-18 Positive Control	DS2301	54	108	212
Molecular Biology Grade Water for PCR	DS0440	162	324	648

*For a 20 µL PCR reaction

Materials needed but not provided

All materials are available through www.himedialabs.com

Product name	Product Code
Real-Time PCR Instrument and equipment	
Insta Q96® Real time PCR System, 96 well block, 5 channels	LA1012
Insta Q48® M4 Real time PCR System, 96 well block, 4 channels	LA1023
Insta Q96® Plus Real time PCR System, 96 well block, 5 channels	LA1073
Insta Q96® - 6.0 Real time PCR System, 96 well block, 6 channels	LA1074
TabSpin™ Microcentrifuge	LA1089/LA1090
Automated nucleic acid extraction system and materials	
Insta NX® Instrument - fully automated nucleic acid purification system utilizing the Innovative Super -S membrane column method	LA1056
Insta NX® Mag16, Insta NX® Mag16 ^{Plus}	LA1118, MBLA018
Insta NX® Mag32, Insta NX® Mag32 ^{Plus}	LA1096, MBLA019
Insta NX® Mag96	LA1097
Extraction Kits	
HiPurA® Viral DNA Purification Kit	MB575
HiPurA® Viral DNA/ RNA Purification Kit	MB582
HiPurA® Pre- filled Cartridges for Paraffin-Embedded Tissue Nucleic Acid Purification	MB530MPF16
HiPurA® Pre-filled Cartridges for Paraffin-Embedded Tissue Nucleic Acid Purification	MB530PC16
HiPurA® Pre- filled Clinical Multi- purpose Nucleic Acid Purification Kit (Cartridges)	MB583PC16200

HiPurA® Pre-filled Clinical Multi-purpose Magnetic Nucleic Acid Purification kit (Plates)	MB583MPF16200
HiPurA® Pre-filled Clinical Multi-purpose Magnetic Nucleic Acid Purification kit	MB583MPF32200
HiPurA® Pre-filled Clinical Multi-purpose Magnetic Nucleic Acid Purification kit	MB583MPF96200
Tubes, plates and other consumables	
Varivol II Micropipette-10 (Capacity: 0.5 to 10 µl)	LA611
Varivol II Micropipette-100 (Capacity: 10 to 100 µl)	LA615
Varivol II Micropipette-1000 (Capacity: 200 to 1000 µl)	LA614
Barrier Tips, Maximum capacity 10 µl	LA749/LA749A
Barrier Tips, Maximum capacity 200 µl	LA751/LA751A
Barrier Tips, Maximum capacity 1000 µl	LA859/LA859A
8-strip tubes & optically clear flat caps for PCR	PR17, PR22, PR23
PCR Tubes, 0.2 mL; PCR Plates	PW1255/ PR2/PR3/PR19
Optical Sealing film	PR18

Kit compatibility with Real-Time PCR Systems

The Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit contains fluorophores compatible with following PCR systems. However, its performance has been specifically validated on the Bio-Rad (CFX Opus 96), Applied Biosystems (QuantStudio 5), Insta Q96® AG and Insta Q96® Plus Real-Time PCR Systems. Below is the fluorophore compatibility matrix for commonly used instruments:

Real-Time PCR system	Company	Dye 1	Dye 2	Dye 3	Dye 4
Insta Q96®AG/ Insta Q96®AG 6.0/Insta Q96® - 6.0/Insta Q96® Plus/Insta Q48® M4	HiMedia Laboratories Pvt. Ltd.	FAM	JOE/HEX	Texas Red	Cy5
BioRad CFX Opus 96/CFX96 Touch/ CFX384 Touch	Bio-Rad Laboratories, Inc.	FAM	JOE/HEX	Texas Red	Cy5
QuantStudio™ 5	Applied Biosystems	FAM	JOE/HEX/VIC	Texas Red/ ROX	Cy5
ABI® Prism SDS 7500	Applied Biosystems	FAM	JOE/HEX/VIC	Texas Red/ ROX	Cy5
QIAquant 96 & 384 5plex	QIAGEN	FAM	JOE/HEX	Texas Red	Cy5
Rotor-Gene®6000 & Q	QIAGEN	Green	Yellow	Orange	Red
LightCycler® 96 /LightCycler® 480	Roche	FAM	JOE/HEX/VIC	ROX/ Texas Red	Cy5
qTOWER ³	Analytik Jena	FAM	JOE/HEX/VIC	ROX/ Texas Red	Cy5

Note: Ensure that your instrument is calibrated for the specific dyes listed above and is maintained according to the manufacturer's guidelines. Inconsistent calibration or dye incompatibility may result in abnormal amplification plots or failed detections.

Warning and Precautions

Not for Medicinal Use. Read the procedure carefully before beginning the protocol. Wear protective gloves/protective clothing/eye protection/face protection. Follow good clinical laboratory practices while handling clinical samples. Standard precautions should be followed as per established guidelines. Safety guidelines may be referred in safety data sheets of the product.

General Preparation Instructions

- Before use, all PCR components should be completely thawed on ice (4°C).
- Perform the amplification reactions in a clean area, preferably in a biosafety cabinet.
- Use of aerosol barrier pipette tips is recommended to reduce contamination risks from extraneous DNA templates.
- Extract and store positive control sample (if used) separately from all other reagents to avoid contamination and add it to the reaction mix in a separate area.

Protocol for PCR Master Mix Preparation (For one reaction)

1. In the “Master mix Preparation” area, thaw all components from the kit on ice, mix by inverting the tubes and centrifuge the reagents for several seconds. Keep on ice for later use.
2. Based on the number of specimens to be tested (N), including the PTC and NC, calculate the volume of the components to be added as N* volume of 1X.
3. Use 1.5 mL Nuclease free centrifuge tube(s) for the preparation of the reaction system. After all the reagents are added, mix them thoroughly and centrifuge for several seconds.

Components	Product code	Volume for “1X” (one reaction)
Preparation of PCR Reaction Mix		
HPV16-18 Master Mix	DS2285	5.0 µL
HPV16-18 Primer Probe Mix	DS2286	3.0 µL
Water (PCR grade)	DS0440	2.0 µL
Total PCR Reaction Mix	-	10.0 µL
Template addition		
Template/ Purified Viral DNA	-	10.0 µL
Total reaction volume	-	20.0 µL

4. Load 10µL of master mix into the 0.1/0.2 mL PCR reaction tube/plate/strips, compatible to the instrument to be used; add 10 µL of master mixture with 10 µL Nuclease free water to the negative control.
5. In the “Nucleic acid handling” area add 10 µL HPV Multiplex Positive Control and extracted test DNA into the plate/strip.
6. For positive control and no template control tube, Template Bacterial DNA is replaced by MTB/NTM Positive Control and by PCR grade water respectively. Refer the following table.

Set up of controls for the PCR run		
Components	Product code	Volume for “1X” (one reaction)
Total PCR Reaction Mix	-	10.0 µL
Positive Control		
HPV16-18 Positive Control	DS2301	10.0 µL
No Template Control		
Water (PCR grade)	DS0440	10.0 µL
Total reaction volume	-	20.0 µL

7. Tightly cap the tubes/strips or seal the plate using an optically clear adhesive film.
 8. Briefly, spin the strips/tubes to settle the reagent to the bottom of the tube.
- Place the plate/strips/tubes in Real-time PCR machine and set the PCR program. Interpret the data from the amplification plot (observe the Ct values).

Sr. No	Step	Temperature	Time	Sampling	No. of cycles
1.	Initial denaturation	95°C	10 minutes	---	1
2.	Denaturation	95°C	15 seconds	---	45
3.	Annealing	60°C	30 seconds	YES	
4.	Extension	72°C	15 seconds	--	

Selection of channels

Target	Channels	Quencher
HPV16	FAM	None
HPV18	Texas Red/ROX	None
Internal Control (IC)	Cy5	None

Please select 'Passive reference dye' as 'None' wherever applicable


Data Analysis

The following conditions should be met for a valid test:

Control	Target		
	HPV-16 (FAM)	HPV-18 (Texas Red)	Internal Control (Cy5)
Positive Control (PC)	+	+	+
No Template Control (NTC)	-	-	-

Target	Ct value	Result
HPV16	≤ 34	Detected (+)
HPV18	≤ 34	Detected (+)

Amplification Data



Sr. No	PC	Ct value	
		PC	NTC
1	HPV16	24	-
2	HPV18	23	-
3	IC	26	-

Note: Representative image showing amplification plot of HPV-16, HPV-18 and Internal Control with Ct values, using Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit on Applied Biosystems QuantStudio™ 5 System. The results completely depend upon sample types.

Data Interpretation:

Detection Channel			Result Interpretation
HPV-16 (FAM)	HPV-18 (Texas Red/ROX))	Internal Control (Cy5)	
+	-	+/-*	Positive for HPV16
-	+	+/-*	Positive for HPV18
+	+	+/-*	HPV16 and HPV18 Coinfection
-	-	+	Negative for HPV16 & HPV18
-	-	-	Inconclusive test**

*The presence or absence of a signal in the Cy5 channel is not relevant for the validity of the test run due to competition between the test template and Internal Control template.

** When an inconclusive result is obtained, repeat PCR or re-test the extracted DNA or re-extract the specimen and test the newly extracted DNA

Note:

- 1) Kindly correlate the results with clinical findings. A positive result should be interpreted in conjunction with clinical findings, such as cytology (Pap smear) or histopathology, to guide further diagnostic evaluation and patient management.
- 2) If sample is positive for two genotypes, re-test the specimen. If sample is repetitively positive for both the genotypes, the result may be indicative of a “*co-infection*”.
- 3) A negative result (no detection of HPV16 and HPV 18 genotypes) does not exclude HPV infection, as the assay is limited to these genotypes. The patient may still be infected with other common high-risk or low-risk HPV types, which this kit does not detect. Therefore, ***clinical findings and cytology must be considered***, and additional HPV testing may be recommended where broader genotyping is clinically indicated.
- 4) Visually inspect all amplification plots. Do not rely on Ct values alone. Please note that amplification curves exhibit a characteristic sigmoidal shape. Curves that are non-sigmoidal or/and show a sudden rise in fluorescence at very early Ct values (e.g., Ct <6) may be due to non-specific amplification, background noise, or instrument artifacts. Such results should be repeated or confirmed before making any clinical interpretation.

Performance Evaluation

Limit of Detection (LoD) - Analytical Sensitivity

The Limit of Detection (LoD) is defined as the concentration (copies per μl of the eluate) of target molecule that can be detected at 95% or greater probability according to CLSI EP17-A2. Sensitivity for the Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit was conducted using 20 replicates of dilution series of known copy number of WHO International Standard Collection of Human Papillomavirus (HPV) DNA Genotypes HPV16 (NIBSC code: 06/202) and HPV18 (NIBSC code: 06/206), each on Biorad CFX Opus 96, Applied Biosystems QuantStudio 5, Insta Q96[®] AG and Insta Q96[®] Plus Real Time PCR Systems. The detectable limit of the Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit on Real Time instrument was determined to be 2 IU/ μl for both HPV16 and HPV 18.

Inclusivity - Analytical Sensitivity

In silico analysis for the assessment of inclusivity for the Hi-PCR[®] Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit was conducted by mapping the primers and probes against all the

available Human Papilloma Virus (HPV) sequences in GenBank. The Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit targets 100% of the known Hi-PCR® Human Papilloma Virus (HPV) Genotyping 16 and 18.

Specificity

Wet testing analysis was performed against the following pathogens mentioned below in the table. No cross-reaction was observed with any strains.

Pathogen name	Name of Genomic standards used for the study	Catalog No.
<i>Neisseria gonorrhoeae</i>	Quantitative Genomic DNA from <i>Neisseria gonorrhoeae</i>	ATCC 700825DQ
<i>Trichomonas vaginalis</i>	Quantitative Genomic DNA from <i>Trichomonas vaginalis</i> strain C-1: NIH	ATCC 30001DQ
<i>Mycoplasma genitalium</i>	Quantitative Genomic DNA from <i>Mycoplasma genitalium</i> strain G37	ATCC33530DQ
<i>Chlamydia trachomatis</i>	Quantitative Genomic DNA from <i>Chlamydia trachomatis</i> Serovar D strain UW-3/Cx	ATCCVR-885DQ
<i>Treponema pallidum</i>	Quantitative Synthetic <i>Treponema pallidum</i> DNA	ATCCBAA-2642SD
<i>Gardnerella vaginalis</i>	Quantitative Genomic DNA from <i>Gardnerella vaginalis</i>	ATCC 14019DQ
Human adenovirus	Quantitative Genomic DNA from Human adenovirus 1 strain Adenoid 71	ATCCVR-1DQ
<i>Candida dubliniensis</i>	Quantitative Genomic DNA from <i>Candida dubliniensis</i> strain CBS 7987	ATCC 646DQ
<i>Neisseria meningitidis</i>	Quantitative Genomic DNA from <i>Neisseria meningitidis</i> strain FAM18	ATCC 700532DQ
<i>Candida albicans</i>	Quantitative genomic DNA from <i>Candida albicans</i> strain 3147	ATCC10231DQ
	Quantitative genomic DNA from <i>Candida albicans</i> strain SC5314	ATCC 2876DQ
Human papillomavirus 31	WHO International Standards for Human Papillomavirus (HPV) DNA genotypes HPV31, HPV33, HPV45, HPV52, HPV58	NIBSC 19/226
Human papillomavirus 45		
Human papillomavirus 33		
Human papillomavirus 52		
Human papillomavirus 58		
Human herpesvirus 3	Quantitative Genomic DNA from Human herpesvirus 3 (HHV-3) strain Ellen	ATCC 1367DQ
Human herpesvirus 5	Quantitative Genomic DNA from Human herpesvirus 5 (HHV-5) strain AD-169	ATCC 538DQ
Hepatitis B virus	Quantitative Synthetic DNA from Hepatitis B virus	ATCC 3232SD
Zika virus	Quantitative Synthetic RNA from Zika virus	ATCC 3252SD
Human immunodeficiency virus 1	Quantitative Synthetic Human immunodeficiency virus 1 (HIV-1) RNA	ATCC 3245SD
Human gammaherpesvirus 4	Quantitative Synthetic DNA from Human gammaherpesvirus 4	ATCC 3247SD
Human herpesvirus 1	Quantitative Genomic DNA from Human herpesvirus 1	ATCC 539DQ
Human herpesvirus 2	Quantitative Genomic DNA from Human herpesvirus 2	ATCC 540DQ

Hepatitis C virus	Quantitative Synthetic RNA from Hepatitis C virus	ATCC 3233SD
Ebolavirus	RNA from Bundibugyo Ebolavirus, Prototype Isolate #811250 (200706291 Uganda)	BEI NR-31812

Cross-Reactivity Analysis – *in silico*

The oligonucleotide sequences (primers and probes) used in the Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit were subjected to BLAST (Basic Local Alignment Search Tool) analysis against the organisms listed in the table below. No significant cross-reactivity was observed for any of the sequences during the in-silico BLAST analysis.

Pathogen	Taxonomy	Accession No.
<i>Haemophilus ducreyi</i>	[Haemophilus] ducreyi strain VAN2 chromosome, complete genome	NZ_CP015425.1
	[Haemophilus] ducreyi strain FDAARGOS_297 chromosome, complete genome	NZ_CP022037.2
Human Simplex virus 1	Human herpesvirus 1 strain 17, complete genome	NC_001806.2
Human Simplex virus 2	Human herpesvirus 2 strain HG52, complete genome	NC_001798.2
Human herpesvirus 6	Human betaherpesvirus 6A, variant A DNA, complete virion genome, isolate U1102	NC_001664.4
Group B Streptococcus (GBS)	Streptococcus sp. 'group B' strain FDAARGOS_229 chromosome, complete genome	NZ_CP020432.1
<i>Ureaplasma urealyticum</i>	Ureaplasma urealyticum strain 132 chromosome, complete genome	NZ_CP041200.1
	Ureaplasma urealyticum strain MIN-201 chromosome, complete genome	NZ_CP086123.1
<i>Ureaplasma parvum</i>	Ureaplasma parvum strain hebnu3h04 chromosome, complete genome	NZ_CP021987.1
	Ureaplasma parvum serovar 3 str. ATCC 700970, complete sequence	NC_002162.1
<i>Mycoplasma hominis</i>	Metamycoplasma hominis ATCC 23114, complete sequence	NC_013511.1
	Metamycoplasma hominis ATCC 27545 strain LBD-4 chromosome, complete genome	NZ_CP009652.1
HPV 6	Human papillomavirus type 6b, complete genome	NC_001355.1
HPV 11	Human papillomavirus type 11 (HPV-11) complete genome	NC_075235.1
HPV 40	Human papillomavirus type 7 genomic DNA	NC_001595.1
HPV 42	Human papillomavirus type 42 isolate TJ43-42, complete genome	GQ472847.1
HPV 43	Human papillomavirus type 43 complete genome	AJ620205.1
HPV 44	Human papillomavirus type 44, complete genome	NC_075252.1
HPV 54	Human papillomavirus 54, complete genome	NC_001676.1

HPV 61	Human papillomavirus type 61 isolate Qv3721, complete genome	KF436857.1
HPV 30	Human papillomavirus type 30 genomic DNA	NC_038889.1
Human papillomavirus	Human papillomavirus isolate SE379, complete genome	NC_027779.1
HPV 39	Human papillomavirus ORFs	NC_075238.1
HPV 35	Human papillomavirus type 35 isolate QV29782, complete genome	HQ537729.1
HPV 51	Human papillomavirus 51 isolate BF315, complete genome	KF436887.1
	Human papillomavirus 51 isolate 51YN25, complete genome	OR997941.1
HPV 56	Human papillomavirus type 56 genomic DNA	NC_075271.1
HPV 59	Human papilloma virus type 59, complete viral genome	NC_075272.1
HPV 66	Human papillomavirus type 66, complete genome	NC_075253.1
HPV 68	Human papillomavirus type 68a, complete genome	NC_075157.1
HPV 31	Human papillomavirus type 31 (HPV-31) complete genome	NC_075191.1
HPV 33	Human papillomavirus type 33, complete genome	NC_075233.1
HPV 45	Human papillomavirus type 45 genomic DNA	NC_075269.1
HPV 52	Human papillomavirus type 52 genomic DNA	NC_075270.1
HPV 58	Human papillomavirus type 58 complete genome	D90400.1
	Human papillomavirus 58 isolate ZWE054176, complete genome	KY225961.1

Limitations

- Strict compliance with the Instructions for Use is required for optimal results and the use of the kit is limited to staff qualified clinical laboratory personnel trained in the techniques of real-time PCR and in vitro diagnostic procedures.
- Appropriate specimen collection, transport, storage and processing procedures are required for the optimal performance of this test.
- This assay must not be used on the specimen directly. Viral DNA should be extracted from clinical sample using appropriate nucleic acid extraction method.
- Presence of PCR inhibitors and other interferences may lead to false negative or invalid results.
- Although rare, mutations within the highly conserved regions of the targets genes covered by the kit's primers and/or probe may result in failure to detect the presence of pathogen.
- As with any diagnostic test, results of the Hi-PCR® Human Papilloma Virus (HPV) HR Genotyping Probe PCR Kit need to be interpreted in consideration of all clinical and laboratory findings.
- Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit cannot detect other common high-risk HPV genotypes and all low-risk HPV genotypes.
- Performance of the kit in monitoring treatment of HPV infection has not been evaluated.

Evaluation

Each lot of Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is tested against predetermined specifications to ensure consistent product quality.

Quality Control

Each lot of Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is assayed for contaminating endonuclease, exonuclease and non-specific DNase activities. Functionally tested in amplification.

Troubleshooting Guide

Sr. No.	Problem	Possible Cause	Solution
1.	No amplification in test and/or control wells	Degraded samples or poor-quality DNA template	Use freshly extracted, high-quality DNA. Avoid repeated freeze-thaw cycles for stored DNA.
		Missing or incorrect addition of reagents	Verify all reagents were added in the correct volumes and order. Recheck reaction setup steps. Use a master mix to reduce pipetting errors.
		Incorrect thermal cycling conditions	Cross-check the PCR cycling profile with the IFU.
		Incorrect fluorophore/dye setting	Refer the selection of channels table and set the correct fluorophore in the real-time PCR instrument
		Expired or improperly stored reagents	Confirm the expiry date and proper storage. Avoid using reagents that have undergone multiple freeze-thaw cycles
		Instrument malfunction	Verify that the real-time PCR instrument is functioning properly and calibrated.
2.	Variability between replicates	Inconsistent pipetting or error in reaction setup	Prepare a single master mix for all replicates to minimize variation. Vortex thoroughly and aliquot carefully. Use calibrated pipettes and consistent technique.
		Air bubbles in reaction mix	Briefly centrifuge PCR tubes or plate before placing in the instrument to eliminate air bubbles.
		Uneven mixing of reagents	Ensure all reagents are fully thawed and mixed by vortexing before use. Spin down before pipetting.
		Edge effect (thermal variation across plate)	Avoid using outer wells in PCR plates if not temperature-uniform; use a plate seal and consistent plate layout.
3.	Amplification in negative control	Reagent contamination	Replace all critical reagents. Clean workspace, pipettes, and repeat analysis with fresh aliquots. Use filter tips and maintain a unidirectional workflow.
		Cross-contamination or aerosol contamination during reaction setup	Maintain strict unidirectional workflow. Set up reactions in a PCR hood or clean bench. Use aerosol-resistant filter tips for all pipetting steps. Avoid fast, forceful pipetting.

			Regularly clean work surfaces and equipment with DNA-decontaminating agents. Minimize opening of positive control tubes and avoid splashing.
		Template contamination in workspace or pipettes	Decontaminate work surfaces, pipettes, and equipment using DNA/RNA decontamination solutions. Perform regular cleaning.
		Improper sealing of PCR plate/tubes	Ensure plates/tubes are properly sealed to prevent cross-well contamination during thermal cycling. Use optical-grade seals if required.
4.	No signal with positive control or partial target amplification	Degradation of the positive control material due to improper storage or repeated freeze-thaw cycles	Use a fresh aliquot of positive control. Ensure storage conditions follow IFU and avoid repeated freeze-thaw. Discard expired or compromised controls.
		Incorrect thermal cycling conditions	Cross-verify cycler settings with IFU. Ensure annealing/extension temperature and time are as specified.
		Improper mixing of PC before use	Vortex and briefly spin down the positive control before adding it to the PCR mix. Mix well to ensure homogeneity.
5.	Early Ct value with non-sigmoidal amplification	Instrumental noise, background fluorescence, or non-specific amplification	Visually inspect amplification plots. Disregard flat or atypical (non-sigmoidal) curves that show a sudden rise in fluorescence with very early Ct values (e.g., <6). These may not indicate true amplification. Repeat the assay if necessary. Adjust threshold setting in software within exponential phase. Refer to instrument-specific guidance.
6.	Non-specific amplification of targets in samples and/or negative control	Improper threshold setting in the real-time PCR software	Visually inspect all amplification plots. Do not rely on Ct values alone. Manually adjust the threshold as recommended in the threshold setting table Refer to the instrument-specific guidelines or user manual for proper threshold adjustment.
		Contamination of reagents or workspace	Use fresh aliquots of reagents. Prepare reactions in a contamination-free environment using separate areas for pre- and post-PCR steps. Use aerosol-resistant filter tips and routinely decontaminate surfaces and pipettes.
7.	No Internal Control (IC) amplification	High HPV viral load may suppress IC amplification	If valid Ct for HPV target is observed, it is still considered a valid result. Test repetition with freshly extracted sample is recommended.
		Improper / insufficient sample	If no signals from any of the HPV targets and also from the IC was observed, this indicates improper or insufficient sample collection. Repeat the sample collection.

		Presence of PCR inhibitor	Dilute the extracted nucleic acid in molecular biology grade water and repeat the PCR assay.
		Sample inhibition, extraction failure, or expired reagents	Repeat extraction or test with fresh sample. Check IC amplification to rule out PCR inhibition. Use validated extraction method and fresh reagents.
		Instrument malfunction	Verify that the real-time PCR instrument is functioning properly and calibrated.
8.	High Ct values in positive samples (late amplification)	Low viral load or suboptimal sample	Results near the cut-off should be interpreted with caution. Repeat test using freshly extracted DNA. Confirm with orthogonal test if clinically significant.
9.	Signal in only one replicate (of duplicate or triplicate reactions)	Pipetting error or borderline positivity	Repeat the test. If consistent upon retesting, interpret cautiously in context of clinical findings. Borderline cases (close to Ct cutoff) may require repeat sampling or orthogonal testing (e.g., NGS). Use calibrated pipettes and proper technique.

Safety Information

Hi-PCR® Human Papilloma Virus (HPV) Genotyping (Multiplex) Probe PCR Kit is for laboratory use only, not for drug, household or other uses. Take appropriate laboratory safety measures and wear gloves when handling.










Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures while disposing the infectious materials. Material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Technical Assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com

Symbols

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number
	In vitro diagnostic medical device		

Identification No.: PIMBPCR105

Rev.No.:15

Date of Issue: 2025-09

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.