

MBPCR264

Hi-PCR[®] *Vibrio parahaemolyticus* Semi-Q PCR Kit

Description

Vibrio parahaemolyticus is a Gram-negative halophilic bacterium belonging to the genus *Vibrio*. It is ordinarily found in estuarine, marine and coastal environments, but can also be found in livestock and poultry meat, freshwater fish, preserved eggs, and pickles as a contaminant. It is a major food-borne pathogen causing stomachache, diarrhea, vomiting, dehydration, chills, and fever after the consumption of contaminated raw or undercooked seafoods including crab, shrimp, shellfish, lobster, fish, and oysters. Recently, *V. parahaemolyticus* was reported as a major causative agent of acute hepatopancreatic necrosis syndrome troubling penaeid shrimp, thereby inflicting major losses to the shrimp aquaculture industry. An early detection of *V. parahaemolyticus* contamination is therefore a prerequisite. To meet enhanced diagnostic accuracy within a shorter time as compared to traditional methods, Semi-Q PCR method is developed for detection of *V. parahaemolyticus* in food samples.

NOTE: The Hi-PCR[®] *Vibrio parahaemolyticus* Semi-Q PCR Kit is for research use only.

Intended Use

Recommended for sensitive and specific detection of *Vibrio parahaemolyticus* in food samples.

Principle

Polymerase chain reaction is used to amplify a targeted DNA sequence by use of specific primers. HiMedia's Hi-PCR[®] *Vibrio parahaemolyticus* Semi-Q PCR Kit is a qualitative conventional PCR kit which allows amplification of a *V. parahaemolyticus* specific target gene along with the internal control in a single tube reaction. The amplified target of 341 bp product size is confirmed by using agarose gel electrophoresis. This kit also contains an **Internal control** and a **Positive control**.

Positive control

This is a control reaction using a known template. A positive control is usually used to check that the primers have been designed properly and the PCR conditions have been set up correctly.

Internal control

This is a control sequence which is amplified in the same reaction tube along with the target sequence (target pathogen) but detected with a different primer (i.e., Multiplex PCR). An internal control is often used to detect the failure of amplification in cases where the target sequence is not amplified.

Features

- Fast and simple
- Good sensitivity and specific results
- Guaranteed reproducible results

Sample Source: Food sample

Storage and Shelf life

The provided kit has a shelf-life of 12 months when stored between -15°C to -25°C. Repeated thawing and freezing of PCR reagents should be avoided, as this may reduce the sensitivity. If the reagents are to be used multiple times, we recommend storing reagents as aliquots to avoid repeated freeze and thaw. Degradation of sample DNA specimens can also reduce the sensitivity of the assay. HiMedia Laboratories does not recommend using the kit after the expiry date stated on pack.

Kit Contents: The provided PCR contains reagents for 25 µL/reaction volume:

Components	Product Code	Reagents provided for (reactions)*	
		25R	50R
2X PCR TaqMixture	MBT061	338	675
V. parahaemolyticus primer mix	DS1661	68	135
V. parahaemolyticus Positive Control	DS1663	25	50
Water	DS0440	135	293
6X Gel Loading Buffer	ML015	54	108
100 bp DNA Ladder	MBT049	20	40

Specimen collection and Handling

Follow appropriate techniques for handling contaminated food specimens; after use, contaminated materials must be sterilized by autoclaving before discarding. Standard precautions as per established guidelines should be followed while handling contaminated food specimens. Safety guidelines may be referred to in individual safety data sheets.

Sample Preparation

For preparation of bacterial DNA, perform nucleic acid purification using HiPurA® Multi-Sample DNA Purification Kit (MB554) as described in the protocol.

Materials needed but not provided

- PCR tubes (Product code PW1255) or PCR Strips (Product code: PR17) or PCR Plates (Product code: PR2/ PR3/ PR19) & Sealing film (PR18)
- Thermal Cycler (Product Code: LA948/ LA949/ LA950/ LA1006/ LA1015/ LA1059/ LA1060/ LA1066)
- Barrier Micropipette Tips (Product Code: LA749/ LA749A/ LA751/ LA751A/ LA750/ LA750A/ LA859/ LA859A)
- Micropipettes

General Preparation Instructions

- Before use all PCR components should be completely thawed on ice (4°C).
- Perform the amplification reactions in a clean area, preferably in a biosafety cabinet.
- Use of aerosol barrier pipette tips is recommended to reduce contamination risks from extraneous DNA templates.
- Extract and store positive control sample (if used) separately from all other reagents to avoid contamination and add it to the reaction mix in a separate area.
- Centrifuge the components briefly once thawed.

A. Protocol for PCR Master Mix Preparation

Perform PCR reactions for each DNA sample as per the following table:

Components	Product Code	Recommended volume to be added per reaction (µL)
2X PCR TaqMixture	MBT061	12.5
V. parahaemolyticus primer mix	DS1661	2.5
Template DNA / Positive Control / Negative	-	5
Water	DS0440	5
Total volume		25

NOTE: (Optional) – The user can also set up an additional PCR reaction containing 1µL of Positive control DNA (provided) in a separate tube.

Centrifuge the tube briefly at 6000 rpm for about 10 seconds and place the tubes in the PCR machine and set the recommended PCR program (mentioned below). Interpret the data using Agarose Gel Electrophoresis.

B. Recommended PCR program


- | | | |
|-------------------------|-----------------------|---------------------|
| 1. Initial denaturation | : 95°C for 5 minutes | } No. of cycles: 35 |
| 2. Denaturation | : 95°C for 30 seconds | |
| 3. Annealing | : 63°C for 20 seconds | |
| 4. Extension | : 72°C for 30 seconds | |
| 5. Final extension | : 72°C for 10 minutes | |

After amplification, the products can be kept at 4°C overnight or frozen at -20°C for long-term storage.

C. Visualizing PCR Assay Results

- For analysis of the PCR results, load 5 µL of amplicon on a 1.5% agarose gel along with 1 µL of 6X Gel loading buffer (ML015). Load 3 µL of 50 bp DNA ladder (MBT084) in a separate well.
- Incorporate EtBr in the agarose gel or stain the agarose gel after electrophoresis using EtBr solution for 10-15 mins. Confirm the expected amplicon size by comparing with 50 bp DNA marker.

Data Interpretation:



LANES	CONTENT
1	100 bp Ladder
2-3	Negative control
4-5	Positive control

Warning and Precautions

Read the procedure carefully before beginning the protocol. Wear protective gloves/protective clothing/eye protection/face protection. Follow good clinical laboratory practices while handling contaminated food samples. Standard precautions should be followed as per established guidelines. Safety guidelines may be referred in safety data sheets of the product.

Performance and Evaluation

Each lot of HiMedia's Hi-PCR® *Vibrio parahaemolyticus* Semi-Q PCR Kit is tested against predetermined specifications to ensure consistent product quality.

Quality Control

Each lot of HiMedia's Hi-PCR® *Vibrio parahaemolyticus* Semi-Q PCR Kit is assayed for contaminating endonuclease, exonuclease and non-specific DNase activities. It has been functionally tested in DNA amplification assays.

Troubleshooting Guide

Sr. No.	Problem	Cause	Solution
1.	No amplification	Degraded samples	1. Check the integrity of DNA using agarose gel electrophoresis. 2. Use freshly prepared DNA to ensure the availability of intact template sequence for efficient amplification.
		Error in protocol setup	Verify that the correct reagent volumes, dilutions and storage conditions have been used.
2.	Variability between replicates	Error in reaction set-up	Prepare a large volume master mix, vortex thoroughly and aliquot into reaction tubes.
		Air bubbles in reaction mix	Briefly centrifuge reaction samples/plate prior to running on a real-time PCR instrument.
		Pipetting error	Replicates can show increased variation due to poor laboratory technique or imprecise pipettes.
3.	Amplification in negative control	Reagents contaminated	1. Replace all critical solutions. 2. Repeat the analysis of all tests with fresh aliquots of critical reagents.
4.	No signal with positive controls	Incorrect programming of the temperature profile of the thermocycler	Compare the temperature profile to the manual.

Safety Information

HiMedia's Hi-PCR® *Vibrio parahaemolyticus* Semi-Q PCR Kit is for laboratory use only, not for drug, household or other uses. Take appropriate laboratory safety measures and wear gloves when handling.









Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with contaminated food must be decontaminated and disposed of in accordance with current laboratory techniques.

Technical Assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIMBPCR264

Rev.No.:02

Date of Issue: 2025-07

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA. Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: mb@himedialabs.com Website: www.himedialabs.com