



Tetrathionate Brilliant Green Bile Broth

O3477

Intended Use:

Recommended for isolation and identification of *Salmonellae*.

Composition**

Component	Quantity
Peptone	8.600
Bile	8.000
Sodium chloride	6.400
Calcium carbonate	20.000
Potassium tetrathionate	20.000
Brilliant green	0.070
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Preparation

Suspend 63.07 grams in 1000 ml purified / distilled water. Heat just to boiling. DO NOT AUTOCLAVE OR REHEAT. Dispense into sterile tubes or flasks as desired.

Note: Due to the presence of calcium carbonate, the prepared medium forms opalescent solution with white precipitate.

Applications

Salmonella are gram-negative, facultatively anaerobic, non-sporulating, non-motile rods in the family *Enterobacteriaceae*. They are widely distributed in animals affecting mainly the stomach and the intestines. These organisms are difficult to differentiate biochemically from *Escherichia coli*. Tetrathionate Broth was originally described by Mueller (1) and later modified by Kauffman (2,3). Tetrathionate Brilliant Green Bile Broth is used as an enrichment medium for *Salmonella*.

Enrichment broth is usually recommended to facilitate the recovery of small numbers of *Salmonella* species (4). Tetrathionate Brilliant Green Bile Broth is also mentioned in I.P. (5) for isolation and identification of *Salmonella* species from foods, water and other materials of sanitary importance.

Peptone in the medium provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and nutrients for growth of *Salmonellae*. Brilliant green and bile inhibit both gram-positive as well as some selected gram-negative organisms. Potassium tetrathionate inhibits normal flora of faecal specimens. Sodium chloride helps in maintaining osmotic equilibrium.

After incubation, streak the culture from Tetrathionate Brilliant Green Bile Broth (M1255) onto differential medium for isolation and identification. Tetrathionate Brilliant Green Bile Broth is not suitable for growth of *Salmonella* Typhi and *Salmonella* Paratyphi (6).

Uses

Food samples

Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (1,11).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. For further confirmation, streak the enriched cultures after incubation, on plates of Brilliant Green Agar (M016), MacConkey Agar (M081) and Bismuth Sulphite Agar (M027).

Swenk{"Eqvtq

Crgctpeg

Light yellow to pale green homogeneous free flowing powder

Eqqwt"cpf"Enctkv{"qh"rtgrctgf"ogfkw o

Bluish green coloured opalescent solution with white precipitate.

Tgcevkqp

Reaction of 6.3% w/v aqueous solution at 25°C. pH : 7.0±0.2

r J

6.80-7.20

Ewnwtcn" Tgurqpug

Cultural characteristics observed when subcultured on MacConkey Agar (M082) after an incubation at 35-37°C for 18-24 hours.

Qti cpluo	Kpqewnw o *EHW+	I tqyv j	Tgeqxtg{	Eqqwt"qh eqqp{
<i>Escherichia coli</i> ATCC 25922	50-100	fair	20-30%	pink to red with bile precipitate
<i>Salmonella Typhi</i> ATCC 6539	50-100	luxuriant	>=50%	colourless
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	>=50%	colourless
<i>Salmonella Enteritidis</i> ATCC 13076	50-100	luxuriant	>=50%	colourless
<i>Staphylococcus aureus</i> ATCC 25923	>=10 ⁴	inhibited	0%	
<i>Staphylococcus aureus</i> ATCC 6538	>=10 ⁴	inhibited	0%	
<i>Escherichia coli</i> ATCC 8739	50-100	fair	20-30%	pink to red with bile precipitate
<i>Escherichia coli</i> NCTC 9002	50-100	fair	20-30%	Pink to red with bile precipitate
<i>Staphylococcus aureus</i> NCIMB 9518	>=10 ⁴	inhibited	0%	

Tghgtpeg

- Mueller L., 1923, C. R. Soc. Biol., (Paris), 89, 434.
- Kauffman F., 1930, Hyg. Abt. I. Orig., 113, 148.
- Kauffman F., 1935, Z. Hyg. Infektionskr., 117, 26.
- Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Eds.). 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- Indian Pharmacopoeia, 1996, Ministry of Health and Family Welfare, Govt. of India,
- MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

Revision : 02/ 2016

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

