



GBS Medium Base

M1073

GBS Medium Base is used for the isolation and rapid detection of Group B Streptococci (GBS) in clinical specimens

Composition**

Ingredients	Gms / Litre
Proteose peptone	23.000
Sodium dihydrogen phosphate	1.500
Disodium hydrogen phosphate	5.750
Starch, soluble	80.000
Final pH (at 25°C)	7.5±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 55.12 grams in 475 ml distilled water. Dissolve completely by gently heating to boiling for 15-20 minutes. Sterilize by autoclaving at 15 lbs (121°C) for 15 minutes. Cool to 60°C and aseptically add 25 ml sterile inactivated Horse serum (RM1239) and sterile rehydrated contents of 1 vial of GBS Supplement (FD054). Mix well and dispense into sterile test tubes. For the formation of gel, tubes are to be refrigerated (2-8°C) over-night before use.

Principle And Interpretation

Beta-haemolytic Streptococci with Lancefield group B antigen (*Streptococcus agalactiae*) are an important cause of serious neonatal infection characterized by sepsis and meningitis. Heavy colonization of the maternal genital tract is associated with colonization of infants and risk of neonatal disease (1). GBS Medium, formulated by Islam (2-4) is recommended for the isolation and detection of group B Streptococci (GBS) from clinical specimens. GBS Medium is designed to exploit the ability of most Group B Streptococci (GBS) to produce orange /red pigmented colonies when incubated under anaerobic conditions. The orange red pigment of group B Streptococci also has the characteristic of a carotenoid (3). GBS Medium Base also supports growth of other genital bacteria that cause perinatal infections (1), e.g. anaerobic *Streptococcus*, *Bacteroides* and *Clostridium* species.

Proteose peptone provides the necessary nutrients for the growth of Group B Streptococci. The phosphate salts buffer the medium. The antibiotic supplement (FD054) makes the medium selective for Group B Streptococci, while the horse serum enriches the media. Colonies of Group B Streptococci are 0.5 to 1 mm in diameter, round, entire and give pigmented growth (orange/red) after 24-48 hours anaerobic incubation. Other organisms that can grow on this medium do not produce the orange/red pigment.

Vaginal or rectal swabs should be inserted vertically into the medium. Incubation is carried out at 35-37°C. Pigment production is observed at hourly interval. Colour change (due to pigment production) of the butt occurs gradually, starting from the bottom of the tube towards the upper end. Presence of blood in the specimen may give false positive results. Presumptively positive tubes should be further confirmed by biochemical analysis to identify Group B Streptococci.

Quality Control

Appearance

Cream to beige homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent solution without any precipitate in tubes

Reaction

Reaction of 11.02% w/v aqueous solution at 25°C. pH : 7.5±0.2

pH

7.30-7.70

Cultural Response

Cultural characteristics observed with added inactivated Horse serum (RM1239) and GBS Supplement (FD054) after an incubation at 35-37°C for 24-48 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth	Pigmentation
Cultural Response			
<i>Bacteroides fragilis</i> ATCC 25285	50-100	fair to good	no pigmentation
<i>Streptococcus agalactiae</i> ATCC 13813	50-100	good-luxuriant	orange/red
<i>Enterococcus faecalis</i> ATCC 29212	50-100	good-luxuriant	no pigmentation

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

Reference

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Islam A. K. M. S., 1977, Lancet i : 256-7 (letter).
3. Merrit K. and Jacobs N. J. 1978, J. Clin. Microbiol. 8, 105-7.
4. Atlas R. M. 2004, Handbook of Microbiology Media, 3rd Edition, CRC Press, 704-705

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.