



HiCrome™ MRSA Agar Base, Modified

M1953

Intended Use:

Recommended for the differentiation and identification of MRSA and MRSE *Staphylococcus* species from clinical samples.

Composition**

Ingredients	Gms / Litre
Peptone	23.000
Sodium chloride	10.000
Sodium pyruvate	5.000
Chromogenic substrate	0.770
Inhibitor mixture	7.000
Agar	15.000
Final pH (at 25°C)	7.2±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 30.38 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Aseptically add sterile rehydrated contents of 1 vial of MeReSa Selective Supplement (FD229) or Cefoxitin supplement (FD259) or both in combination for more selectivity as desired. Mix well and pour into sterile Petri plates.

Principle And Interpretation

MRSA is a resistant variation of the common bacterium *Staphylococcus aureus* and MRSE is a resistant variation of the common bacterium *Staphylococcus epidermidis*. *Staphylococcus aureus* is an invasive pathogen that can cause disease in almost any tissue or organ in the human body, primarily in compromised individuals (1). Staphylococcal infections were earlier treated using Penicillin. But over the years resistance to this drug developed. Methicillin was the next drug of choice. While methicillin is very effective in treating most *Staphylococcus* infections, some strains have developed resistance to methicillin and can no longer be killed by this antibiotic. These resistant bacteria are called Methicillin Resistant *Staphylococcus aureus* (MRSA) (5).

Patients with breaks in their skin due to wounds, indwelling catheters or burns are those with certain risk of developing MRSA infection (5). Spread of MRSA infections can be controlled to a great extent by maintaining personal hygiene after interaction with an MRSA infected person (5).

Peptone provide the essential nutrients along with carbonaceous, nitrogenous and Vitamin B complex nutrients. The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* to give green coloured colonies. Sodium pyruvate enhances the growth of *Staphylococcus* species. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting the accompanying microflora. Cefoxitin is recommended to use for selective isolation of MRSA. The medium is made selective for MRSA by the addition of MeReSa Selective Supplement (FD229) or Cefoxitin supplement (FD259) or both in combination.

Type of specimen

Clinical samples - skin, blood, body fluids

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (3,4). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Some intermediate strains may show poor growth due to nutritional variations and resistance to methicillin/cefoxitin.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation must be carried out by sensitivity testing.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to beige homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light purple, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.08% w/v aqueous solution 25°C. pH : 7.2±0.2

pH

7.00-7.40

Cultural Response

Cultural characteristics observed with added MeReSa Selective Supplement (FD229) or Cefoxitin Supplement (FD259) or both after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Growth w/ FD229 or FD259 or both	Recovery w/ FD229 or FD259 or both	Colour of Colony
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	poor-good	30-40%	inhibited	0%	purple
<i>Staphylococcus aureus</i> , MRSA ATCC 43300	50-100	luxuriant	≥50%	luxuriant	≥50%	green
<i>Staphylococcus epidermidis</i> , MRSE	50-100	luxuriant	≥50%	luxuriant	≥50%	blue
<i>Staphylococcus xylosus</i> ATCC 29971	50-100	luxuriant	≥50%	inhibited	0%	blue
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	none-poor	≤10%	inhibited	0%	

Key : (*) Corresponding WDCM numbers

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

Reference

1. DWorkin M et. al 2006. The Prokaryotes (a Handbook on the Biology of Bacteria) 3rd ed, Vol. 2, page 345.
2. Dr. Alan Johnson, methicillin resistant staphylococcus aureus (MRSA) infection. The Support group for MSRA sufferers and Dependents, Aug 1st, 2005.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
5. Methicillin Resistant Staphylococcus aureus Copyright a 1997-2005 Canadian Centre for Occupational Health and Safety, Sept 19th, 2005.


Revision : 02 / 2020



In vitro diagnostic medical device



CE Marking



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Pvt. Limited,
23 Vadhani Industrial Estate,
LBS Marg, Mumbai-86, MS, India



CE Partner 4U ,Esdoornlaan 13, 3951
DB Maarn The Netherlands,
www.cepartner4u.eu

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.