



## OF Basal Medium

M

### Intended use

OF (Oxidation Fermentation) Basal medium is used for the determination of oxidative and fermentative metabolism of carbohydrates by gram-negative bacteria.

### Composition

| Ingredients           | Gms / Litre |
|-----------------------|-------------|
| Tryptone              | 2.000       |
| Sodium chloride       | 5.000       |
| Dipotassium phosphate | 0.300       |
| Bromo thymol blue     | 0.080       |
| Agar                  | 2.000       |
| Final pH ( at 25°C)   | 6.8±0.2     |

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 9.38 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense in 100 ml amounts and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. To first 100 ml of sterile basal medium, aseptically add 10 ml of sterile 10% dextrose solution. To second 100 ml add 10 ml sterile 10% lactose solution. To third 100 ml add 10 ml sterile 10% saccharose solution. Mix and dispense aseptically in 5 ml amounts in sterile tubes in duplicate for aerobic and anaerobic fermentation.

### Principle And Interpretation

Hugh and Leifson developed OF Medium to study oxidative and fermentative metabolism of carbohydrates by gram-negative bacteria. This criterion is used during taxonomic studies of *Enterobacteriaceae* (2).

Tryptone in the medium provides the necessary carbon and nitrogen, vitamins etc required for bacterial growth. A carbohydrate whose fermentation reaction is to be studied is added separately. Phosphate buffers the medium and the low agar concentration determines motility and dispersion of the acid produced on the surface ( ). Bromothymol blue acts as the pH indicator. The low concentration of agar permits the determination of motility and aids in the even distribution of any acid produced at the surface of the medium ( ). Motility is observed as diffused zone of flaring out from the line of inoculation. Non-motile organisms grow along the line of inoculation.

Dextrose is the most important carbohydrate for use in OF Basal Medium. However, certain organisms may metabolize other carbohydrates even if they are unable to utilize dextrose. Degradation of the carbohydrate to acid is indicated by the pH indicator bromothymol blue which changes its colour to yellow. Oxidative utilization takes place when the medium is exposed to air while fermentative utilization occurs under exclusion of air.

The authors Hugh and Leifson showed that when a gram-negative organism is inoculated in this medium containing a carbohydrate in duplicate, of which one tube is covered with mineral oil to exclude oxygen and the second tube is uncovered; reactions of differential value can be observed. Fermentative organisms will produce an acid reaction in both the covered and uncovered medium. Oxidative organisms will produce an acid reaction in the uncovered medium and give slight growth without change in the covered medium. Organisms which are not classified either as oxidative or fermentative show no change in the covered medium and an alkaline reaction in the uncovered medium ( ).

Prepare the medium with 1% dextrose and without 1% dextrose. Two tubes of each carbohydrate are used per organism and inoculated by stabbing. One of the inoculated tubes of each carbohydrate medium is covered with 2 ml of sterile mineral oil and the other is left uncovered. The tubes are incubated at 35-37°C for 18-48 hours or longer. The results are read after 48 hours.

The acidic reaction of oxidative organisms is more apparent at the surface of the medium that gradually spreads throughout the medium. If the oxidation reaction is weak or slow, an initial alkaline reaction at the surface of the uncovered tube may persist for several days and eventually convert to an acid reaction.

OF Basal Medium can be supplemented with 2% serum or yeast extract (0.1%) to make the medium more nutritious for the growth of bacteria ( ).

## Type of specimen

Pure isolate from clinical and non clinical samples.

## Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4, ).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (1).

After use, contaminated materials must be sterilized by autoclaving before discarding.

## Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations

1. Since the medium is pH dependent, pH must be monitored.
2. Certain fastidious organisms may not grow and require supplementation of 2% serum or yeast extract (0.1%) (4).

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

Cream to greenish yellow homogeneous free flowing powder

### Gelling

Semisolid, comparable with 0.2% Agar gel.

### Colour and Clarity of Prepared medium

Green coloured clear to slightly opalescent gel forms in tubes.

### Reaction

Reaction of 0.94% w/v aqueous solution at 25°C. pH : 6.8±0.2

### pH

6.60-7.00

### Cultural Response

M395: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

| Organism                                  | Inoculum (CFU) | Only Basal Medium (aerobic)                   | Only Basal Medium (overlaid with mineral oil) | w/ Dextrose (aerobic)                                       | w/Dextrose (overlaid with mineral oil)                      |
|---|----------------|---|---|---|---|
| <b>Cultural Response</b>                  |                |   |   |   |   |
| <i>Acinetobacter baumannii</i> ATCC 1     | 50-100         | alkaline reaction, green colour of the medium | alkaline reaction, green colour of the medium | acidic reaction, yellowing of the medium                    | alkaline reaction, green colour of the medium               |
| <i>Alcaligenes faecalis</i> ATCC 875      | 50-100         | alkaline reaction, green colour of the medium | alkaline reaction, green colour of the medium | alkaline reaction, green colour of the medium               | alkaline reaction, green colour of the medium               |
| <i>Escherichia coli</i> ATCC 2 22 (0001 ) | 50-100         | alkaline reaction, green colour of the medium | alkaline reaction, green colour of the medium | acidic reaction, yellowing of the medium with gas formation | acidic reaction, yellowing of the medium with gas formation |

|  |        |  |  |   |   |
|--|--------|--|--|---|---|
| # <i>Klebsiella aerogenes</i><br>ATCC 1 48           | 50-100 | alkaline<br>reaction, green<br>colour of the<br>medium | alkaline<br>reaction, green<br>colour of the<br>medium | acidic reaction,<br>yellowing of<br>the medium<br>with gas<br>formation | acidic reaction,<br>yellowing of<br>the medium<br>with gas<br>formation |
| <i>Pseudomonas aeruginosa</i><br>ATCC 9027 (00026*)  | 50-100 | alkaline<br>reaction, green<br>colour of the<br>medium | alkaline<br>reaction, green<br>colour of the<br>medium | acidic reaction,<br>yellowing of<br>the medium                          | alkaline<br>reaction, green<br>colour of the<br>medium                  |
| <i>Salmonella Enteritidis</i><br>ATCC 13076 (00030*) | 50-100 | alkaline<br>reaction, green<br>colour of the<br>medium | alkaline<br>reaction, green<br>colour of the<br>medium | acidic reaction,<br>yellowing of<br>the medium<br>with gas<br>formation | acidic reaction,<br>yellowing of<br>the medium<br>with gas<br>formation |
| <i>Shigella flexneri</i> ATCC<br>12022 (0012 )       | 50-100 | alkaline<br>reaction, green<br>colour of the<br>medium | alkaline<br>reaction, green<br>colour of the<br>medium | acidic reaction,<br>yellowing of<br>the medium                          | acidic reaction,<br>yellowing of<br>the medium                          |
| <i>Vibrio cholerae</i> ATCC<br>15748                 | 50-100 | alkaline<br>reaction, green<br>colour of the<br>medium | alkaline<br>reaction, green<br>colour of the<br>medium | acidic reaction,<br>yellowing of<br>the medium                          | acidic reaction,<br>yellowing of<br>the medium                          |

*Corresponding C numbers*

# Formerly known as *Enterobacter aerogenes*

## Storage and Shelf Life

Store between 10 – 30 °C in a tightly closed container and the prepared medium at 2 – 8 °C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry, ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

## Disposal

4,5

## Reference

1. , , , 1 , 14 ,
2. Hugh R. and Leifson E., 1953, J. Bacteriol. 66:24.
3. Cowan, 1974, Cowan and Steeles Manual for the Identification of Medical Bacteria, 2nd Ed., Cambridge University Press, Cambridge, Mass.
4. , , , , 11 , , 1 , , , , , , , , 15
5. MacFaddin J. F., 2000, Biochemical Tests for Identification of Medical Bacteria, 3rd Ed., Lippincott , Williams & Wilkins, Baltimore, Md .
6. MacFaddin J. F., 1985, Media for the Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore , 5 , , 15, , ,


Revision : 02 / 201



In vitro diagnostic medical device



CE Marking



Storage temperature



Do not use if package is damaged



CE Partner 4U, Esdoornlaan 13, 3951 DB Maarn The Netherlands, [www.cepartner4u.eu](http://www.cepartner4u.eu)



CE Partner 4U, Esdoornlaan 13, 3951 DB Maarn The Netherlands, [www.cepartner4u.eu](http://www.cepartner4u.eu)

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.