



Technical Data

C.L.E.D HiCynth™ Agar w/ Andrade Indicator

MCD352

C.L.E.D. HiCynth™ Agar w/ Andrade Indicator is recommended for isolation and differentiation of urinary pathogens on the basis of lactose fermentation.

Composition**

Ingredients	Gms / Litre
HiCynth™ Peptone No.3*	11.000
Lactose	10.000
L-Cystine	0.128
Bromothymol blue	0.020
Andrade indicator	0.100
Agar	15.000
Final pH (at 25°C)	7.5±0.2

**Formula adjusted, standardized to suit performance parameters

*- Chemically defined peptone

Directions

Suspend 36.25 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Sandys reported a new technique where the swarming of *Proteus* on an agar medium could be prevented by restricting the electrolyte content in the culture medium (1). Sandys Medium was modified by Mackey and Sandys (2), by replacing mannitol with lactose and sucrose and elevating the concentration of agar and bromothymol blue. The same authors further modified this medium by retaining the lactose (deleting sucrose) and by including L-cystine for promoting the growth of cystine-dependent dwarf coliform colony (3). This later modified medium was designated as C.L.E.D. (Cystine- Lactose- Electrolyte-Deficient) Medium. This medium is recommended for use in urinary bacteriology, promoting the growth of all urinary pathogens. C.L.E.D. Medium is also recommended for dip stick procedures and as dip inoculum transport medium for urine specimens (2, 3, 4).

C.L.E.D. Medium was further modified by Bevis (5) by incorporation of Andrades indicator. C.L.E.D. HiCynth™ Agar w/ Andrade Indicator is the modification of C.L.E.D. Agar w/ Andrade Indicator using chemically defined peptone free from animal and vegetable peptones, to avoid BSE/TSE risks associated with animal peptones. This medium provides sharper differentiation between lactose-fermenters (LF) and lactose-non-fermenters (NLF) (5). Addition of Andrades indicator enhances the appearance of colony and aids in the identification of microorganisms.

At different pH values, the colour of the medium varies from the standard medium, which is well documented by Bevis (5).

pH	Colour of C.L.E.D. medium
7.4	deep blue
7.0	bluish grey
6.8	pale grey
6.6	pinkish grey
6.4	bright red with whitish tinge
6.0	bright red

For better results, the medium should not be incubated for more than 24 hours because if lactose fermenters predominate, the entire medium may turn pink masking the presence of non-lactose fermenters. Inoculate the medium immediately after urine collection. *Shigella* species may not grow on this medium. Prior initiation of antibiotic therapy, low urine pH (less than 5) etc. may result in low urine count from infected patients.

Please refer disclaimer Overleaf.

Quality Control

Appearance

Light yellow to greyish yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Greenish blue clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.62% w/v aqueous solution at 25°C. pH : 7.5±0.2

pH

7.30-7.70

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours

Cultural Response

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
Cultural Response				
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	good-luxuriant	≥70%	greyish green, mucoid
<i>Enterococcus faecalis</i> ATCC 29212	50-100	good-luxuriant	≥70%	orange-yellow or greenish
<i>Escherichia coli</i> ATCC 25922	50-100	good-luxuriant	≥70%	bright pink with pink halo
<i>Proteus mirabilis</i> ATCC 25933	50-100	good-luxuriant	≥70%	blue-green
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant	≥70%	pale pink to golden-yellow
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	good-luxuriant	≥70%	greyish green

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Sandys, 1960, J. Med. Lab. Technol., 17:224.
2. Mackey and Sandys, 1965, Br. Med. J., 2:1286.
3. MacKey and Sandys, 1966, Br. Med. J., 1:1173.
4. Dixson J. M. S. and Clark M. A., 1968, Conc. Med. Assoc. J., 99 (15)
5. Bevis T. D., 1968, J. Med. Lab. Technol., 25:38.

Revision : 01 / 2015



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.