

WORLD CLASS QUALITY

HiCrome™

Single Streak Rapid Differentiation Series

Clinical Testing



COLOURS that

Identify the Pathogens

HIMEDIA®
For Life is Precious

HiMediaLaboratories™
himedialabs.com

Clinical Testing

Urinary Tract Infection

■ HiCrome™ UTI Agar, Modified / HiCrome™ UTI Selective Agar


M1418 / M1505

A chromogenic differential medium for identification, differentiation and confirmation of enteric bacteria from specimens such as urine which may contain large number of *Proteus* species as well as potentially pathogenic gram-positive organisms.

- Chromogenic mixture to detect presence of β -glucosidase and β -D-galactosidase enzymes.
- One chromogenic substrate is cleaved by β -glucosidase possessed by Enterococci resulting in formation of blue colonies.
- *Escherichia coli* produce purple-magenta colonies due to β -D-galactosidase which cleaves the other chromogenic substrate.
- Rich source of phenylalanine and tryptophan provides an indication of tryptophan deaminase activity by *Proteus* species, *Morganella* species and *Providencia* species.

<i>Escherichia coli</i>	purple-magenta
<i>Enterococcus faecalis</i>	blue-blue green (small)
<i>Klebsiella pneumonia</i>	blue to purple, mucoid
<i>Proteus mirabilis</i>	light brown
<i>Pseudomonas aeruginosa</i>	colourless (greenish pigment may be observed)
<i>Staphylococcus aureus</i> *	golden yellow

* inhibited on selective media (M1505)




■ HiCrome™ UTI Agar

M1353 / M1353R

Recommended for presumptive identification and confirmation of microorganisms mainly causing urinary tract infections, can also be used for testing water, food, environmental and other clinical samples.

- It has broader application as a general nutrient agar for isolation of various microorganisms.
- It facilitates the identification on basis of different contrasted colony colours produced by enzymes with two chromogenic substrate.
- Enterococci possess β -glucosidase giving blue coloured colonies
- *E.coli* produce purple-magenta colonies due to β -D-galactosidase which cleaves the other chromogenic substrate.
- M1353R opaque background helps in better visibility

<i>Escherichia coli</i>	pink-purple
<i>Enterococcus faecalis</i>	blue-blue green (small)
<i>Klebsiella pneumonia</i>	blue to purple, mucoid
<i>Proteus mirabilis</i>	light brown
<i>Pseudomonas aeruginosa</i>	colourless (greenish pigment may be observed)
<i>Staphylococcus aureus</i>	golden yellow




ESBL/Carbapenem Resistant *Enterobacteriaceae*

HiCrome™ ESBL Agar Base

M1829

Recommended for selective isolation of Extended Spectrum β -lactamase producing *Enterobacteriaceae*.

- Medium supports rapid growth with detection in 18-24 hours.
- Selective supplement helps in selective growth of ESBL producing *Enterobacteriaceae*
- Distinct colours allows proper differentiation
Escherichia coli - pink to purple
Resistant KESC group - bluish green
Proteus, *Morganella* and *Providencia* - colourless to light brown




HiCrome™ KPC Agar Base

M1831

Recommended for detection of gram negative bacteria with a reduced susceptibility to carbapenem agents.

- Medium supports rapid growth with detection in 18-24 hours.
- Selective supplement helps in selective growth of Carbapenem resistant strains.
- Distinct colours allows proper differentiation
Escherichia coli - pink to purple (confirmation with indole)
Citrobacter freundii - pink to magenta
Carbapenem resistant *Klebsiella*, *Enterobacter*, *Serratia* spp. - bluish green
Acinetobacter and *Salmonella*- smooth colourless colonies
Pseudomonas species - colourless to light yellowish green translucent colonies with wrinkled edges.



Clinical Testing

Yeast and moulds

HiCrome™ Candida Differential Agar/ Base Modified

M1297A/M1456A

Recommended for rapid isolation and identification of *Candida* species from mixed cultures.

- Medium supports growth with detection in 40-48 hours.
- Chromogenic mixture helps differentiate between *Candida* species
- Chloramphenicol and Selective supplement (in M1456A) helps in inhibiting bacterial growth.


Candida albicans - green

Candida tropicalis - blue-purple

Candida krusei - purple, fuzzy

Candida glabrata - cream to white

Candida parapsilosis - cream to white (may have mauve centre)



HiCrome™ Candida Differential Agar Base

M1297AR

Recommended for rapid isolation and identification of *Candida* species from mixed cultures.

- Medium supports growth with detection in 40-48 hours.
- Chromogenic mixture helps differentiate between *Candida* species
- Chromogenic mixture contains X-NAG which detects hexosaminidase activity and BCIP which detects phosphatase activity.
- Opaque background for better visibility.
- Selective supplement helps in inhibiting bacterial growth.

Candida albicans - light green smooth


Candida tropicalis - blue - metallic blue raised

Candida krusei - pink -purple, fuzzy

Candida glabrata - cream to white

Candida kefyr - cream to white

Candida parapsilosis - cream to white (may have mauve centre)



HiCrome™ Malassezia Agar (Twin Pack)

M1985

Recommended for isolation, cultivation and identification of *Malassezia furfur*.

- Medium supports growth with detection in 48-72 hours.
- Tween 80, Glycerol monooleate and fatty acids helps support luxuriant growth of *Malassezia furfur*.


Malassezia furfur - mauve, small

Candida albicans - light green

Candida tropicalis - metallic blue raised

Candida krusei - pink -purple, fuzzy

Candida glabrata - cream to white



Clinical Testing


Pseudomonas

HiFluoro™ *Pseudomonas* Agar Base

M1469

Recommended for selective isolation of *Pseudomonas aeruginosa* from non-clinical and clinical samples by fluorogenic method.

- Rapid detection in 18-24 hours
- Fluorogenic compound is specifically cleaved by *Pseudomonas* to give fluorescence under uv
- Cefrimide - inhibits accompanying microflora other than *Pseudomonas*
- Salts impart pigmentation




Klebsiella

HiCrome™ *Klebsiella* Selective Agar Base

M1573

Recommended for detection and isolation of *Klebsiella* from water and clinical samples.

- Selective medium for the growth of *Klebsiella*
- Chromogenic substrate is cleaved by *Klebsiella* to give purple-magenta coloured mucoid colonies
- Sodium lauryl sulphate and Bile salts mixture - inhibits gram positive organisms.
- Selective supplement (carbenicillin) inhibits other accompanying microflora.




Acinetobacter

HiCrome™ *Acinetobacter* Agar Base

M1938

Recommended for detection of *Acinetobacter* species from clinical specimen.

- Medium supports growth with detection in 24-48 hours.
- Chromogenic mixture imparts light purple colour to *Acinetobacter* (MDR)
- Selective supplement helps in selective growth.



Clinical Testing


Group B Streptococci

■ HiCrome™ Strep B Selective Agar Base

M1840

Recommended for the selective isolation of Group B *Streptococcus* from clinical samples.

- Medium to support growth in 18-24 hours
- Chromogenic substrate is cleaved by β -glucosidase enzyme that imparts blue colour to Group B *Streptococcus*
- Selective supplement - inhibits accompanying microflora
- Opaque background for better visibility



■ HiCrome™ Strep B Selective Agar Base, Modified

M1966

Recommended for the selective isolation of Group B *Streptococcus* from clinical samples.

- Chromogenic substrate is cleaved by group B-Streptococci resulting in purple coloured colonies.
- Other Streptococci give bluish green to green coloured colonies with yellow background due to fermentation indicated by phenol red.
- Selective supplements - inhibits other microorganisms.



Clinical Testing

Staphylococcus aureus (MRSA/MRSE)

HiCrome™ MeReSa Agar Base

M1674

Recommended for the isolation and selective identification of Methicillin Resistant *Staphylococcus aureus* (MRSA) from clinical isolates.

- The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* (MRSA) to give bluish green coloured colonies.
- *Staphylococcus aureus* (MSSA) and other accompanying flora - inhibited by selective supplement (Methicillin and Cefoxitin)




HiCrome™ MRSA Agar Base, Modified

M1953

Recommended for the differentiation and identification of MRSA and MRSE *Staphylococcus* species from clinical isolates.

- The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus*
- *Staphylococcus aureus* (MRSA) - green
- *Staphylococcus epidermidis* (MRSE)-blue
- *Staphylococcus aureus* (MSSA) and other accompanying flora - inhibited by Inhibitor mixture and selective supplement (Methicillin and Cefoxitin)



HiCrome™ Rapid MRSA Agar Base

M1974

It is recommended for rapid isolation and identification of Methicillin Resistant *Staphylococcus aureus* (MRSA) from clinical isolates

- Rapid detection in 18-24 hours
- The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus*. Carbohydrate fermentation is detected by phenol red indicator
- Methicillin Resistant *Staphylococcus aureus* (MRSA) - greenish yellow (Note: Green colour may develop after 48 hours)
- Methicillin Resistant *Staphylococcus epidermidis* (MRSE) - blue
- *Staphylococcus aureus* (MSSA) and other accompanying flora - inhibited by selective supplement.



Clinical Testing


Vancomycin Resistant Enterococci

HiCrome™ VRE Agar Base

M1830

Recommended for the identification and differentiation of VRE from clinical samples.

- Medium to support growth in 24-48 hours
- Chromogenic substrate is cleaved by β -glucosidase enzyme which imparts bluish green colour to VRE species
- Selective supplement - inhibits VSE and other gram positive organisms
- Opaque background for better visibility




HiCrome™ VRE Agar Base, Modified

M1925

Recommended for the identification and differentiation of Vancomycin Resistant *Enterococcus faecalis* and *Enterococcus faecium* from clinical samples.

- Medium to support growth in 24-48 hours
- Chromogenic substrate is cleaved by β -glucosidase, enzyme which imparts blue green colour to *Enterococcus* species
- Presence of arabinose and phenol red aids to differentiate between *Enterococcus faecalis* and *Enterococcus faecium*
- *E. faecalis*(VRE) - blue colour and *E. faecium* (VRE) - green w/yellow background
- Selective supplement - inhibits VSE and other gram positive organisms



Clinical Testing


HiCrome™ Mueller Hinton Agar

M2010

Recommended for differentiation of organisms from clinical samples based on chromogenic identification and determination of susceptibility of microorganisms to antimicrobial agents.

- Traditional method takes 48 hours for organism identification and antimicrobial susceptibility
- This medium gives rapid and reliable results in 24 hours
- Chromogenic differentiation of various Urinary Tract pathogens
- Simultaneous detection of Antimicrobial susceptibility
- Can be employed in clinical testing of urinary tract infection

<i>Escherichia coli</i>	purple-magenta
<i>Enterococcus faecalis</i>	blue-blue green
<i>Klebsiella or Enterobacter</i>	metallic blue
<i>Proteus mirabilis</i>	light brown
<i>Pseudomonas aeruginosa</i>	colourless (greenish pigment may be observed)
<i>Staphylococcus aureus</i>	golden yellow



HiCrome™ Mueller Hinton Agar (For Antifungal Testing)

M2067

Recommended for the chromogenic differentiation of yeast from clinical samples and determination of susceptibility to antifungal agents.

- Identification of organisms and antifungal susceptibility takes 96 hrs with traditional method
- The medium gives rapid and reliable results in 48 hours.
- Chromogenic differentiation of yeast cells along with antifungal susceptibility
- Can be employed in clinical testing

<i>Candida albicans</i>	Green
<i>Candida tropicalis</i>	Blue
<i>Candida krusei</i>	Purple
<i>Candida parapsilosis</i>	Cream
<i>Saccharomyces cerevisiae</i>	Cream



Asia

Africa

Australia

HiMedia Laboratories Pvt. Ltd.
www.himedialabs.com



HiMedia is a global brand with network reach spanning over 150 countries across the world.

North America

South America

Europe

OVERSEAS OFFICE

- North America -

HiMedia Laboratories LLC,
507 School House Rd,
Kennet Square 19348,
PA, USA.

Tel : +1-484-734-4401

Fax : +1-484-734-4402

Email : infous@himedialabs.com

- CORPORATE OFFICE -

Plot No. C-40,
Road No. 21Y, MIDC, Wagle Industrial Estate,
Thane (West) - 400604,
Maharashtra, INDIA.

Tel : +91-22-6147 1919 / 6116 9797 / 6903 4800

Fax : +91-22-6147 1920

Email : info@himedialabs.com

OVERSEAS OFFICE

- Europe -

HiMedia Laboratories GmbH,
Marie-Curie-Str. 3, 64683,
Einhausen, Germany.

Tel : +49 6251 989 24 26

Fax : +49 6251 989 24 27

Email : infoeu@himedialabs.com



For iOS



For Android



For Website

HIMEDIA®

For Life is Precious