



VF Selective Supplement

FD277

Intended use:

Recommended for selective isolation of Vancomycin Resistant Enterococci (VRE).

Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Vancomycin	4mg
Fluconazole	5mg

Directions:

1. This product is available in one size containing 5 vials.
2. The bottles when supplied are intact. Ensure that bottles do not have any cracks or defects.
3. User may remove the desired number of bottles from the box as per their requirement.
4. It should be handled by trained person wearing appropriate personal protective equipment (PPE) and sterile gloves.
5. Place the bottles on sterile surfaces such as laminar air flow or sterile working bench.
6. Label them accordingly.
7. Disinfect the outer surface of cap or closures with suitable disinfectant example 70% IPA.
8. Rehydrate the contents of 1 vial aseptically with 5 ml of 50% v/v aqueous ethanol.
9. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) HiCrome™ VRE Agar Base M1830/ HiCrome™ VRE Agar Base, Modified M1925. Mix well and pour into sterile Petri plates.
10. Follow good lab practices for procedures and disposal.

Type of specimen

Clinical samples - faecal sample

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (1,2). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning & Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1. Isenberg (Ed.),2004, Clinical Microbiology Procedures Handbook, Vol.3, American Society for Microbiology, Washington. D.C.
2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology,11th Edition. Vol. 1.

* Not For Medicinal Use

Revision : 00/2024

Packaging

FD277-5VL - VF Selective Supplement



HiMedia Laboratories Pvt. Limited,
Plot No.C-40, Road No.21Y,
MIDC,WagleIndustrial Area,
Thane (W) -400604, MS, India




CEpartner4U,Esdoornlaan13,
3951DB Maarn, NL
www.cepartner4u.eu



**In vitro diagnostic
medical device**



CE Marking



Storage temperature



**Do not use if
package is damaged**



**Single sterile
barrier system**

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.