

## Mucormycosis Detection Kit

K144

### Intended Use:

Recommended for sample collection, transportation, detection of suspected *Mucor* species (after sporulation appear as black fungi) and *Candida* species (appears white in colour) causing Mucormycosis.

### Kit contents :

	<b>2 kits</b>	<b>10 kits</b>
	<b>Sufficient for</b>	<b>Sufficient for</b>
	<b>2 samples</b>	<b>10 samples</b>
1. MP5476 : Mucormycosis Selective Agar Plate	2 plates	10 plates
2. MP5477 : Candida Selective Agar Plate	2 plates	10 plates
3. MS5478 : HiFungal™ Transport medium w/swab (Sterile flocked nylon swab for specimen collection).	2 nos.	10 nos.

Refer individual technical data sheets for more details.

### Advantages :

1. A complete one stop solution for the collection and transport of sample and selective isolation of *Mucor* species (black fungus and *Candida* species (white fungus).
2. MS5478 : Collection and transport of sample of fungi.
3. MP5476 : Selective isolation of *Mucor* fungi inhibiting bacterial contamination and yeasts.
4. MS5477 : Selective isolation of *Candida* species inhibiting bacterial contamination and mycelial fungi.
5. Lockable plates to ensure clinician safety.

### Directions

1. Cut open the pouch containing the swab for specimen collection
2. Collect the specimen using standard techniques.
3. After sample collection, insert the swab into the HiFungal™ Transport medium and break the swab at the breakpoint.
4. Tighten the cap firmly.
5. The specimen will be preserved during transportation and also the viability of the organisms will be maintained. After the transportation, the specimen should be inoculated in respective medium as soon as possible.
6. MP5476- Mucormycosis Selective Agar Plate - selective cultivation of *Mucor* species (after sporulation appear as black fungi).
7. MP5477- Candida Selective Agar Plate -selective isolation of *Candida* species (appears as white fungi)
8. For cultivation of other fungi any general fungal medium as Sabouraud Dextrose Agar (GM063/MP063) or Potato Dextrose Agar (M096/ MP096) etc.

### Principle And Interpretation

Mucormycosis (previously called zygomycosis) is a serious but rare fungal infection caused by a group of molds called mucormycetes (1). In most cases it is due to an invasion of the genera *Rhizopus* and *Mucor*, common bread molds (2). Effective recovery of microorganisms and its identification is dependent on a number of factors such as collection and transportation to the laboratory under conditions which allow maintenance of viability. This medium gives improved recovery of fungal cultures including yeasts and moulds. The recovery of mycelial fungi including *Aspergillus*, *Mucor*, *Rhizopus* is also well supported. The fungal cultures can be safely transported. Antibacterial components have been added to the transport medium to avoid bacterial contamination.

Mucormycosis selective Agar Plate is specifically designed to selectively promote the growth of *Mucor* species and inhibit the contaminating bacterial flora and yeast.

Candida Selective Agar Plate is specifically designed to selectively promote the growth of *Candida* species and inhibit the contaminating bacterial flora and mycelial fungi.

## Type of specimen

Clinical samples : eye lesion, nasal swabs, other sites of infection

## Specimen Collection and Handling:

HiFungal™ Transport medium w/swab (MS5478) which contains the transport medium along with the swab is used for specimen collection from the eye, nose, nasopharynx, other sites of infection and transported to the lab. can be used. (Refer directions).

After use, contaminated materials must be sterilized by autoclaving before discarding.

## Warning and Precautions :

In Vitro Diagnostic Use only. Read the label before opening the kit. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations :

1. Microscopy, morphological tests, staining tests and further biochemical tests are recommended for confirmation.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

MP5476 - Mucormycosis Selective Agar Plate  
 MP5477 - Candida Selective Agar Plate  
 MS5478 - HiFungal™ Transport Medium with sterile flocced nylon swab.

### Sterility test

Passes release criteria

### Cultural response

Viability of following organisms was established for a period of 48 hours. For selective cultivation of Mucormycosis fungi, the swab should be swabbed onto Mucormycosis Selective Agar Plate (MP5476) and for selective isolation of white fungi, *Candida* species, the swab must be swabbed onto Candida Selective Agar Plate (MP5477). For cultivation of fungi any general fungal medium such as Sabouraud Dextrose Agar, Potato Dextrose Agar etc can be used. The plates are then incubated at 25-30°C for 24-48 hours. Overincubation should be avoided.

### Cultural Response

Organism	Growth on Mucormycosis Selective Agar Plate (MP5476)	Growth on Candida Selective Agar Plate (MP5477)
<i>Mucor racemosus isolate</i>	Luxuriant (Black spores observed)	none-poor
<i>Rhizopus oryzae</i> MTCC 1987	Luxuriant (Black spores observed)	none-poor
# <i>Aspergillus brasiliensis</i> ATCC 16404 (00053*)	none-poor	none-poor
<i>Candida albicans</i> ATCC 10231 (00054*)	none-poor	Luxuriant
<i>Trichophyton rubrum</i> ATCC 28191	none-poor	none-poor
<i>Escherichia coli</i> ATCC 25922 (00013*)	inhibited	inhibited
<i>Staphylococcus aureus subsp. aureus</i> ATCC 25923 (00034*)	inhibited	inhibited
<i>Pseudomonas aeruginosa</i> ATCC 27853 (00087*)	inhibited	inhibited

Key : (#) - Formerly known as *Aspergillus niger*, (\*) - corresponding WDCM numbers

## Storage and Shelf Life

Store between 20-30°C. Use before expiry date on the label. Product performance is best if used within stated expiry period.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

## Reference

1. "Mucormycosis". NORD (National Organization for Rare Disorders). Archived from the original on May 26, 2021. Retrieved May 25, 2021
2. Lee, Soo Chan; Idmurm, Alexander (2018). "8. Fungal sex: The Mucoromycota". In Heitman, Joseph; Howlett, Barbara J.; Crous, Pedro W.; Stukenbrock, Eva H.; James, Timothy Yong; Gow, Neil A. R. (eds.). The Fungal Kingdom. Wiley. pp. 177–192. ISBN 978-1-55581-958-3.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.


Revision : 00/ 2021



In vitro diagnostic medical device



CE Marking



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Pvt. Limited,  
23 Vadhani Industrial Estate,  
LBS Marg, Mumbai-86, MS, India



CE Partner 4U ,Esdoornlaan 13, 3951  
DB Maarn The Netherlands,  
[www.cepartner4u.eu](http://www.cepartner4u.eu)

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.