

## HiCombi™ Dual Performance Salmonella Medium - SS

LQ029A

### Intended use

Recommended as a qualitative test for rapid growth and confirmation of *Salmonella*. Combination of solid (20 ml) and liquid (40 ml) media in single bottle.

### Composition\*\*

Ingredients	g/ L
Solid	7.000 ml
Peptone	5.000
HM peptone B #	5.000
Lactose	10.000
Bile salts mixture	8.500
Sodium citrate	10.000
Sodium thiosulphate	8.500
Ferric citrate	1.000
Brilliant green	0.00033
Neutral red	0.025
Agar	15.000
Liquid	20.000 ml

Same as solid media without Agar

\*\*Formula adjusted, standardized to suit performance parameters

# Equivalent to Beef extract

### Directions

Label the ready to use LQ029A bottle. Remove the top seal of the cap. Disinfect the part of the rubber stopper which is now exposed. Transfer the sample immediately into the culture bottle by puncturing the rubber stopper with the needle. Venting: Use sterile venting needle (LA038). Keep the bottle in an upright position preferably in a biological safety cabinet, place an alcohol swab over the rubber stopper and insert the venting needle with filter through it. Insertion and withdrawal of the needle should be done in a straight line. Discard the needle and mix the contents by gently inverting the bottle 2-3 times. Do not vent the bottle for anaerobic cultures. Incubate at 35-37°C for 18-24 hours. Recommended volume of blood to be tested in LQ029A: 8-10 ml (For Adult use).

### Principle And Interpretation

SS Agar medium is recommended as differential and selective medium for the isolation of *Salmonella* and *Shigella* species from pathological specimens (1). SS Agar is a moderately selective medium in which gram-positive bacteria are inhibited by bile salts, brilliant green and sodium citrate. Peptone, HM peptone B provides nitrogen and carbon source, long chain amino acids, vitamins and essential growth nutrients. Lactose is the fermentable carbohydrate. Brilliant green, bile salts and thiosulphate selectively inhibit gram-positive and coliform organisms. Sodium thiosulphate is reduced by certain species of enteric organisms to sulphite and H<sub>2</sub>S gas and this reductive enzyme process is attributed by thiosulphate reductase. Production of H<sub>2</sub>S gas is detected as an insoluble black precipitate of ferrous sulphide, formed upon reaction of H<sub>2</sub>S with ferric ions or ferric citrate, indicated in the center of the colonies.

The high selectivity of Salmonella Shigella Agar allows the use of large inocula directly from faeces, rectal swabs or other materials suspected of containing pathogenic enteric bacilli. On fermentation of lactose by few lactose-fermenting normal intestinal flora, acid is produced which is indicated by change of colour from yellow to red by the pH indicator-neutral red. Thus these organisms grow as red pigmented colonies. Lactose non-fermenting organisms grow as translucent colourless colonies with or without black centers. Growth of *Salmonella* species appears as colourless colonies with black centers resulting from H<sub>2</sub>S production. *Shigella* species also grow as colourless colonies which do not produce H<sub>2</sub>S.

This combination is recommended for paediatric use.

### Type of specimen

Clinical: blood

## Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (2,3). After use, contaminated materials must be sterilized by autoclaving before discarding.

## Warning and Precautions

InVitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/ eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.
2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.
3. Recommended for adult use only.
4. Further biochemical test must be carried out for confirmation.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

In a sterile glass bottle combination of broth and one agar coated surface.

<b>Colour of Agar medium</b>	<b>Colour of liquid medium</b>
Reddish orange coloured medium	Red coloured medium

### Quantity of medium

20ml of solid medium in glass bottle 40ml of liquid medium in glass bottle

<b>pH of Agar medium</b>	<b>pH of liquid medium</b>
6.80-7.20	6.80-7.20

### Sterility test Passes

release criteria

### Cultural response

Cultural characteristics was observed after incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth on agar medium	Growth on liquid medium	Colour of colony
<i>Escherichia coli</i> ATCC 25922 (00054*)	50-100	Fair	Fair	Pink with bile precipitate
# <i>Klebsiella aerogenes</i> ATCC 13048 (00175*)	50-100	Fair	Fair	Cream-Pink
<i>Salmonella</i> Typhimurium ATCC 14028 (00031*)	50-100	Luxuriant	Luxuriant	Colourless colonies with black center(H <sub>2</sub> S production)
<i>Salmonella</i> Typhi ATCC 6539	50-100	Luxuriant	Luxuriant	Colourless colonies with black center(H <sub>2</sub> S production)
<i>Salmonella</i> Enteritidis ATCC 13076 (00030*)	50-100	Luxuriant	Luxuriant	Colourless colonies with black center(H <sub>2</sub> S production)
<i>Shigella flexneri</i> ATCC 12022	50-100	Luxuriant	Luxuriant	Colourless
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	Luxuriant	Luxuriant	Colourless
<i>Salmonella</i> Choleraesuis ATCC 12011	50-100	Luxuriant	Luxuriant	Colourless colonies with black center(H <sub>2</sub> S production)
<i>Proteus mirabilis</i> ATCC 25933	50-100	Luxuriant	Luxuriant	colourless, may have black centre

Key : (\*) Corresponding WDCM numbers. # formerly known as *Enterobacter aerogenes*

## Storage and Shelf Life

On receipt store between 2-8°C. Use before expiry date on the label. Product performance is best if used within stated expiry period.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2,3).

## Reference

1. Lennette and others (Eds.), 1985, Manual of Clinical Microbiology, 4th ed., ASM, Washington, D.C.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.


Revision : 01/2023



HiMedia Laboratories Pvt. Limited,  
Plot No.C-40, Road No.21Y,  
MIDC, Wagle Industrial Area,  
Thane (W) -400604, MS, India



*In vitro* diagnostic  
medical device



Storage temperature



CEpartner4U, Esdoornlaan 13,  
3951DB Maarn, NL  
www.cepartner4u.eu



CE Marking



Do not use if  
package is damaged

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.