



MeReSa Agar

M1594

Intended use

Recommended for the selection, isolation and identification of Methicillin Resistant *Staphylococcus aureus* from clinical specimens.

Composition**

Ingredients	g / L
Tryptone	10.000
HM peptone B #	5.000
Glycine	10.000
Sodium pyruvate	10.000
Lithium chloride	5.000
Mannitol	10.000
Sodium chloride	10.000
Indicator mixture	0.130
Agar	20.000
Final pH (at 25°C)	7.1±0.2

**Formula adjusted, standardized to suit performance parameters

Equivalent to beef extract

Directions

Suspend 40.06 grams in 500 ml purified/ distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add sterile rehydrated contents of 1 vial of MeRS Selective Supplement (FD229) and CF Selective Supplement II (FD259) both in combination for more selectivity as desired. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Staphylococcus aureus sometimes referred to as “Staph” is a common bacterium found on the skin of healthy people. It is responsible for infections ranging from superficial to systemic (1,2). *Staphylococcus aureus* resistant to the antibiotic methicillin are referred to as Methicillin Resistant *Staphylococcus aureus* (MRSA) (3). Initially staphylococcal infections were treated using penicillin. But over the years, resistance to penicillin developed, so methicillin was the next drug of choice. Unfortunately certain strains (MRSA) have now developed resistance to methicillin also. Patients with breaks in their skin due to wounds, indwelling catheters or burns are those with certain risk of developing MRSA infection (4). Symptoms in serious cases may include fever, lethargy and headache. MRSA can cause UTI, pneumonia, toxic shock syndrome and even death, which can be controlled to a great extent by maintaining personal hygiene after interaction with an MRSA infected person (3). MRSA were recognized in 1980s as a major clinical and epidemiological problem. MRSA strains were heterogeneous in their expression of resistance to β -lactam agents, in that large differences in the degree of resistance were seen among the individual cells in a population. The basis of methicillin-resistance is the production of an additional penicillin-binding protein mediated by the *mec A* gene, an additional gene found in methicillin-resistant *Staphylococci*. MeReSa Agar Base was developed to detect the presence of the *mec A* gene in *S.aureus* i.e. methicillin-resistant *S.aureus*.

Tryptone and HM peptone B provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients. Lithium chloride and methicillin inhibit most of the contaminating microflora except methicillin-resistant *S.aureus* (MRSA). Glycine and sodium pyruvate enhance the growth of *Staphylococcus* species. Colour of the colonies is due to the indicator mixture and mannitol used in the medium. Sodium chloride maintains the osmotic equilibrium of the medium as well as supports the growth of *Staphylococcus* species.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (5,6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. MIC must be carried out for confirmation.
2. Further biochemical tests must be carried out for further identification

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% Agar gel.

Colour and Clarity of prepared medium

Pale pink coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 8.01% w/v aqueous solution at 25°C. pH : 7.1±0.2

pH

6.90-7.30

Cultural Response

Cultural characteristics observed with added MeRS Selective Supplement (FD229) and CF Selective Supplement II (FD259) both in combination after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth w/ FD229 & FD259	Recovery w/ FD229 & FD259	Colour of Colony
<i>Escherichia coli</i> ATCC 25922 (00013*)	≥10 ⁴	inhibited	0%	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	≥10 ⁴	inhibited	0%	
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	50-100	good-luxuriant	≥50%	light pink
<i>Staphylococcus epidermidis</i> ATCC 12228 (00036*)	50-100	Inhibited	0%	
<i>Staphylococcus gallinarum</i> MTCC 2992	50-100	Inhibited	0%	
<i>Staphylococcus saprophyticus</i> subsp. <i>saprophyticus</i> ATCC 15305 (00159*)	50-100	Inhibited	0%	

Key : *Corresponding WDCM numbers.

Storage and Shelf life

Store dehydrated powder and prepared medium at 2-8°C. Use before expiry period on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

- 1.Doyle, Beuchat and Montville, (Eds.), 1997, Food Microbiology Fundamentals and Frontiers. American Society for Microbiology, Washington, D.C.
- 2.Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- 3.Methicillin Resistant *Staphylococcus aureus*, Copyright © 1997-2005, Canadian Centre for Occupational Health and Safety, Sept 19th, 2005.
- 4.Dr. Alan Johnson, Methicillin resistant *Staphylococcus aureus* (MRSA) infection, The support group for MRSA sufferers and Dependents, AUG 1st, 2005.
- 5.Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- 6.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.


Revision : 06/2024




HiMedia Laboratories Pvt. Limited,
Plot No.C-40, Road No.21Y,
MIDC, Wagle Industrial Area,
Thane (W) -400604, MS, India



**In vitro diagnostic
medical device**



Storage temperature



CEpartner4U, Esdoornlaan 13,
3951DB Maarn, NL
www.cepartner4u.eu



CE Marking



**Do not use if
package is damaged**

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.