



## BHI w/PABA

M212

### Intended Use:

Recommended for examination of blood from patients under Sulphonamide therapy.

### Composition\*\*

Ingredients	g / L
HM infusion from ##	12.50
BHI Powder#	250.000
Peptone	10.000
Sodium chloride	5.000
Dextrose (Glucose)	2.000
Disodium hydrogen phosphate	2.500
p-Amino benzoic acid (PABA)	0.050
Final pH ( at 25°C)	7.4±0.2

## Equivalent to Calf brain infusion from

# Equivalent to Beef heart infusion from

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 37.05 grams in 1000 ml purified / distilled water. Heat if necessary to dissolve the medium completely. Dispense as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Brain Heart Infusion w/ PABA is a highly nutritious medium which can support luxuriant growth of wide variety of microorganisms including bacteria, yeasts and moulds (3) and is often used for isolation of pathogens from clinical specimens especially blood (5).

Para amino benzoic acid is an active inhibitor of the bacteriostasis produced by the sulfonamide drugs; also it serves as an accessory growth factor for several species of bacteria (4). Therefore para amino benzoic acid incorporated in the medium helps to neutralize the effect of antimicrobials present in the blood of patients under sulphonamide therapy making isolation of organisms from blood easier.

Peptone and HM infusion from and BHI Powder provides carbon, nitrogen, amino acids and vitamins. Dextrose serves as a source of energy. Sodium chloride helps in maintaining the osmotic equilibrium.

### Type of specimen

Clinical samples -blood.

### Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5). After use, contaminated materials must be sterilized by autoclaving before discarding.

### Warning and Precautions

In vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves protective clothing eye protection face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

### Limitations

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.
2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.
3. Further biochemical tests must be carried out for complete identification.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

Cream to yellow homogeneous free flowing powder

### Colour and Clarity of prepared medium

Light amber coloured, clear to very slightly opalescent solution without any precipitate

### Reaction

Reaction of 3.7% w/v aqueous solution at 25°C. pH : 7.4±0.2

### pH

7.20-7.60

### Cultural Response

Cultural characteristics observed with added 0.5 grams of sulphadiazine per litre after an incubation i) Bacteria at 35-37°C ii) Bacteroides species anaerobically for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth
<i>Bacteroides fragilis</i> ATCC 25285	50-100	good-luxuriant
<i>Candida albicans</i> ATCC 10231 (00054*)	50-100	good-luxuriant
<i>Neisseria meningitidis</i> ATCC 13090	50-100	good-luxuriant
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	good-luxuriant
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	good-luxuriant

Key : \*Corresponding WDCM numbers.

## Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-3 °C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

## Reference

1. MacFaddin J. F., 1985, Media for the Isolation-Cultivation-Identification- Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore
2. Murray P. R., Baron E. J., Jorgensen J. H., Pfaller M. A., Tenover F. C., Tenover J. C., (Eds.), 8th (Eds.), 2003, Manual of Clinical Microbiology, ASM, Washington, D.C.
3. Mirick G. S., 1943, Exp. Med., 78:255
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



HiMedia Laboratories Pvt. Limited,  
Plot No.C-40, Road No.21Y,  
MIDC, Wagle Industrial Area,  
Thane (W) -400604, MS, India




CEpartner4U, Esdoornlaan 13,  
3951DB Maarn, NL  
[www.cepartner4u.eu](http://www.cepartner4u.eu)



**In vitro diagnostic  
medical device**



**CE Marking**



**Storage temperature**



**Do not use if  
package is damaged**

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.