



Technical Data

Stuart Transport Medium (Transport Medium, Stuart)

M306

Intended use

Recommended for the preservation and transportation of *Neisseria* species and other fastidious organisms from the clinic to laboratory.

Composition**

Ingredients	g / L
Sodium glycerophosphate	10.000
Sodium thioglycollate	1.000
Calcium chloride anhydrous	0.100
Methylene blue	0.002
Agar	3.000
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 14.1 grams in 1000 ml double purified / distilled water. Heat to boiling to dissolve the medium completely. Dispense into tubes with screw caps to give a depth of approximately 7 cm. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes and after sterilization tighten the caps. Cool the tubes immediately in an upright position. Care should be taken that the water is free from chlorine.

Principle And Interpretation

Stuart Transport media were originally designed by Stuart while studying Gonococci (1). Stuart et al (2) later on modified the Stuart Medium for the transportation of gonococcal specimens for culturing. Ringertz included thioglycollate in the Stuart Medium and omitted charcoal (3). The medium may be used for the transportation of many fastidious organisms including anaerobes by maintaining the organism's viability without significant multiplication (4). Crooks and Stuart (5) suggested the addition of Polymyxin B sulphate which facilitates the recovery of *Neisseria gonorrhoeae*.

This medium is a chemically defined, semisolid, non-nutrient medium which prevent microbial proliferation. Because of this composition the medium ensures that microorganisms present are able to survive for a sufficiently long period of time. The medium provides an adequate degree of anaerobiosis which can be monitored by means of the redox indicator methylene blue. Prepared sterile medium will undergo a slight degree of oxidation at the upper periphery of the medium, however, in the tube or vial exhibits a distinct blue colour throughout the medium, it should be discarded. Calcium chloride along with sodium glycerophosphate act as good buffering agent and also maintains osmotic equilibrium in the medium.

Type of specimen

Clinical samples - Gonococcal specimens.

Specimen Collection and Handling:

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (2,3).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Due to nutritional variations, some strains may show poor growth.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Please refer disclaimer Overleaf.

Quality Control

Appearance

White to light blue coloured homogeneous free flowing powder

Gelling

Semisolid, comparable with 0.3% Agar gel.

Colour and Clarity of prepared medium

Colourless to whitish coloured slightly opalescent butt with upper 10% or less portion blue on standing.

Reaction

Reaction of 1.41% w/v aqueous solutions at 25°C. pH : 7.4±0.2

pH

7.20-7.60

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 72 hours when subcultured from Stuart Transport Medium.

Organism	Growth	Subculture Medium
<i>Haemophilus influenzae</i> ATCC 49247	good	Chocolate Agar (incubated in CO ₂ atmosphere)
<i>Neisseria gonorrhoeae</i> ATCC 19424	good	Chocolate Agar (incubated in CO ₂ atmosphere)
<i>Streptococcus pneumoniae</i> ATCC 6303	good	Tryptone Soya Agar with 5% sheep blood

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 5-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (6,7).


Reference

1. Stuart, 1946, Glasgow Med. J. 27:131.
2. Stuart, Toshach and Patsula, 1954, Can. J. Public Health, 45:73.
3. Ringertz, 1960, Acta Pathol. Microbiol. Scand., 48:105.
4. Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.
5. Crookes E.M.L. and Stuart R.D., 1959, J. Path. Bact., 78:283.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 05/2024



HiMedia Laboratories Pvt. Limited,
Plot No.C-40, Road No.21Y,
MIDC, Wagle Industrial Area,
Thane (W) -400604, MS, India




CEpartner4U, Esdoornlaan 13,
3951DB Maarn, NL
www.cepartner4u.eu



**In vitro diagnostic
medical device**



CE Marking



Storage temperature



**Do not use if
package is damaged**

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.