

MB568

**HiPurA® Food Pathogen DNA Purification Kit
[(Salmonella) Provided with Culture Media]**

Kit Contents

| Product Code | Reagents provided | MB568 | |
|--------------|---|----------|----------|
| | | 10 Preps | 25 Preps |
| DS0015 | Lysis Solution (AL) | 3 ml | 7.5 ml |
| DS0010 | Lysis Solution (C1) | 3 ml | 7.5 ml |
| DS0031 | Prewash Solution (PWB) | 6 ml | 15 ml |
| DS0012 | Wash Solution Concentrate (WS) | 2 ml | 5 ml |
| DS0040 | Elution Buffer (ET) [10mM Tris- Cl, pH8.5] | 3 ml | 7.5 ml |
| MBO86 | Proteinase K | 5 mg | 12.5 mg |
| DS0003 | RNase A Solution (20mg/ml) | 250 µl | 625 µl |
| M1494I | Buffered Peptone Water | 56 g | 140 g |
| DSCA02 | HiShredder (in DBCA016 Collection tube) | 10 nos. | 25 nos. |
| DBCA03 | HiElute Miniprep Spin Column (Capped) [in DBCA016 Collection Tube] | 10 nos. | 25 nos. |
| DBCA016 | Collection Tube (Uncapped), Polypropylene (2.0ml) | 10 nos. | 25 nos. |
| PW1139 | Collection Tube, Polypropylene (2.0ml) | 20 nos. | 50 nos. |

Intended Use

Recommended for isolation of DNA from Food samples enriched in culture media.

Introduction:

Salmonellosis, a widely distributed foodborne disease, is a major public health concern worldwide resulting in thousands of deaths. This food-borne disease is caused by *Salmonella*, a rod shaped, gram-negative non-spore forming bacterium, a member of Enterobacteriaceae family. The mode of infection is through contaminated food. The traditional culture based detection takes longer for detection of *Salmonella*, hence, need for new, quick and sensitive methods to detect *Salmonella spp.* is a major concern for food industry.

HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media] is developed to provide a fast and easy method for purification of total DNA from variety of food samples (after enrichment according to AOAC guidelines) for reliable applications in PCR of *Salmonella spp.* The DNA purification procedure using the miniprep spin column comprises of three steps viz. adsorption of DNA to the membrane, removal of residual contaminants and elution of pure genomic DNA. HiMedia's HiElute Miniprep Spin Column (Capped) format allows rapid processing of multiple samples. The columns have a high binding capacity and high quality genomic DNA is obtained from various species. The DNA obtained is compatible with downstream applications like PCR. HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media]

HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media]

This kit simplifies isolation of DNA from food samples by the spin-column procedure. The food sample is enriched in an enrichment medium (provided) as per AOAC guidelines and the cells are harvested by centrifugation. After harvesting, the bacterial cell wall is degraded by Proteinase K. Following lysis, the DNA is bound to the silica-gel membrane of the HiElute Miniprep Spin Column

(Capped) to yield pure DNA. Two rapid wash steps remove trace amount of salt and protein contaminants resulting in the elution of pure DNA in the Elution Buffer provided with the kit.

HiElute Miniprep Spin Column (Capped) [DBCA03]

HiElute Miniprep Spin Column (Capped) is based on the advanced silica binding principle presented in a microspin format. The system efficiently couples the reversible nucleic acid-binding properties of the advanced silica gel membrane and the speed plus versatility of spin column technology to yield high quantity of DNA. The use of spin column facilitates the binding, washing and elution steps thus enabling multiple samples to be processed simultaneously. This column eliminates the need for alcohol precipitation, expensive resins and harmful organic compounds such as phenol and chloroform, otherwise employed in traditional DNA isolation techniques. DNA binds specifically to the advanced silica-gel membrane while contaminants pass through. PCR inhibitors such as divalent cations and proteins are completely removed in two efficient wash steps, leaving pure nucleic acid to be eluted in the buffer provided with the kit. The purified DNA can be used for further downstream applications.

Elution

The yield of genomic DNA depends on the sample type and the number of cells in the sample. A single elution with 200 µl of Elution Buffer (ET) will provide sufficient DNA to carry out multiple amplification reaction. Elution with volume less than 200 µl will increase the final DNA concentration, but will reduce the overall DNA yield.

Materials needed but not provided

- 55°C water bath or heating block
- Tabletop Micro centrifuge (with rotor for 2.0 ml tubes)
- Ethanol (96 - 100%)
- Micropipettes and Tips
- Distilled water
- Autoclave
- Incubator Shaker

Storage

Store the HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media] between 15-25°C except certain components as specified on each labels. Under recommended condition kit is stable for 1 year.

General Preparation Instructions

1. Preheat a water bath or heating block to 55°C.
2. **M1494I Media preparation:** Suspend 20.07 grams (the equivalent weight of dehydrated medium per liter) in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.
3. **Thoroughly mix reagents**
Examine the reagents for precipitation. If any kit reagent forms a precipitate (other than enzymes), warm at 55-65°C until the precipitate dissolves and allow cooling to room temperature (15-25°C) before use.
4. Ensure that clean & dry tubes and tips are used for the procedure.

5. **Dilute Wash Solution Concentrate (WS) (DS0012) as follows:**

| Number of Preps | Wash Solution Concentrate (WS) | Ethanol (96-100%) |
|-----------------|--------------------------------|-------------------|
| 10 | 2 ml | 8 ml |
| 25 | 5 ml | 20 ml |

6. **Reconstitute Proteinase K (MB086)**

The HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media] contains Proteinase K. It is completely free of DNase and RNase activity. The specific activity of Proteinase K is 33.5 units/mg dry weight.

Re-suspend the Proteinase K (MB086) powder in Molecular Biology Grade Water (ML024) to obtain a 20 mg/ml stock solution.

| Number of Preps | Proteinase K | Molecular Biology Grade Water |
|-----------------|--------------|-------------------------------|
| 10 | 5 mg | 0.25 ml |
| 25 | 12.5 mg | 0.62 ml |

The product as supplied is stable at room temperature (15-25°C); upon receipt, store at 2-8°C & upon reconstitution store at -20°C.

NOTE: The Proteinase K solution must be added directly to each sample preparation every time. Do not combine the Proteinase K and Lysis Solution for storage.

RNase A enzyme treatment

RNase A is a type of RNase that is commonly used in research. RNase A (e.g. bovine pancreatic ribonuclease A) is one of the sturdiest enzymes in common laboratory usage. It cleaves 3' end of unpaired C and U residues.

Unit Definition for RNase A

One unit of the enzyme causes an increase in absorbance of 1.0 at 260 nm when yeast RNA is hydrolyzed at 37°C and pH 5.0. Fifty units are approximately equivalent to 1 Kunitz unit. It is completely free of DNases and proteases. The specific activity is 90 U/mg.

The product as supplied is stable at room temperature (15-25°C).

Centrifugation

All centrifugation steps are carried out in conventional laboratory e.g. Beckman CS-6KR, Heraeus Varifuge 3.0R, or Sigma 6k10 with fixed angle rotor. The tubes provided with the kit are compatible with almost all laboratory centrifuges and rotors. All centrifugation steps are performed at room temperature and are given in g, the correct rpm can be calculated using the formula:

$$RPM = \sqrt{RCF/1.118} \times 10^5 r$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary g-force.

Specimen Collection and Handling

For Food

Collect food sample in a sterile container and freeze the sample at -20°C for short term storage or -80°C for long term storage.

Types of Specimen

Samples: Food

Enrichment of pathogens:

In order to ensure sensitive detection of pathogens from different variety of food products by PCR, the pathogens need to be enriched in broth.

Procedure:

- Weigh 25g of food sample and add it to autoclaved 225 ml of the Buffered peptone water (M1494I) (Refer General Preparation Instructions).
- Incubate at 37°C for 16-18 hours with shaking at 100 rpm.
- Take 1.5 ml of above sample and proceed as follows:

1. Harvesting of cells

Pellet 1.5 ml of an overnight enriched sample in provided 2ml capped collection tube by centrifuging for 2 minutes at 12,000-16,000 x g (approximately 13,000-16,000 rpm). Remove the culture medium and discard.

2. Re-suspend the pellet

Re-suspend the pellet thoroughly in 180 µl of Lysis Solution (AL) (DS0015).

3. Prepare for cell lysis

Add 20 µl of the Proteinase K solution (20mg/ml) (**Refer to General Preparation Instructions**) to the sample. Mix and incubate for 5 minutes at 55°C. If residual RNA is not a concern, continue with step 4.

Optional RNase A treatment

If RNA-free genomic DNA is required, add 20 µl of RNase A Solution (DS0003), mix and incubate for 5 minutes at room temperature (15-25°C), then continue with step 4.

4. Lyse cells

Add 200 µl of Lysis Solution (C1) (DS0010), vortex thoroughly (about 15 seconds) and incubate at 55°C for 10 minutes.

NOTE: A homogeneous mixture is essential for efficient lysis.

5. Load sample in HiShredder (DSCA02)

Add the entire sample to the HiShredder placed in a 2.0 ml collection tube (uncapped) and centrifuge for 1 minute at a maximum speed (≈14,000 rpm). Transfer the flow-through fraction to a 2.0 ml collection tube (not provided) without disturbing the cell debris pellet.

6. Prepare for binding

Add 200 µl of ethanol (95-100%) to the lysate and mix thoroughly by vortexing for few seconds.

NOTE: A homogenous mixture is essential. A white precipitate may form on addition of ethanol. It is essential to apply all of the precipitate to the HiElute Miniprep Spin Column. This precipitate does not interfere with the DNA isolation procedure or with any subsequent application. Do not use alcohols other than ethanol because this may result in reduced yields.

7. Load lysate onto HiElute Miniprep Spin Column (Capped) [DBCA03]

Transfer the lysate obtained from step 6 onto HiElute Miniprep Spin Column (Capped) provided. Centrifuge at $\geq 6,500$ x g (≈10,000 rpm) for 1 minute. Discard the flow-through liquid and place the spin column in same 2.0 ml collection tube.

NOTE: Use a wide bore pipette tip to reduce shearing of the DNA while transferring contents

onto the column. It is essential to apply all of the precipitate to the column. If the solution has not completely passed through the membrane, spin at a higher speed until all the solution has passed through. Centrifugation at full speed will not affect the yield or purity of the DNA.

8. Prewash

Add 500 µl of Prewash Solution to the column and centrifuge at $\geq 6,500 \times g$ ($\approx 10,000$ rpm) for 1 minute. Discard the flow-through liquid and re-use the same collection tube with the column.

9. Wash

(Prepare Wash Solution Concentrate (WS) (DS0012) as indicated in General Preparation Instructions)

Add 500 µl of diluted Wash Solution (WS) to the column and centrifuge for 3 minutes at maximum speed $12,000-16,000 \times g$ ($\approx 13,000-16,000$ rpm). Discard the flow through and spin again at same speed for the additional 1 minute to dry the column. The column must be free of ethanol before eluting the DNA.

10. DNA Elution

Transfer the HiElute Miniprep Spin Column (Capped) to fresh uncapped collection tube. Pipette 200 µl of the Elution Buffer (ET) (DS0040) directly into the column without spilling to the sides. Incubate for 1 minute at room temperature. Centrifuge at $\geq 6,500 \times g$ ($\approx 10,000$ rpm) for 1 minute to elute the DNA.

NOTE: To increase the elution efficiency, incubate for 5 minutes at room temperature after adding the Elution Buffer (ET), then centrifuge. Elution with volumes less than 200 µl increases the final DNA concentration in the eluate significantly, but slightly reduces the overall DNA yield. Storing DNA in water can cause acid hydrolysis.

11. Transfer the eluate to a fresh capped 2ml collection tube for longer storage.

Storage of the eluate with purified DNA: The eluate contains pure genomic DNA. For short-term storage (24-48 hrs) of the DNA, 2-8°C is recommended. For long-term storage, -20°C or lower temperature (-80°C) is recommended. Avoid repeated freezing and thawing of the sample which may cause denaturing of DNA. The Elution Buffer will help to stabilize the DNA at these temperatures.

Warning and Precautions

Not for Medicinal Use. Read the procedure carefully before beginning the protocol. Wear protective gloves/protective clothing/eye protection/face protection. Follow good laboratory practices while handling samples. Standard precautions should be followed as per established guidelines. Safety guidelines may be referred in safety data sheets of the product.

Limitations

1. The yield of DNA depends upon the type and the volume of starting material used.

Performance and Evaluation

Each lot of HiMedia's HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media] is tested against predetermined specifications to ensure consistent product quality.

Quality Control

| Type of Sample | DNA Purity |
|----------------|------------|
| Fruit juice | 1.6-1.9 |

Troubleshooting Guide:

| Sr. No. | Problem | Possible Cause | Solution |
|---------|---|---|--|
| 1. | HiElute Miniprep Spin Capped Column is clogged | Sample volume is large | Use smaller quantity of sample; to salvage the current preparation, clogging can be alleviated by increasing the g force or spinning for longer time until the lysate passes through the spin column. |
| 2. | Lysate appears to be gelatinous prior to loading onto the column | Sample volume is large | Use fewer cells ($\leq 1 \times 10^{10}$ cells/ml). The incubation time and or the amount of Proteinase K solution can be increased Eg: The incubation time and the amount of enzyme can be doubled. |
| 3. | Poor / Lower yield of genomic DNA | Sample is old | Yield of genomic DNA varies from different species and strains of bacteria. It is necessary to use cells before they reach their maximum density or they become confluent. |
| | | Incomplete lysis of cells | The incubation time and or the amount of Proteinase K solution can be increased.eg: The incubation time and the amount of enzyme can be doubled. |
| | | Lysate/Ethanol mixture is not homogenous | Vortex the tubes for at least 5-10 seconds in order to obtain a homogenous solution before applying it to the column. |
| | | DNA elution is incomplete. Eluate contains residual ethanol from the wash | DNA Yield can be improved by incubating the Elution Buffer for 5 minutes at room temperature (15-25°C) after it is added to the column. Ethanol from the final wash should be eliminated completely before eluting DNA. Spin the tubes for longer time to dry the column completely. |
| | | Wash Solution Concentrate was not diluted before use | Check that the Wash Solution Concentrate is properly diluted with ethanol as per instructions. |
| | | Use of water instead of Elution Buffer for elution of DNA. | Elution Buffer is recommended for optimal yields and storage of the genomic DNA. If water is used instead of the Elution Buffer the pH should be at least 7.0 to avoid acidic conditions, which may cause acid hydrolysis of DNA when stored for long periods of time. (NOTE: Only DNase/RNase and Protease free water should be used for eluting DNA) |
| 4. | Purity of the DNA is lower than expected; A_{260}/A_{280} ratio is | Background reading is high due to silica fines | The DNA sample can be centrifuged at maximum speed for 1 minute, the supernatant can be used to repeat the absorbance readings. |

| | | | |
|----|---|--|---|
| | low | Sample diluted in water | Use either Elution Buffer provided, or (10 mM Tris-HCl, 0.5 mM EDTA pH 9.0) or 10 mM Tris-HCl pH 8.0-8.5 as the eluant. |
| 5. | Purity of the DNA is lower than expected; A_{260}/A_{280} ratio is too high. | RNA contamination | RNase A treatment should be included in future isolations or the final product can be treated with RNase A and repurified. |
| 6. | DNA is sheared | Improper handling of genomic DNA | All pipetting steps should be executed as gently as possible. Wide orifice pipette tips are recommended to eliminate shearing of the DNA to a large extent. If the isolated DNA is to be used for PCR, mix with gentle pipetting or invert until homogenous, instead of vortexing as it reduces shearing of DNA considerably. |
| | | Cells are old | Cells grown for a longer time period may lyse prematurely when subjected to cell wall lysing enzymes, which may result in the release of endogenous nucleases and subsequent DNA degradation. |
| 7. | Downstream applications are inhibited | Traces of ethanol present in the final genomic DNA preparation | After the washing steps, the eluate should not come in contact with the column. Spin the column for 1 minute at maximum speed (12,000-16,000 x g) if necessary, after emptying the collection tube. |
| | | Salt is carried over in the final genomic DNA preparation | The HiElute Miniprep Spin Capped Column should be transferred to a new 2.0 ml collection tube before adding the wash solution. |

Safety Information

HiPurA® Food Pathogen DNA Purification Kit [(Salmonella) provided with culture media] is for laboratory use only; not for drug, household or other uses. Take appropriate laboratory safety measures and wear gloves when handling. Avoid contact with skin, and use eye protection. In case of contact, wash with large amount of water. Seek medical attention. Not compatible with disinfecting agents containing bleach. Please refer the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.


Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed off in accordance with current laboratory techniques.

Please refer disclaimer Overleaf.

Technical Assistance:

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No.
21Y, MIDC, Wagle Industrial Estate,
Thane, (West) 400604, Maharashtra,
INDIA. Web: www.himedialabs.com



08/2027

PIMB568_0/0824

MB568-05

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com