

## CNA Agar Plate with 5% sheep blood

MP5208

### Intended Use:

Recommended for selective isolation of pathogenic gram positive cocci from clinical and non clinical specimens.

### Composition\*\*

Ingredients	g / L
Biopeptone	20.000
Tryptose B #	3.000
Corn starch	1.000
Sodium chloride	5.000
Colistin sulphate	0.010
Nalidixic acid	0.015
Agar	15.000
Sheep Blood	50.000

\*\*Formula adjusted, standardized to suit performance parameters

# Equivalent to Tryptic digest of beef heart

### Directions

Either streak, inoculate or surface spread the test inoculum (50-100CFU) aseptically on the plate.

### Principle And Interpretation

CNA Agar is a nutritionally rich formula containing 5% defibrinated blood, which provides more nutrients and capability of displaying haemolytic reactions. Columbia Blood Agar Base is utilized as a base for preparation of media containing blood and in selective media preparations where various combinations of antimicrobial agents are used as additives. Ellner et al formulated the medium (1) and found that the combination of peptones used gave more rapid and abundant growth of Streptococci, Staphylococci, *Neisseria* and *Haemophilus* with better-defined haemolytic reactions. Columbia C.N.A. Agar Base is prepared with the same formula as Columbia Agar Base with the addition of 10 mg/litre of colistin and 15 mg/litre of nalidixic acid to inhibit the growth of gram-negative bacteria and to support the growth of Staphylococci, haemolytic Streptococci and Enterococci when supplemented with 5% blood.

Biopeptone and tryptose B supports luxuriant growth of microorganisms and visualization of good haemolytic reactions. Sheep blood allows detection of haemolytic reactions and supplies X-factor necessary for the growth of many bacterial species. Horse blood supplies X-factor and V-factor, therefore is mostly preferred in most laboratories. Yeast extract and cornstarch serve as energy source and neutralizer respectively.

It should be noted that this medium has relatively high carbohydrate content and, therefore, beta-hemolytic streptococci may produce a greenish hemolytic reaction that may be mistaken for alpha haemolysis. The addition of the antimicrobial agents, colistin (or polymyxin B) and nalidixic acid, renders the medium selective for gram-positive microorganisms (2). Colistin and nalidixic acid disrupt the cell membrane of gram-negative organisms, whereas nalidixic acid blocks DNA replication in susceptible gram-negative bacteria (3).

C.N.A. Agar with addition of blood gives selective isolation of gram-positive cocci, Staphylococci and Streptococci, particularly when gram-negative bacilli are present and tend to overgrow on conventional blood agar plates.

### Type of specimen

Clinical samples : Pure cultures isolated from urine, stool etc.

### Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5). After use, contaminated materials must be sterilized by autoclaving before discarding.

### Warning and Precautions

In Vitro diagnostic use only. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium
2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended o users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.
3. Further serological and biochemical testing is required for complete identification.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

Sterile Columbia C.N.A Agar w/5% sheep blood in 90 mm disposable plates with smooth surface and absence of black particles/cracks/bubbles

### Colour of medium

Red coloured medium

### Quantity of medium

25ml of medium in 90mm disposable plate

### pH

7.10- 7.50

### Sterility Check

Passes release criteria.

### Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 40-48 hours .

Organism	Inoculum (CFU)	Growth	Recovery	Haemolysis
<i>Escherichia coli</i> ATCC 25922 (00013*)	$\geq 10^3$	inhibited	0%	
<i>Neisseria meningitidis</i> ATCC 13090	$\geq 10^3$	inhibited	0%	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	luxuriant	$\geq 50\%$	beta/gamma
<i>Staphylococcus epidermidis</i> ATCC 12228 (00036*)	50-100	luxuriant	$\geq 50\%$	gamma
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	luxuriant	$\geq 50\%$	alpha
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	luxuriant	$\geq 50\%$	beta

Key : (\*) Corresponding WDCM numbers

## Storage and Shelf Life

On receipt store between 2-8°C. Use before expiry date on the label. Product performance is best if used within stated expiry period.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

## Reference


1. Ellner et al, 1966, Am. J. Clin. Path., 45:502.
2. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
3. Estevez, 1984, Lab. Med., 15:258.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



HiMedia Laboratories Pvt. Limited,  
Plot No.C-40, Road No.21Y,  
MIDC, Wagle Industrial Area,  
Thane (W) -400604, MS, India



*In vitro* diagnostic  
medical device



Storage temperature



CEpartner4U, Esdoornlaan 13,  
3951DB Maarn, NL  
[www.cepartner4u.eu](http://www.cepartner4u.eu)



CE Marking



Do not use if  
package is damaged

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.