

BHI HiVeg™ w/PABA

MV212

Intended Use:

Recommended for examination of blood from patients under sulphonamide therapy.

Composition**

| Ingredients | Gms / Litre |
|-----------------------------|-------------|
| HiVeg™ infusion | 10.000 |
| HiVeg™ special infusion | 7.50 |
| HiVeg™ peptone | 10.000 |
| Sodium chloride | 5.000 |
| Dextrose (Glucose) | 2.000 |
| Disodium hydrogen phosphate | 2.500 |
| p-Amino benzoic acid (PABA) | 0.050 |
| Final pH (at 25°C) | 7.4±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 37.05 grams in 1000 ml purified/ distilled water. Heat if necessary to dissolve the medium completely. Dispense into tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

BHI HiVeg™ w/PABA is a highly nutritious medium which can support luxuriant growth of wide variety of microorganisms including bacteria, yeasts and moulds (1) and is often used for isolation of pathogens from clinical specimens especially blood (2). Para amino benzoic acid is an active inhibitor of the bacteriostasis produced by the sulfonamide drugs; also it serves as an accessory growth factor for several species of bacteria (3). Therefore para amino benzoic acid incorporated in the medium helps to neutralize the effect of antimicrobials present in the blood of patients under sulphonamide therapy making isolation of organisms from blood easier. BHI HiVeg™ w/PABA is prepared by using vegetable peptones in place of animal based peptones which make the media free of BSE/TSE risks. HiVeg™ infusion and HiVeg™ special infusion and HiVeg™ peptone provides carbon, nitrogen, amino acids and vitamins. Dextrose serves as a source of energy. Sodium chloride helps in maintaining the osmotic equilibrium.

Type of specimen

Please add specimens

Specimen Collection and Handling

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Further biochemical tests must be carried out for complete identification.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured, clear to very slightly opalescent solution without any precipitate

Reaction

Reaction of 3.7% w/v aqueous solution at 25°C. pH : 7.4±0.2

pH

7.20-7.60

Cultural Response

Cultural characteristics observed with added 0.5 grams of sulphadiazine per litre after an incubation i) Bacteria at 35-37° C for 18-24 hours ii)Fungal at 25-30°C for 24-48 hours. iii)Bacteroides species anaerobically for 24-48 hours.

| Organism | Inoculum (CFU) | Growth |
|---|----------------|----------------|
| <i>Bacteroides fragilis</i> ATCC 25285 | 50-100 | good-luxuriant |
| <i>Candida albicans</i> ATCC 10231 (00054*) | 50-100 | good-luxuriant |
| <i>Neisseria meningitidis</i> ATCC 13090 | 50-100 | good-luxuriant |
| <i>Streptococcus pneumoniae</i> ATCC 6303 | 50-100 | good-luxuriant |
| <i>Streptococcus pyogenes</i> ATCC 19615 | 50-100 | good-luxuriant |

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

Reference

1. MacFaddin J. F., 1985, Media for the Isolation-Cultivation-Identification- Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.
2. Murray P. R., Baron E. J., Jorgensen J. H., Pfaller M. A., Tenover F. C., Tenover F. C., (Eds.), 8th (Eds.), 2003, Manual of Clinical Microbiology, ASM, Washington, D.C.
3. Mirick G. S., 1943, Exp. Med., 78:255
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 02 / 2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.