



Technical Data

XLT4 Supplement

FD152

Recommended for isolation of non-typhi *Salmonella* .

Composition

Per vial sufficient for 1000 ml medium

Ingredients

(40ml per vial contains Tergitol 4 - 27%)

Tergitol 4

Distilled water

Concentration

10.800g

40ml

Directions:

Warm up XLT4 Supplement to room temperature and add 4.6ml to 1000 ml of XLT4 Agar Base [M1147](#) / XLT4 HiVeg™ Agar Base [MV1147](#) . Heat to boiling to dissolve the medium completely. Avoid overheating. DO NOT AUTOCLAVE. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 2 / 2017

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.