



# Technical Data

## Growth Supplement II for MSM

FD288

A growth supplement for the enrichment of *Salmonella* species.

### Composition

Per vial sufficient for 1000 ml medium

### Ingredients

Ethanolamine

### Concentration

2ml

### Directions:

Warm up the refrigerated contents of one vial to room temperature. Aseptically add the contents of one vial of Growth Supplement II for MSM (FD288) to 1000 ml sterile, cooled (45-50°C) MSM Broth Base [M1864](#) along with Growth Supplement I for MSM [FD287](#). Mix well and dispense as desired.

### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.