



Broth Medium D (Lactose Monohydrate Broth)

ME1003

Lactose Monohydrate Broth is used for the detection of coliform bacteria in water, foods, dairy products in accordance with per European Pharmacopoeia 2008.

Composition**

Ingredients	Gms / Litre
Pancreatic digest of gelatin	5.000
Beef extract	3.000
Lactose monohydrate	5.000
pH after sterilization (at 25°C)	6.9±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 12.75 grams (the equivalent weight of dehydrated powder per litre) in 1000 ml purified/distilled water. Heat if necessary to dissolve the medium completely. For larger inocula (10 ml or more), concentrated medium may be prepared to account for medium dilution by the inoculum. Dispense in tubes containing inverted fermentation vial (Durhams tube) as desired. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes and cool immediately.

Principle And Interpretation

Lactose Broth is recommended by European pharmacopoeia(4) for selective pre-enrichment of *Enterobacteriaceae* as well as for *E.coli* and *Salmonella* in water, food and pharmaceutical products. The medium is used for detection of specified microorganisms of non-sterile products according to European Pharmacopoeia (EP). This medium is also recommended by various other pharmacopoeia (5,6,7)

Lactose Broth is also recommended by APHA in the performance and confirmation of the presumptive test for coliform bacteria in water (1), food (2), milk (3) .This medium can be used as an alternate to Lauryl Sulphate Broth in the presumptive test of the MPN of standard coliforms.

Pancreatic digest of gelatin and beef extract supply essential nutrients to the organisms. Lactose monohydrate is a fermentable carbohydrate for the coliforms. Tubes of Lactose Broth are inoculated with dilutions of water or milk, etc. under test, and incubated at 35°C and examined for gas formation after 24 and 48 hours.

Members of the coliform group are defined as aerobic and facultative anaerobic gram-negative and non-sporing bacilli which ferment lactose with gas formation within 48 hours at 35°C. In testing dairy products, Lactose Broth is used only in the completed test (3). Large water samples may require double strength Lactose Broth to minimize the final volume.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and clarity of prepared medium

Light amber coloured clear solution

Reaction

Reaction of 1.27% w/v aqueous solution after sterilization. pH : 6.9±0.2

pH

6.70-7.10

Growth Promotion Test

As per European Pharmacopoeia.

Cultural Response

ME1003: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Gas
Cultural Response			
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	Positive reaction
<i>Enterococcus faecalis</i> ATCC 29212	50-100	luxuriant	Negative reaction
<i>Escherichia coli</i> ATCC 8739	50-100	luxuriant	Positive reaction
<i>Pseudomonas aeruginosa</i> ATCC 9027	50-100	luxuriant	Negative reaction
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	Negative reaction
<i>Salmonella Abony</i> NCTC 6017	50-100	luxuriant	Negative reaction

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2- 8°C. Use before expiry date on the label.

Reference

- 1.Eaton A. D., Clesceri L. S. and Greenberg A W.,(Eds.), 2008, Standard Methods for the Examination of Water and Wastewater, 21st ed., APHA, Washington, D.C.
- 2.Downes F P and Ito K(Eds.), 2001, Compendium of Methods For The Microbiological Examination of Foods, 4th ed., APHA, Washington, D.C
- 3.Wehr H M and Frank J H., 2004, Standard Methods for the Examination of Dairy Products, 17th ed., APHA Inc., Washington, D.C.
- 4.European Pharmacopoeia, 2008, European Department, for the quality of Medicines.
- 5.The United States Pharmacopoeia, 2008, The United States Pharmacopoeial Convention. Rockville, MD.
- 6.British Pharmacopoeia, 2009, The Stationery office British Pharmacopoeia
- 7.The Indian Pharmacopoeia 2007, Govt. of India, 2007. The Controller of Publication, Delhi.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.