

HiMedia's TabSpin Microcentrifuge LA1001



Outstanding Productivity at an exceptional quality

- German brushless motor with excellent performance & long life
- High speed with low sound
- 14,500 rpm with accuracy and precision
- Timer controlled
- Ideal for applications across pharmaceutical, biotechnology and academic research



Salient Features...

Specifications

Max Capacity	12 x 0.2 / 0.5 / 1.5 or 2 ml tubes (0.2ml & 0.5ml Tube Adapter)
Max Speed	14,500
Min Speed	500
Max RCF	14,000
Min RCF	100
Drive System	Brushless motor
Sound Level	56 dBA
Max Timer Range	20 seconds to 99 minutes + "Hold" (continuous mode)
Time Increments	20 second
Dimensions (H xWx D)	138 x 212 x 240
Weight	4.1 kg/ 9.3 lb

