

Adulteration Testing of Milk

Kit for Adulteration Testing of Milk

Milk is a part of the diet of every individual, its benefits being enumerable. It is a complete health drink. However the increased rate of milk adulteration has been a matter of concern, posing health hazards thus masking its advantages. Despite the checks on purity, adulteration continues to be on a rising stand.

The Kit for Adulteration Testing of Milk is a kit designed for the detection of common seventeen milk adulterants in the milk sample. The kit is available in two parts

K088A Adulteration Testing of Milk


Part A : Reagents for detection of 11 adulterants

1. Alizarine
2. Formaline
3. Urea
4. Starch
5. Neutralizers
6. Detergent
7. Sodium chloride
8. Skim Milk powder
9. Sugar
10. Glucose
11. Hydrogen Peroxide



Main features

- Kit comprising of 50 tests
- Simple, Easy and Rapid detection
- Most economical
- Can be performed even outside the laboratory
- Suitable for large scale as well as household use



K088B Adulteration Testing of Milk

Part B : Reagents for detection of 6 adulterants

1. Cellulose
2. Maltodextrin / Maltose
3. Ammonium Sulphate
4. Proteins
5. Boric Acid
6. Pond Water (Nitrate-Nitrogen)

