

AS115 Dichloromethane, HPLC**Product Number Packing**

AS115 : 2.5L
AS115 : 1L

Product Information

Product Code : AS115
Product Name : Dichloromethane, HPLC
Synonym : Methylene chloride
Molecular Formula : CH₂Cl₂
Molecular Weight : 84.93
CAS No. : 75-09-2
EC No. : 200-838-9
HS Code : 2903 12 00
EC Index No. : 602-004-00-3
Shelf Life : 3 years

Technical Specification

Appearance : Clear colorless liquid
Density (d 20/4) : 1.324 - 1.326 g/mL
Free acid (as HCl) : <= 0.001%
Water (K.F.) : <= 0.05%
Non-volatile matter : <= 0.001%
UV transmittance 1 : >= 20% (at 235 nm)
UV transmittance 2 : >= 80% (at 245 nm)
UV transmittance 3 : >= 98% (at 260 nm)
Assay (GC) : min. 99.80%

GHS Safety Information

Hazard Statement(s) : H351
Precautionary Statement(s) : P281
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 40
S-Phrase(s) : 23-24/25-36/37
WGK : 2
RTECS : PA8050000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein



Product Information

| | | |
|------------------|---|---------------|
| UN No. | : | 1593 |
| Class | : | 6.1 |
| Packaging Group | : | 3 |
| Marine Pollutant | : | YES |
| ADR/RID | : | 1593 6.1/PG 3 |
| IMDG | : | 1593 6.1/PG 3 |
| IATA | : | 1593 6.1/PG 3 |

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein