



Polymyxin B Selective Supplement

FD003

An antibiotic supplement recommended for the selective isolation of various microorganisms.

Composition

Per vial sufficient for 500/ 1000 ml medium

*Ingredients

Polymyxin B sulphate

Concentration

50000Unit

Directions:

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 475 ml of sterile, molten Bacillus Cereus Agar Base [M833](#) / Bacillus Cereus HiVeg™ Agar Base [MV833](#) / or to 450 ml of KG Agar Base [M658](#) / KG HiVeg™ Agar Base [MV658](#) / MYP Agar Base [M636](#) / [M636S](#) / MYP HiVeg™ Agar Base [MV636](#) / Modified MYP Agar Base [M1139](#) / Modified MYP HiVeg™ Agar Base [MV1139](#) along with 25 ml / 50 ml Egg Yolk Emulsion [FD045](#) to make a total volume of 500 ml or to 500 ml of SDS Agar [M1155](#) / SDS HiVeg™ Agar [MV1155](#) / Salt Polymyxin Broth Base [M821](#) / [M821I](#) / Salt Polymyxin HiVeg™ Broth Base [MV821](#) / HiCrome™ Staph Agar Base, Modified [M1837](#) / Soyabean Casein Digest Medium Base [M011F](#) or to 1000 ml of Hicrome™ Bacillus Agar [M1651](#) . Mix well and pour into sterile petri plates / tubes.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.