



# Technical Data

## Bordetella Selective Supplement

FD004

An antibiotic supplement for the selective isolation of *Bordetella pertussis* .

### Composition

Per vial sufficient for 500 ml/ 1000 ml medium

### Ingredients

Cephalexin

### Concentration

20mg

### Directions:

Rehydrate the content of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Bordet Gengou Agar Base [M175](#) / [M175A](#) / Bordet Gengou HiVeg™ Agar Base [MV175](#) / [MV175A](#) or 1000 ml of sterile, molten Charcoal Agar Base w/Niacin [M1\\_053](#) / Charcoal HiVeg™ Agar Base w/Niacin [MV1053](#) together with 10% v/v defibrinated horse blood. Mix well and pour into sterile petri plates.

The vial content may be added to 500 ml of sterile half strength Charcoal Agar Base [M344](#) / Charcoal HiVeg™ Agar Base [MV344](#) with 10% v/v defibrinated horse blood for use as a transport medium for *Bordetella pertussis* .

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 1 / 2012

\* Not For Medicinal Use

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.