



# Technical Data

## Clostridium Difficile Supplement

FD010

An antibiotic supplement recommended for the selective isolation of *Clostridium difficile* .

### Composition

Per vial sufficient for 500 ml medium

#### \*Ingredients

D-Cycloserine

Cefoxitin

#### Concentration

250mg

8mg

### Directions:

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add along with 7% v/v defibrinated horse blood to 500 ml sterile, molten, cooled (45-50°C) Clostridium Difficile Agar Base [M836](#) / Clostridium Difficile HiVeg™ Agar Base [MV836](#) / Clostridium Brazier Agar Base [M1803](#) . Mix well and pour into sterile petri plates. Sheep blood may be used in place of horse blood but some strains of the organism will show a slightly reduced growth.

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label

Revision : 1 / 2012

\* Not For Medicinal Use

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.