



Technical Data

Listeria Selective Supplement II

FD063I

An antimicrobial supplement recommended by ISO Committee for the selective cultivation of *Listeria* species.

Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Acriflavin hydrochloride (Trypaflavin)	5mg
Nalidixic acid	20mg
Cycloheximide	25mg

Directions:

Rehydrate the contents of one vial with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml sterile, cooled (45-50°C) Listeria Selective Broth Base [M889](#) / Listeria Selective HiVeg™ Broth Base [MV889](#) / Buffered Listeria Enrichment Broth Base (BLE Broth Base) [M1578](#) / Buffered Listeria Enrichment HiVeg™ Broth Base (BLE HiVeg™ Broth Base) [MV1578](#) or to 500 ml sterile, molten, cooled (45-50°C) Listeria Selective Agar Base [M1474](#) Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.