



Technical Data

Tributyryn (10 ml per vial)

FD081

Recommended for detection of lipolytic microorganisms.

Composition

Per vial sufficient for 990 ml medium

Ingredients

Tributyryn

Concentration

10ml

Directions:

Warm up the refrigerated solution to room temperature. Add contents of one vial aseptically to 990 ml of Tributyrin Agar Base w/o Tributyrin [M157](#) / Tributyrin HiVeg™ Agar Base w/o Tributyrin [MV157](#) and then sterilize by autoclaving at 15 lbs pressure for 15 minutes. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.