



Technical Data

Tellurite - Cefixime Supplement

FD147

A selective supplement recommended by ISO Committee for the isolation of *Escherichia coli* 0157:H7 .

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Potassium tellurite
Cefixime

Concentration

1.250mg
0.025mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the contents to 495 ml of sterile, molten, cooled (45-50°C) MacConkey Sorbitol Agar Base [M2981](#) / HiCrome MacConkey Sorbitol Agar Base [M1340](#) . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.