



Technical Data

HiCrome EC 0157:H7 Selective Supplement

FD187

Recommended for selective isolation and easy detection of *Escherichia coli* O157:H7 from food samples.

Composition

Per vial sufficient for 1000 ml medium

*Ingredients

Novobiocin

Potassium tellurite

Concentration

10mg

1mg

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and add aseptically to 990 ml sterile, molten, cooled (45-50°C) HiCrome™ EC O157:H7 Selective Agar Base [M1575](#) / HiCrome™ EC O157:H7 Selective HiVeg™ Agar Base [MV1575](#) . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.