



Technical Data

Novobiocin Selective Supplement

FD290

An antibiotic supplement recommended for the selective isolation of Escherichia coli O157:H7 from foods.

Formula: (per vial, sufficient for 1000 ml medium)

Novobiocin	20.00 mg
------------	----------

Directions:

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water and mix well to dissolve. Aseptically add the rehydrated contents to 1000 ml of sterile, molten, cooled (45-50°C) Soyabean Bile Broth Base (M1286) / Soyabean HiVeg™ Broth Base (MV1286) / Modified Soyabean Bile Broth Base (M1286I) / Modified EC Broth Base (M1285). Mix well and pour into sterile petri plates.

* Not for medicinal use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.