



# Technical Data

## Lecithin solution

FD332

Lecithin solution is recommended for selective differentiation of *L.monocytogenes* and *Listeria* species based on PCPLC activity.

### Composition

Per vial sufficient for 480 ml medium

#### Ingredients

Soya lecithin

#### Concentration/vial

20ml

### Directions:

Thaw the content of one vial at room temperature. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) HiCrome L.mono Differential Agar Base (M2009) along with sterile contents of 1 vial of Modified L.mono Selective Supplement (FD333). Mix well and pour into sterile Petri plates.

### Storage and Shelf Life

On receipt product should be stored at -20°C. Use before the expiry date on the label.

Revision : 00 / 2016

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.