

## Calcium Caseinate Agar, Granulated

GM1309

Calcium Caseinate Agar, granulated is used for the detection and enumeration of proteolytic microorganisms in foodstuffs and other materials.

### Composition\*\*

<b>Ingredients</b>	<b>Gms / Litre</b>
Peptic digest of animal tissue	4.000
Meat extract	2.000
Casein enzymic hydrolysate	2.000
Calcium caseinate	3.500
Calcium chloride.2H <sub>2</sub> O	0.200
Tri-potassium citrate, H <sub>2</sub> O	0.350
Disodium hydrogen phosphate	0.105
Potassium dihydrogen phosphate	0.035
Sodium chloride	5.000
Agar	13.000
Final pH ( at 25°C)	7.0±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 30.19 grams in 1000 ml distilled water. Heat gently to boiling with frequent shaking for 10 minutes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix thoroughly while pouring into sterile Petri plates to suspend the precipitate. If desired, to increase turbidity, add 5-10 grams of skim milk powder before heating.

### Principle And Interpretation

Protein hydrolysis by microorganisms in foods may produce a variety of odour and flavour defects. On the other hand, microbial proteolytic activity may be desirable in certain foods such as in the ripening of cheese, where it contributes to the development of flavour, body and texture. In some foods the level of proteolytic microorganisms may be of value to predict refrigerated storage life and to assess processing methods (1, 2). Calcium Caseinate Agar is a modification of the original formulation of Frazier and Rupp (3) and is used for the detection and enumeration of proteolytic microorganisms in foodstuffs and other materials. Peptic digest of animal tissue, casein enzymic hydrolysate, meat extract provides nitrogenous, carbonaceous nutrients along with vitamins and amino acids. Phosphates are added to buffer the medium. Sodium chloride maintains the osmotic equilibrium. Casein in the medium is degraded by the proteolytic organisms. This results in formation of clear zones around the proteolytic colonies, in the otherwise opaque medium. Agar acts as a solidifying agent.

The test sample can be directly surface inoculated or the inoculation can be carried out by the pour plate technique. After an incubation for 24-48 hours, proteolytic organisms, if present will form clear zones on the medium. For better visualization of the zones, the plates can be flooded with 5-10% acetic acid.

### Quality Control

#### Appearance

Cream to yellow coloured granular medium

#### Gelling

Firm, comparable with 1.3% Agar gel.

#### Colour and Clarity of prepared medium

Whitish coloured, turbid gel forms in Petri plates

#### Reaction

Reaction of 3.02 % w/v aqueous solution at 25°C. pH : 7.0±0.2

#### pH

6.80-7.20

#### Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Proteolytic activity
<i>Bacillus cereus</i> ATCC 14579	50-100	good-luxuriant	≥70%	positive, clear zone surrounding colonies
<i>Escherichia coli</i> ATCC 25922	50-100	good-luxuriant	≥70%	negative, no clear zone surrounding colonies
<i>Proteus vulgaris</i> ATCC 13315	50-100	good-luxuriant	≥70%	negative, no clear zone surrounding colonies
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	good-luxuriant	≥70%	positive, clear zone surrounding colonies

### Storage and Shelf Life

Store below 30°C in tightly container and prepared medium at 2-8°C. Use before expiry period on the label.

### Reference

1. Chai T., Chen C., Rosen A. and Levin R. E., 1968, Appl. Microbiol., 16: 1738.
2. Martely F. G., Jayashankar S. R. and Lawrence, R. C., 1970, J. Appl. Bacteriol., 33: 363.
3. Frazier W. C. and Rupp P., 1928, J. Bacteriol., 16, 57-63

Revision : 00 / 2014



#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.