

**GRM3060 Rubidium chloride, A.R.**

**Product Number Packing**

---

GRM3060 : 25G  
GRM3060 : 5G

**Product Information**

---

Product Code : GRM3060  
Product Name : Rubidium chloride, A.R.  
Molecular Formula : R b C l  
Molecular Weight : 120.92  
CAS No. : 7791-11-9  
EC No. : 232-240-9  
HS Code : 2827 39 85  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : White crystals or powder  
Solubility : 100 mg soluble in 1 mL of water  
Magnesium (Mg) : <= 0.0005%  
Heavy metal (as Pb) : <= 0.0003%  
Iron (Fe) : <= 0.0003%  
Sulfate (SO4) : <= 0.005%  
Assay (AT) : min. 99.50 %

**Risk and Safety Information**

---

S-Phrase(s) : 22-24/25  
WGK : 2  
RTECS : VL8575000VL8575000  
Storage Temperature(°C) : Store below 30°C

**Transport Information**

---

Marine Pollutant : -  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods