

RM321 Canada balsam, Natural

Product Number Packing

RM321 : 500G

Product Information

Product Code : RM321
Product Name : Canada balsam, Natural
Synonym : Balsam canada , Natural
CAS No. : 8007-47-4
EC No. : 232-362-2
HS Code : 1301 90 00
Shelf Life : 4 years

Technical Specification

Appearance : Clear colourless to yellowish to greenish yellow viscous liquid.
Solubility : Miscible with alcohol, chloroform.
Density : 0.987 to 0.994
Refractive index (n 20/D) : 1.52 to 1.54

Risk and Safety Information

R-Phrase(s) : 10
S-Phrase(s) : 16
WGK : 2
RTECS : CP2352500
Flash Point(°F) : 143.6 °F
Flash Point(°C) : 62 °C
Storage Temperature(°C) : Store below 30°C.
Hazard Symbol(s)



Transport Information

UN No. : 1307
Class : 3
Packaging Group : 3
Marine Pollutant : No
ADR/RID : 1307 3/PG3
IMDG : 1307 3/PG3
IATA : 1307 3/PG3