

---

## Aluminium oxide active, neutral

**GRM4014**

### For column chromatography

#### Product Identifier

---

CAS No.	:	1344-28-1
EC No.	:	215-691-6
Molecular Formula	:	Al <sub>2</sub> O <sub>3</sub>
Molecular Weight	:	101.96
HS Code	:	2818 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

#### Technical Specification

---

Appearance	:	White to off-white powder
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	5.00 - 8.00
Particle size 1	:	min. 80.00% passing through 70 mesh (210 µm)
Particle size 2	:	max. 20.00% passing through 230 mesh (63 µm)

#### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	BD1200000
WGK	:	NWG