
Magnesium metal ribbon

GRM4339


Product Identifier

CAS No.	:	7439-95-4
EC No.	:	231-104-6
Atomic Formula	:	Mg
Atomic Weight	:	24.31
HS Code	:	8104 30 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Silvery white metal ribbon
Solubility	:	Reacts readily with dilute acids
Assay (EDTA Titration)	:	99.00 - 103.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H228- H260
Precautionary Statement(s)	:	P210- P223- P231+P232- P370 + P378- P422
UN No.	:	1869
Class	:	4.1
Packing Group	:	III
RTECS	:	OM2100000
WGK	:	NWG