

GRM6231 di-Sodium tartrate dihydrate**Product Number Packing**

GRM6231 : 500G

Product Information

Product Code : GRM6231
Product Name : di-Sodium tartrate dihydrate
Synonym : Sodium tartrate dibasic dihydrate; L-(+)-Tartaric acid disodium salt
Molecular Formula : $C_4H_4Na_2O_6 \cdot 2H_2O$
Molecular Weight : 230.08
CAS No. : 6106-24-7
EC No. : 212-773-3
HS Code : 2918 13 00
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder or needles
Solubility : 100 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 7.00 - 9.00
Chloride (Cl) : $\leq 0.005\%$
Iron (Fe) : $\leq 0.001\%$
Loss on drying (at 150°C, 4hr) : 15.61 - 15.71 %
Assay (NT, on dry basis) : 99.00 - 101.00%

Risk and Safety Information

WGK : 3
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods