



Product Information

LA016 MetaLoop SL

Product Number

Packing

LA016	: 1 X 8 No.
LA016	: 1 X 24 No.

Product Information

Product Code	: LA016
Product Name	: MetaLoop SL

Technical Specification

Description	: Fixed straight nichrome loop embedded in SS Rod, with heat resistant handle. Useful for stab culturing.
-------------	---

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods