

Metaloops

LA650

An inoculation loop, simple tool used for streaking. The loop of wire at the tip is made of nichrome and handle with stainless steel.

Application : Microbiology laboratory, Molecular biology laboratory, Culture preparation and Various other laboratories.

Product Name	Product Code	Description	Dia (mm)
Metaloop -SS-2	LA650	Fixed nichrome loop embedded in SS rod with heat resistant handle. Calibrated to 0.005 ml.	2.20
Metaloop -SS-4	LA016	Fixed nichrome loop embedded in SS rod with heat resistant handle. Calibrated to 0.01 ml.	4.00

Product Features

- Heat Resistant handle.
- Pack Size : 1X8 Nos. & 1X24 Nos.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.