



## Violet Red Bile Agar

M049A

Violet Red Bile Agar is recommended for selective isolation and enumeration of coli-aerogenes bacteria in water, milk, other dairy food products and also from clinical samples.

### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	7.000
Yeast extract	3.000
Lactose	10.000
Bile salts mixture	1.500
Sodium chloride	5.000
Neutral red	0.030
Crystal violet	0.002
Agar	12.000
Final pH ( at 25°C)	7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 38.53 grams in 1000 ml distilled water. Heat with stirring to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45°C and pour into sterile Petri plates containing the test sample.

### Principle And Interpretation

Violet Red Bile Agar, a modification of MacConkeys original formulation (1) is used for the enumeration of coli-aerogenes bacterial group. It relies on the use of the selective inhibitory components crystals violet and bile salts and the indicator system lactose, and neutral red. Thus, the growth of many unwanted organisms is suppressed, while tentative identification of sought bacteria can be made. Organisms, which rapidly attack lactose, produce purple colonies surrounded by purple halos. Non-fermenters or late lactose-fermenters produce pale colonies with greenish zones (2). VRBA is recommended by APHA (3, 4). Violet Red Bile Agar (1.2 % Agar) (M049A) is prepared, in accordance with the ISO Committee (9).

Selectivity of VRBA can be increased by incubation under anaerobic conditions and/ or at elevated temperature, i.e. equal to or above 42°C (5-7).

Peptic digest of animal tissue and yeast extract serve as sources of carbon, nitrogen, vitamins and other essential growth nutrients. Lactose is the fermentable carbohydrate, utilization of which leads to the production of acids. Neutral red indicator detects the acidity so formed. Crystal violet and bile salts mixture help to inhibit the accompanying gram-positive and unrelated flora. Sodium chloride maintains the osmotic equilibrium. Violet Red Bile Agar is not completely specific for enterics; other accompanying bacteria may give the same reaction. Further biochemical tests are necessary for positive identification (8).

### Quality Control

#### Appearance

Light yellow to pink homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.2% Agar gel.

#### Colour and Clarity of prepared medium

Reddish purple coloured clear to slightly opalescent gel forms in Petri plates.

#### Reaction

Reaction of 3.85% w/v aqueous solution at 25°C. pH : 7.4±0.2

#### pH

7.20-7.60

#### Cultural Response

M049A: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<b>Cultural Response</b>				
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	>=50%	pink to pinkish red
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	>=50%	pinkish red with bile precipitation
<i>Salmonella Enteritidis</i> ATCC 13076	50-100	luxuriant	>=50%	Colourless to orangish yellow
<i>Staphylococcus aureus</i> ATCC 25923	>=10 <sup>3</sup>	inhibited	0%	

## Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

## Reference

1. MacConkey A., 1905, J. Hyg., 5, 333-379
2. Corry J. E. L., Curtis G. D. W. and Baird R. M., (Ed.), 1995, Culture Media for Food Microbiology, Vol. 34, Progress in Industrial Microbiology, Elsevier, Amsterdam.
3. Downes F. P. and Ito K., (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
4. Marshall R. T., (Ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th Ed., APHA, Washington, D. C.
5. Mossel D. A. A. and Vega C. L., 1973, Hlth. Lab. Sci., 11:303
6. Mossel D. A. A., Eclerink I., Koopmans M. and Van Rossem F., 1979, Food Protect., 42 : 470
7. Mossel D. A. A. et al, 1986, J. Appl. Bacteriol., 60:289.
8. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.
9. International Organization for Standardization (ISO), 1991, Draft ISO/DIS 4382

Revision : 02 / 2015

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.