



Technical Data

Potato Dextrose Agar w/ 2% Agar

M096F

Intended Use:

Recommended for isolation and enumeration of yeasts and moulds from dairy and other food products in accordance with FDA BAM, 1998.

Composition**

Ingredients	Gms / Litre
Potatoes, infusion from	200.000
Dextrose (Glucose)	20.000
Agar	20.000
Final pH (at 25°C)	5.6±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 44.0 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well before dispensing. In specific work, when pH 3.5 is required, acidify the medium with sterile 10% tartaric acid. The amount of acid required for 100 ml. of sterile, cooled medium is approximately 1 ml. Do not heat the medium after addition of the acid.

Principle And Interpretation

A large group of fungi found in the nature can spoil the foods and at times turns to be pathogenic to humans. This can also lead to economic losses at producer and consumer levels, apart from its toxic effects. Dilution plating and direct plating are the most common methods used in the isolation and enumeration of fungi. The direct plating is more efficient than the dilution plating method for detecting individual mold species including most of the toxin producers, but it is less effective in detecting yeasts. Potato Dextrose agar w/2% agar is recommended by FDA BAM (1) for plate counts of yeasts and moulds in the examination of foods and dairy products.

The media has also been used for stimulating sporulation, for maintaining stock cultures and for differentiation of typical varieties of dermatophytes on the basis of pigment production (2). Potato infusion and dextrose promote luxuriant fungal growth. Adjusting the pH of the medium by tartaric acid to 3.5 inhibits the bacterial growth. Heating the medium after acidification should be avoided as it may hydrolyse the agar which can render the agar unable to solidify.

Add 25-50g of the sample into 0.1% peptone water to achieve 10-1 dilution and homogenize in a stomacher for 2 min or blend it for 30-60 sec. Make appropriate dilutions in 0.1% peptone water. In normal conditions, dilutions up to 10⁻⁶ shall satisfy the needs. Add 0.1 ml into Dichloran Glycerol Medium Base w/Rose Bengal/ Dichloran (M1000) during spread plate technique; while add 1ml into Dichloran Glycerol Medium Base (M1129) during pour plate technique. Incubate at 25°C and count the colonies. Individual colonies from these plates can be isolated using Potato Dextrose agar w/2% agar (M096F) and species identification can be done if necessary(1).

Type of specimen

Hqqf"uc o rngu"

Specimen Collection and Handling:

Hqt"hqf" uc o rngu."hqmqy"crtrtrtkcvg"vgejpkswgu"ht"uc o rng"eqnngewkqp"cpf"rtqeguukpi"cu"rgt"i wkgfngkpgu"*1,2,8+0
Chvgt"wg."eqpvc o kpcvgf" o cvgtkcu" o wuv"dg"uvgtknk |gf"d{"cwwqencxkpi"dghqtg"fkuectfkpi 0"

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Heating the medium after acidification should be avoided as it may hydrolyse the agar which can render the agar unable to solidify.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the gzkrt{"rgtkqf"when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% Agar gel

Colour and clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

pH of 4.4% w/v aqueous solution at 25°C. pH : 5.6±0.2

pH

5.40-5.80

Cultural Response

Cultural characteristics observed after an incubation at 25°C for 5-7 days.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Candida albicans</i> ATCC 10231 (00054*)	50 -100	luxuriant	≥70 %
<i>Saccharomyces cerevisiae</i> ATCC 9763 (00058*)	50 -100	luxuriant	≥70 %

Mg{"<,Eqttgurqpfkpi"Y FEO"pw o dgtu0

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

Reference

1. FDA, U.S. 1998. Bacteriological Analytical Manual. 8 ed. Gaithersburg, MD: AOAC International.
2. MacFaddin, J. F. 1985. Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria vol. 1. Baltimore: Williams and Wilkins.

Revision : 1 / 2011

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.