



Iron Sulphite Agar Modified

M1852I

Iron Sulphite Agar Modified is used for the enumeration of sulphite-reducing bacteria growing under anaerobic conditions.

Composition**

| Ingredients | Gms / Litre |
|------------------------------------|-------------|
| Enzymatic digest of casein | 15.000 |
| Pancreatic digest of soyabean meal | 5.000 |
| Yeast extract | 5.000 |
| Disodium disulfite | 1.000 |
| Iron (III) Ammonium citrate | 1.000 |
| Agar | 15.000 |
| Final pH (at 25°C) | 7.6±0.2 |

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 42 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Dispense as desired.

Principle And Interpretation

Iron Sulphite Agar, Modified is recommended by ISO for the enumeration of sulphite reducing bacteria.(1). Most Clostridia possess sulfite reductase in their cytoplasm but they are unable to expel them to the exterior. So when H₂S is produced from sulfite, the colony becomes dark due to the formation of precipitates of iron sulfide from citrate.

Enzymatic Digest of Casein and pancreatic digest of soyabean meal provide nitrogen, vitamins, minerals and amino acids necessary for the growth of organism. Yeast extract serves as a rich reservoir of vitamins especially B-complex vitamins. Ferric citrate ammonium citrate and Disodium sulfite serves as are H₂S indicators, wherein *Clostridium perfringens* reduces the sulfite to sulfide which in turn reacts with the iron and forms a black iron sulfide precipitate, seen as black colonies. Agar is the solidifying agent.

Enumeration with this medium can be performed using either tubes or plates. In case of tubes distribute 20-25 ml of the medium in tubes and inoculate 1 ml of test sample or 1 ml of serial dilutions of 10⁻¹ and 10⁻² in molten state. Allow to solidify, and pour 2-3 ml of the same medium in each tube to overlay. In case of Petri plates, transfer 1 ml of test sample or initial dilution. Further dilution can be carried out and 1 ml of each dilution (10⁻¹ and 10⁻²) is transferred to an empty Petri plate. Cool the medium to 44-47°C and pour 15-20 ml of the medium to the Petri plate containing the inoculum. Mix the inoculum and allow the medium to solidify. Overlay the medium with 5-10 ml of the same medium.

After solidification, incubate the medium at 36-38°C for 24-48 hours. If thermophilic bacteria are suspected, a second of tubes is incubated at 49-51°C for 24-48 hours. After incubation, black coloured colonies, possibly surrounded by a black zone are counted as sulphite reducing bacteria.

Quality Control

Appearance

Light yellow to brownish yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured, slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.2% w/v aqueous solution at 25°C. pH : 7.6±0.2

pH

7.40-7.80

Cultural Response

Cultural characteristics observed under anaerobic conditions, after an incubation at 36-38°C for 24-48 hours.

Cultural Response

| Organism | Inoculum | Growth | Recovery | Colour of colony |
|--|----------|-----------|----------|------------------|
| Cultural Response <i>Clostridium botulinum</i> ATCC 25763 | 50-100 | luxuriant | >=50% | black |
| <i>Clostridium butyricum</i> ATCC 13732 | 50-100 | luxuriant | >=50% | black |
| <i>Clostridium sporogenes</i> ATCC 19404 | 50-100 | luxuriant | >=50% | black |
| <i>Desulfotomaculum</i> <i>nigrificans</i> ATCC 19998 | 50-100 | luxuriant | >=50% | black |
| <i>Escherichia coli</i> ATCC 25922 | 50-100 | good | 40-50% | no blackening |

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

Reference

1. Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of sulphite reducing bacteria growing under anaerobic conditions, ISO 15213.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.