



Technical Data

Rapid Perfringens Medium Base (Twin pack)

M1898

For rapid detection of *Clostridium perfringens* in food.

Composition**

Ingredients	Gms / Litre
Part A	-
Litmus milk powder	70.000
Part B	-
Pancreatic digest of casein	15.13
Yeast extract	8.04
Glucose	10.55
Sodium chloride	4.02
L-Cystine	0.51
Sodium thioglycollate	0.51
Resazurin sodium	0.001
Gelatin	60.000
Peptone	5.000
Dipotassium hydrogen orthophosphate	5.000
Iron (II) sulphate	0.500
Agar	0.755
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 70 grams of Part A in 500 ml distilled water. Mix well and adjust the pH to 6.8. Sterilize by autoclaving at 15 lbs pressure (121°C) for 5 minutes. Cool to 45-50°C and aseptically add sterile rehydrated contents of one vial of Perfringens Selective Supplement (FD307).

Suspend 110 grams of Part B in 500 ml distilled water. Heat if necessary to dissolve the medium completely. Adjust the pH to 7.1. Dispense 5 ml amount in screw-capped glass tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add 5 ml of previously cooled Part A solution to Part B. Mix well and store at 2-8°C. Before use, liquefy the medium by placing the tubes in a water bath at 45-50°C for 30 minutes.

Principle And Interpretation

Rapid Perfringens Medium Base is formulated by Erickson & Deibel (1). The Mesophilic spore forming anaerobes belonging to the genus *Clostridia* of food concern are Gram-positive, catalase negative, rods of varying sizes.

The medium can be used to initiate growth from small inocula and to obtain the highest viable count of *Clostridia*. Rapid Perfringens Medium Base is a liquid medium with a litmus milk base and is prepared in tubes. Selectivity is provided by the antibiotics Polymyxin B sulfate and neomycin sulfate, coupled with an incubation temperature of 46°C(2).

Quality Control

Appearance

homogeneous free flowing powder

Colour and Clarity of prepared medium

Light Brown coloured Light brown opaque solution in tubes.

Reaction

Reaction of 7.0% w/v of Part A + 11.0% w/v of Part B at 25°C. pH : 6.80±0.2

Cultural Response

Cultural characteristics observed in an anaerobic atmosphere after an incubation at 46°C for 48 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth
Cultural Response <i>Clostridium perfringens</i> ATCC 13124	50-100	good - luxuriant with stormy fermentation
<i>Proteus mirabilis</i> ATCC 25933	50-100	good - luxuriant

Storage and Shelf Life

Reference

1. Erickson, J.E. and Deibel, R.H. (1978) New medium for rapid screening and enumeration of *Clostridium perfringens* in foods. *Appl. Environ. Microbiol.* 36, 567-571. 2. Handbook of Culture Media for Food and Water Microbiology <(,<)> 3rd edition. Edited by Janet E.L. Corry, Gordon D.W. Curtis and Rosamund M. Baird.

Revision : 1 / 2011



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.