



## Phenylalanine Agar

M281

### Classification

EXb` ` XaVWV Ye WXXeTgba` bY *Proteus` TaW Providencia` Zehc` bY beZTaV` f` Yb` ` bgXe` ` X UX` bY Enterobacteriaceae` ba`g XUIff` bYg XeTU`g` g`Ye` `c[Xal`\_cl` di` VTWVb` `c[Xal`\_TTa`aX`*

### Composition\*\*

Ingredients	Gms / Litre
Yeast extract	3.000
Sodium chloride	5.000
DL-Phenylalanine	2.000
7`fbWh` `1`VbZa`c` bfc` TgX`	1.000
Agar	15.000
Final pH ( at 25°C)	7.3±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 26 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense in tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. 6bb`g` (Z`#6!` Allow the tubed medium to cool in a slanting position.

### Principle And Interpretation

The ability of *Proteus* species to convert phenylalanine to phenylpyruvic acid is an important reaction in the differentiation of *Enterobacteriaceae* (1, 2). Based on this criterion, Buttiaux developed Phenylalanine Agar for differentiation of *Proteus* and *Providencia* group from other members of *Enterobacteriaceae* (3, 4) by the ability of organism in the genera within *Proteus* to deaminate phenylalanine. Phenylalanine Agar is the modification of the medium originally developed by Ewing et al (5).

Yeast extract in the medium supports the growth of the organisms. Sodium chloride maintains osmotic equilibrium. The phenylalanine serves as the substrate for enzymes, which are able to deaminate it to form phenylpyruvic acid. A recommended technique is to inoculate the slant surface with plenty of inoculum and incubate it for 12-16 hours. After incubation, add 0.2 ml of 10% ferric chloride solution so that the solution floods all over the growth. The addition of (0.2 ml 3-5 drops) of a 10% aqueous ferric chloride solution (or a 12% aqueous ferric chloride solution acidified with 2.5 ml of concentrated HCl per 100 ml of reagent) to the cultures following incubation results in the appearance of a light to deep green color (positive reaction) or no color change (negative reaction). In a positive reaction, any phenylpyruvic acid present will react with the ferric salt in the reagent to give a green color. Interpret the results within 5 minutes upon addition of reagent as the green colour fades quickly (1, 4).

### Gl`cXbYfcX` Xa`

Rwtg`kuqncvgu`

### Specimen Collection and Handling:

A recommended technique is to inoculate the slant surface with plenty of inoculum and incubate it for 12-16 hours. After incubation, add 0.2 ml of 10% ferric chloride solution so that the solution floods all over the growth. The addition of (0.2 ml 3-5 drops) of a 10% aqueous ferric chloride solution (or a 12% aqueous ferric chloride solution acidified with 2.5 ml of concentrated HCl per 100 ml of reagent) to the cultures following incubation results in the appearance of a light to deep green color (positive reaction) or no color change (negative reaction). In a positive reaction, any phenylpyruvic acid present will react with the ferric salt in the reagent to give a green color. Interpret the results within 5 minutes upon addition of reagent as the green colour fades quickly (1, 4).

Chvgt`wug.`eqpvc`o`kpcv`gf` `o`cvgtk`cnu` `o`wuv`dg`uvgtknk`|`gf`d`{"`cwwqencxkpi`"dghqtg`"fkuectfkpi`0`

## J Te\azTaVceWIngbaF-

EXIWg X\_TUX`UXbeXbcXa\az`g XVbag`aX` J XTē cbgXgi XZ`bi X"cbgXgi XVbg`az`X XcbgXgba" YVX`cbgXgba!  
 9b\_bj` ZbbW` \WbUb\_bZVT`\_JU`cēVgV` j [\X [TaWaZ`fcX\`Xaf` TaW`Vh`gēX` FgTaWceWIngbaF` If`cXē  
 XgUf[XWZh\X\ax`f[bh\_WUX`Y\_bj`XWj`[\X [TaWaZ`ucX\`Xaf`FTYg`Zh\X\ax`` Il`UXēXēXWa`aW`MIL`  
 fTYg`Wg`f[Xg`)

?\` \gTgbaF

Š\Uqog`qticpluo"oc{"ujqy"rqqt"itqyvj"fwg"vq"pwvtkvkqpen"xctkēvkqp0

### Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

### Quality Control

**Appearance**

Cream to yellow homogeneous free flowing powder

**Gelling**

Firm, comparable with 1.5% Agar gel

**Colour and Clarity of prepared medium**

Light amber coloured slightly opalescent gel forms in tubes as slants

**Reaction**

Reaction of 2.6% w/v aqueous solution at 25°C. pH : 7.3±0.2

**pH**

7.10-7.50

**Cultural Response**

Cultural characteristics observed after an incubation at 35-37°C for 12-16 hours

Organism	Inoculum (CFU)	Growth	Phenylalanine deaminase
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	negative reaction
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	negative reaction
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	positive reaction, green colouration after addition of 10% ferric chloride
<i>Proteus vulgaris</i> ATCC 13315	50-100	luxuriant	positive reaction, green colouration after addition of 10% ferric chloride
<i>Providencia alcalifaciens</i> ATCC 9886	50-100	luxuriant	positive reaction, green colouration after addition of 10% ferric chloride

## Reference

1. Singer J. and Volcani B. E., 1955, J. Bacteriol., 69:303.
2. Henrikson S. D., 1950, J. Bacteriol., 60:225.
3. Buttiaux R., Osteux R., Fresnoy R. and Moriamez J., 1954, Ann. Inst. Pasteur Lille., 87:375.
4. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams & Wilkins, Baltimore, Md.
5. Ewing W. H., Davis B. R. and Reavis R. W., 1957, Public Health Lab., 15:153.

Revision : 1 / 2011

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.