

Triple Sugar Iron HiCynth™ Agar

MCD021

Triple Sugar Iron HiCynth™ Agar is used for the identification of gram-negative enteric bacilli on the basis of dextrose, lactose and sucrose fermentation and hydrogen sulphide production.

Composition**

Ingredients	Gms / Litre
HiCynth™ Peptone No.1*	16.000
HiCynth™ Peptone No.6*	3.000
Lactose	10.000
Sucrose	10.000
Dextrose	1.000
Sodium chloride	4.000
Sodium thiosulphate	3.000
Ferric sulphate	0.200
Ferric ammonium citrate	0.300
Histidine	2.000
DL-Methionine	2.000
Arginine	1.000
Phenol red	0.024
Agar	12.000
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

*Chemically defined peptones

Directions

Suspend 64.52 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Mix well and distribute into test tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Allow the medium to set in sloped form with a butt about 1 inch long.

Note: For better results, the medium can be sterilized by autoclaving at 10 lbs pressure (115°C) for 15 minutes.

Principle And Interpretation

Triple Sugar Iron Agar was originally proposed by Sulkin and Willett (1) and modified by Hajna (2) for identifying *Enterobacteriaceae*. This medium complies with the recommendation of APHA, for the examination of meat and food products (3), for the examination of milk and dairy products (4) and for microbial limit test for confirming the presence of *Salmonellae* (5, 6) and in the identification of gram-negative bacilli (6,7). Triple Sugar Iron HiCynth™ Agar is the modification using chemically defined peptones free from animal and vegetable peptones to avoid BSE/TSE risks associated with animal peptones.

HiCynth™ Peptone No.1 and HiCynth™ Peptone No.6 provide nitrogenous and carbonaceous compounds, sulphur, trace elements and vitamin B complex etc. Sodium chloride maintains osmotic equilibrium. Lactose, sucrose and dextrose are the fermentable carbohydrates. Sodium thiosulphate and ferrous ions make H₂S indicator system. Phenol red is the pH indicator. Organisms that ferment glucose produce a variety of acids, turning the colour of the medium from red to yellow. More amount of acids are liberated in butt (fermentation) than in the slant (respiration). Growing bacteria also form alkaline products from the oxidative decarboxylation of peptone and these alkaline products neutralize the large amounts of acid present in the butt. Thus the appearance of an alkaline (red) slant and an acid (yellow) butt after incubation indicates that the organism is a glucose fermenter but is unable to ferment lactose and/or sucrose. Bacteria that ferment lactose or sucrose (or both), in addition to glucose, produce large amounts of acid enables no reversion of pH in that region and thus bacteria exhibit an acid slant and acid butt. Gas production (CO₂) is detected by the presence of cracks or bubbles in the medium, when the accumulated gas escapes. Thiosulphate is reduced to hydrogen sulphide by several species of bacteria and H₂S combines with ferric ions of ferric salts to produce the insoluble black precipitate of ferrous sulphide. Reduction of thiosulphate proceeds only in an acid

environment and blackening usually occurs in the butt of the tube. Triple Sugar Iron HiCynth™ Agar should be used in parallel with Urea Agar / Broth to distinguish between *Salmonella* and *Proteus* species. The reactions can be summarized as follows:

Alkaline slant / acid butt-only glucose fermented

Acid slant / acid butt-glucose and sucrose fermented or glucose and lactose fermented or all the three sugars, glucose, lactose and sucrose fermented.

Bubbles or cracks present-gas production

Black precipitate present-H₂S gas production

Some members of the *Enterobacteriaceae* and H₂S producing *Salmonella* may not be H₂S positive on TSI HiCynth™ Agar. Some bacteria may show H₂S production on Kligler Iron HiCynth™ Agar but not on TSI HiCynth™ Agar. This can happen because utilization of sucrose in TSI HiCynth™ Agar suppresses the enzymic pathway that result in H₂S production.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel.

Colour and Clarity of prepared medium

Pinkish red coloured clear to slightly opalescent gel forms in tubes as slants.

Reaction

Reaction of 6.45% w/v aqueous solution at 25°C. pH : 7.4±0.2

pH

7.20-7.60

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth	Slant	Butt	Gas	H ₂ S
Cultural Response						
<i>Citrobacter freundii</i> ATCC 8090	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, yellowing of the medium	positive reaction	positive, blackening of medium
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, yellowing of the medium	positive reaction	negative, no blackening of medium
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, yellowing of the medium	positive reaction	negative, no blackening of medium
<i>Klebsiella pneumoniae</i> ATCC 13883	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, yellowing of the medium	positive reaction	negative, no blackening of medium
<i>Proteus vulgaris</i> ATCC 13315	50-100	luxuriant	alkaline reaction, red colour of the medium	acidic reaction, yellowing of the medium	negative reaction	positive, blackening of medium
<i>Salmonella Paratyphi A</i> ATCC 9150	50-100	luxuriant	alkaline reaction, red colour of the medium	acidic reaction, yellowing of the medium	positive reaction	negative, no blackening of medium
<i>Salmonella Typhi</i> ATCC 6539	50-100	luxuriant	alkaline reaction, red	acidic reaction, yellowing of the medium	negative reaction	positive, blackening of medium

<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	colour of the medium alkaline reaction, red colour of the medium	acidic reaction, positive yellowing of the medium	positive, blackening of medium
<i>Shigella flexneri</i> ATCC 12022	50-100	luxuriant	alkaline reaction, red colour of the medium	acidic reaction, negative yellowing of the medium	negative, no blackening of medium
<i>Escherichia coli</i> ATCC 8739	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, positive yellowing of the medium	negative, no blackening of medium
<i>Escherichia coli</i> NCTC 9002	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, positive yellowing of the medium	negative, no blackening of medium
<i>Klebsiella pneumoniae</i> ATCC 10031	50-100	luxuriant	acidic reaction, yellowing of the medium	acidic reaction, positive yellowing of the medium	negative, no blackening of medium

Storage and Shelf Life

Store below 30°C in a tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Sulkin E.S. and Willett J.C., 1940, J. Lab. Clin. Med., 25:649.
2. Hajna A.A., 1945, J. Bacteriol, 49:516.
3. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
4. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
5. Finegold S. M. and Baron E. J., 1986, Bailey and Scotts Diagnostic Microbiology, 7th Ed., The C.V. Mosby Co., St. Louis.
6. Eaton A. D., Clesceri L. S. and Greenberg A. W., (Eds.), 2012, Standard Methods for the Examination of Water and Wastewater, 22nd Ed., APHA, Washington, D.C.
7. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.

Revision : 00 / 2015



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.